

# Sampling Accessories

## Integrated Sampling Systems

These Integrated Sampling Systems are direct-attach cuvette holder and light source combinations created specifically for our USB2000+ and USB4000 Spectrometers (page 14). Both systems receive power and control signals through a connector on the spectrometer.

### USB-ISS-UV-VIS

#### Integrated Sampling System

The USB-ISS-UV-VIS Integrated Sampling System is a direct-attach sample holder and deuterium tungsten halogen light source (200-1100 nm) combination for 1-cm square cuvettes. The USB-ISS-UV-VIS allows you to adjust the intensity of the bulb via software. The sampling system has an electronic shutter for taking dark measurements and comes with a 5-volt power supply.

### USB-ISS-VIS

#### Integrated Sampling System

The USB-ISS-VIS Integrated Sampling System has a violet LED-boosted tungsten source (390-900 nm) and a sample holder that bolts to the front of a USB4000 or USB2000+ Spectrometer. The spectrometer provides the power and control signals for the light sources. The USB-ISS-VIS holds 1-cm cuvettes.

### ISS-UV-VIS

#### Integrated Sampling System

The ISS-UV-VIS Integrated Sampling System is a combination RF deuterium source with a tungsten bulb in a housing connected to a holder for 1-cm cuvettes. This sampling system couples to an Ocean Optics spectrometer with optical fiber to create a small-footprint system for relative absorbance. This sampling system is best used with Ocean Optics' 300  $\mu$ m solarization-resistant optical fiber.

### ISS-2

#### Integrated Sampling System

The ISS-2 Integrated Sampling System is a fully integrated 1-cm cuvette holder and tungsten halogen light source for relative absorbance measurements. It couples to Ocean Optics spectrometers with optical fiber to create a small-footprint system for VIS-NIR (~360-1100 nm) applications.



| Specifications              | USB-ISS-UV-VIS             | USB-ISS-VIS                   | ISS-UV-VIS                 | ISS-2                  |
|-----------------------------|----------------------------|-------------------------------|----------------------------|------------------------|
| Dimensions (mm):            | 198 x 105.1 x 40.6         | 40.7 x 88.8 x 34.1            | 198 x 104.9 x 40.9         | 155 x 50 x 53.3        |
| Weight:                     | 200 g                      | 130 g                         | 400 g                      | 240 g                  |
| Power consumption:          | 1.8 A @ 5 VDC              | 160 mA @ 5 VDC                | 420 mA @ 12 VDC            | 600 mA @ 12 VDC        |
| Wavelength range (source):  | ~200-1100 nm (Typical)     | 390-900 nm (Typical)          | ~200-1100 nm (Typical)     | ~360-1100 nm (Typical) |
| Pathlength:                 | 1 cm                       | 1 cm                          | 1 cm                       | 1 cm                   |
| Cuvette shape:              | Square                     | Square                        | Square                     | Square                 |
| Light source:               | Deuterium tungsten         | Tungsten and violet LED       | Deuterium tungsten         | Tungsten               |
| Bulb life (hours):          | 800 (deut.); 2,000 (tung.) | 2000 (tungsten); 45,000 (LED) | 800 (deut.); 2,000 (tung.) | 900                    |
| Time to stabilized output:  | ~30 minutes                | ~5 minutes                    | ~30 minutes                | ~30 minutes            |
| Filter slot:                | None                       | None                          | None                       | 6.35 mm                |
| Recommended optical fibers: | None                       | None                          | QP400-025-SR               | QP400-2-UV-VIS         |
| Spectrometers:              | USB2000+ and USB4000       | USB2000+ and USB4000          | All                        | All                    |
| "Z" dimension:              | 15 mm                      | 15 mm                         | 15 mm                      | 15 mm                  |