IKA® MS 1**IKA® MS 2**

<i>BETRIEBSANLEITUNG</i>	<i>DE</i>	3
<i>OPERATING INSTRUCTIONS</i>	<i>EN</i>	11
<i>MODE D'EMPLOI</i>	<i>FR</i>	19
<i>INDICACIONES DE SEGURIDAD</i>	<i>ES</i>	23
<i>VEILIGHEIDSIINSTRUCTIES</i>	<i>NL</i>	24
<i>NORME DI SICUREZZA</i>	<i>IT</i>	26
<i>SÄKERHETSANVISNINGAR</i>	<i>SV</i>	27
<i>SIKKERHEDSHENVISNINGER</i>	<i>DA</i>	28
<i>SIKKERHEDSHENVISNINGER</i>	<i>NO</i>	29
<i>TURVALLISUUSOHJEET</i>	<i>FI</i>	30
<i>γποδείξεις ασφάλειας</i>	<i>EL</i>	31
<i>INSTRUÇÕES DE SERVIÇO</i>	<i>PT</i>	33

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Contents

	Page
Declaration of CE-conformity	9
Safety instructions	9
Correct use	10
Unpacking	10
General information	11
Commissioning	11
Switching on and setting of the speed	
Operation mode	
Attachments	
Supported weight and speed range	
Maintenance and Cleaning	13
Accessories	13
Warranty	14
Technical data	14
List of spare parts MS 1 and MS 2	21
Spare parts diagram MS 1 and MS 2	22

Declaration of CE-conformity

We declare under our sole responsibility that this product corresponds to the regulations 89/336/EG and 73/23/EG and conforms with the standards or standardized documents DIN EN IEC 61010-1 and DIN EN IEC 61326-1.

Safety instructions

- In order to be able to use the appliance properly and safely, every user must first read the operating instructions and observe the safety instructions contained therein.
- Take care of these operating instructions and keep them in a place

where they can be accessed by everyone.

- Only staff who have been trained accordingly, know the appliance and are authorised to carry out work in this field should use this appliance.
- Please follow the relevant safety instructions and guidelines, and occupational health and safety regulations for use in the laboratory.

! Note! Covering or parts that are capable of being removed from the unit without accessory equipment have to be reattached to the unit for safe operation in order to prevent, for example, the ingress of fluids, foreign matter, etc..

- When working with the shaking unit, the user must select and wear his personal protective equipment according to the mixing hazard category of the medium. Defective or inappropriate protective equipment can expose the user to the risk of spurting liquids, projectile parts or being pulled in at the shaking table or support.

- Never touch moving parts (risk of crushing, impact and cutting, see fig.1: point of danger)!

- Ensure that parts of the body, hair or items of clothing cannot be trapped by the moving parts.

- Always install the appliance on a flat, stable non-slip base.

- Before use, please ensure that the knob for setting the speed is at the left stop as the appliance starts up at the last speed which was set.

- Pay attention to the vessels on the shaking table when setting the speed. This will prevent any of the medium to be shaken from spurting out of the sample vessels.

- All accessories and vessels in place for the shaking process must be firmly secured. Shaking vessels which are not properly secured could get damaged or be projected out, thus causing injury. It is essential to regularly check that the vessels to be shaken and the attachments are firmly secured, especially before using the appliance again.

- If you notice that the device is not running smoothly, the speed


Fig. 1: point of danger

must always be reduced until no more unevenness occurs in the operation.

- Please make certain that individual agitation containers are placed in the middle of the agitation table and that if several containers are placed on the agitation table together, they are evenly distributed and well secured.
- Because of improper loading at the position of the center of gravity, dynamic forces may arise during the shaking process that cause the shaker to move about the table. For limits on the load or the weight placed on the agitation table at high speed, please refer to the table in the chapter „Commissioning“.
- The appliance starts up again automatically following a cut in the power supply during a work process.
- The user may also be at risk if inflammable materials are used in the work process or if the transmitted mechanical power causes any glass equipment in use to break.
- The device may only be used to process those materials or material mixtures that the user knows will not react dangerously to the extra energy produced by processing. This also applies to extra energy generated by means of solar radiation during the work process.
- The shaking unit may not be used in explosive atmospheres, for mixing dangerous substances or under water.
- Accessories may only be assembled once the plug has been disconnected.
- The safety of the user cannot be guaranteed if the appliance is operated with accessories that are not supplied or recommended by the manufacturer or if the appliance is operated improperly contrary to the manufacturer's specifications.

The machine may only be opened by trained specialists - even during repairs. The machine is to be unplugged from the mains before opening. Live parts inside the machine may still be live for some time after unplugging from the mains.

Correct use

The shakers **MS 1** and **MS 2** are suitable for usage in various attachments for mixing liquids in e.g. bottles, flasks, test tubes, microtiter plates (only **MS 1**) and bowls for a maximum supported load of 0.5 kg (incl. the suitable attachment).

They are designed for use in laboratories. The motion of the agitation table or of the containers placed on it is circular.

For usage in accordance with requirements, the machine must be standing on a stable, even surface that is slipfree. In addition, care must be taken that objects in the vicinity are a sufficient distance away from the shaking unit and that they will not come too close to the shaking unit while it is in operation.

Unpacking

Please unpack the equipment carefully and check for any damages. It is important that any damages which may have arisen during transport are ascertained when unpacking. If applicable a fault report must be set immediately (post, rail or forwarder).

The delivery scope covers:

MS 1: A shaker **MS 1**, a standard attachment MS 1.1, a one-hand attachment MS 1.2, a test tube insert MS 1.32, a microtiter plate attachment MS 1.4 (also see p. 6, accessories), a mains cable and an operating instructions.

MS 2: A shaker **MS 2**, a standard attachment MS 1.1, a one-hand attachment MS 1.2 (also see p. 6, accessories), a mains cable and an operating instructions.

General information

With the purchase of this device, you have acquired a high-quality product.

The speed-controlled motor allows for infinity speed adjustment in the range from 0/200 to 2500 rpm. Electronic motor control holds the set speed constant even if the weight of the material on the surface increases.

The shaker **MS 1** and **MS 2** are suitable for continuous operation.

The casing may heat up when the appliance is in use for a prolonged time.

The various attachments can be changed easily by pulling them up off the receptacle on the appliance. Ensure that the attachment is positioned correctly when fitting. see spare parts diagram **MS 1** and **MS 2**, p. 22.

Commissioning

Check whether the voltage specified on the type plate matches the mains voltage available. The power socket used must be earthed (protective earth conductor contact). If these conditions are met, the device is ready to operate after plugging in the mains plug. If these procedures are not followed, safe operation cannot be guaranteed and/or the equipment may be damaged. Observe the ambient conditions (temperature, humidity, etc.) listed under "Technical Data".

Switching on and adjustment of speed

The rocker switch (A) for switching on the shaker is on the left of the appliance.

The appliance is switched off in the centre position "off".

The knob (B) for setting the speed is on the front of the appliance in the centre. The speed can be continuously adjusted within a range of 200 to 2500 rpm.

The shaker cannot be set below 200 rpm as it will be at a standstill. When cold, the lowest possible speed may be slightly higher.


Fig. 2: Operating elements

Adjustment of the operation mode

Touch operation mode

To set the shaker to touch mode, set the rocker switch (A) to "Touch". You must also release the rocker catch (D) on shaker **MS 1**: select the "Touch" position. If the speed is set to 200 rpm or higher, you can start the shaking process by pressing on the attachment.

Continuous operation mode

To set the shaker to continuous mode, set the rocker switch (A) to "On". If the speed is set to 200 rpm or higher, the appliance will begin to shake immediately. When the knob (B) is at the left stop, the shaker will initially stand still, allowing you to set the required speed. You can also set the rocker catch (D) for shaker **MS 1**: select the "platform" position.

Attachments

Attachments recommended by **IKA** must be connected to the agitation table to ensure safe operation. In addition, the usable weight and speed range must not be exceeded.

MS 1.1 Standard attachment

This attachment is used to attach and shake tubes and other small vessels of diameter up to 30 mm. A trough prevents the vessels from slipping. For work in touch mode.

MS 1.2 Test tube insert

This attachment comprises the MS 1.3 test tube attachment and the MS 1.21 one-hand insert.

MS 1.21 One-hand insert

Foam inlay for inserting into MS 1.3 test tube attachment. Flat containers such as Petri dishes or culture bottles can be shaken at a slow shaking rate in continuous mode. Any larger or taller vessels you may wish to shake, for example Erlenmeyer flasks, must be held in place. This insert is also suitable for working in "touch mode".

MS 1.3 Test tube attachment


This attachment can hold various inserts. It is possible to work in both continuous mode and touch mode. The inserts are described below:

MS 1.31 Test tube insert

Foam inlay for inserting into MS 1.3 test tube attachment. This insert can hold up to 14 test tubes with a diameter of 8 mm. For work in continuous mode.

MS 1.32 Test tube insert

Foam inlay for inserting into MS 1.3 test tube attachment. This insert can hold up to 6 test tubes with a diameter of 10 mm. For work in continuous mode.

MS 1.33 Test tube insert

Foam inlay for inserting into MS 1.3 test tube attachment. This insert can hold up to 4 test tubes with a diameter of 15 mm. For work in continuous mode.

MS 1.34 Test tube insert

Foam inlay for inserting into MS 1.3 test tube attachment. You can make as many holes as you like in this insert. For work in continuous mode.

MS 1.4 Microtiter plate attachment

This attachment is used to hold a microtiter plate.

Supported load and speed range

Shakers **MS 1** and **MS 2** are designed for a max. load of 0.5 kg (incl. attachment). The recommended maximum speeds for the respective attachments are specified below. These values are, however, only guidelines. The user must always ensure a safe operation. The user must determine the limit speed for the existing conditions by gradually increasing the speed. If these limits are exceeded, the media to be shaken may escape or be projected out of the vessels. With microtiter applications, excessive speed can cause samples from different chambers to get mixed together. Furthermore, excessive speed can cause vessels to slide off or be projected out of the inserts.

	Max. speed [rpm]	Touch- operation mode	Continuous operation mode
MS 1.1 Standard attachment	2500	yes	no
MS 1.2 One-hand attachment	1400	yes	yes
MS 1.3 Test tube attachment with one-hand insert	1400	no	yes
MS 1.4 Microtiter plate attachment	800	no	yes

If the shaker is screwed down (i.e. built into a system), please note that dynamic forces can be generated during operation as a result of unfavourable loading and a poor centre of gravity position. This could cause the surrounding area to vibrate or induce uncontrollable vibrations in the table. If you notice that the appliance is not running smoothly, reduce the speed until it runs smoothly again.

Maintenance and cleaning

The shaker **MS 1** and **MS 2** are maintenance-free. It is subject only to the natural wear and tear of components and their statistical failure rate.

When ordering spare parts, please give the manufacturing number shown on the type plate, the machine type and the name of the spare part.

Please send in equipment for repair only after it has been cleaned and is free from any materials which may constitute a health hazard. Please tell us if you have used any materials which are dangerous or which may constitute a health hazard.

Use only cleansing agents which have been approved by **IKA** to clean **IKA** devices.

To remove use:

Dyes	isopropyl alcohol
Construction materials	water containing tenside/ isopropyl alcohol
Cosmetics	water containing tenside/ isopropyl alcohol
Foodstuffs	water containing tenside
Fuels	water containing tenside

For materials which are not listed, please request information from **IKA**. Wear the proper protective gloves during cleaning of the devices. Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.

Before using another than the recommended method for cleaning or decontamination, the user must ascertain with the manufacturer that this method does not destroy the instrument.

Accessories

Designation	Description	Id-No.	MS 1	MS 2
MS 1.1 Standard attachment	For small tubes and vessels upto diameter 30 mm	L001260	X	X
MS 1.2 One-hand attachment	Consisting of MS 1.3 and MS 1.21	L001900	X	X
MS 1.3 Test tube attachment	Attachment for inserts MS 1.21, MS 1.31, MS 1.32, MS 1.33, MS 1.34	L001490	X	X
MS 1.4 Microtiter plate attachment	For microtiter plate	L001280	X	-
MS 1.21 One-hand insert	For vessels with flat base	L001540	X	X
MS 1.31 Test tube insert	For 14 vessels with diameter 8 mm	L001840	X	X
MS 1.32 Test tube insert	For 6 vessels with diameter 10 mm	L001850	X	X
MS 1.33 Test tube insert	For 4 vessels with diameter 15 mm	L001860	X	X
MS 1.34 Test tube insert	For any number of bore holes	L001830	X	X

Warranty

You have purchased an original **IKA** laboratory machine which meets the highest engineering and quality standards.

In accordance with **IKA** warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover wearing parts, nor does it apply to faults resulting from improper use or insufficient care and maintenance contrary to the instructions in this operating manual.

If you require servicing, return the appliance in its original packaging. Storage packaging is not sufficient. Please also use suitable transport packaging.

Technical data

Design voltage:	VAC	230±10%
	oder VAC	115±10%
Design frequency:	Hz	50/ 60
Input power:	W	31
Output power	W	10
Shaking rate:	rpm	0 - 2500 (infinitely adjustable from 200)
Speed adjustment:		analog (infinitely)
Power supply acc. to IEC 60320:		cold device plug with fuses
Fuses:	mA	2 x 400T (230V)
		Ident-Nr. L001060
(in cold device plug)		2 x 630T (115V)
		Ident-Nr. 2008600
Mains cable:		with ground wire and in compliance with country specification
Shaking rate display:		Scale 0 - 2500
Drive:		Speed controlled DC-motor, with ball bearings
Protection at overloaded:		Temperature sensor in motor winding
Agitation stroke:	mm	4.5
Shaking motion:		circular
Perm. duration of operation:	%	100
Perm. ambient temperature:	°C	+5 bis +40
Perm. relative humidity:	%	80
Protection type acc. to DIN EN 60 529:		IP 21
Protection class:		I
Overvoltage category:		II
Contamination level:		2
Operation at a terrestrial attitude:	m	max. 2000
Max. supported load incl. attachment:	kg	0,5
Dimernsions: (W x D x H)	mm	115 x 225 x 60
Weight without supported load:	kg	3.5

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List of spare parts

Item	Designation
10	Potentiometer
20	Transformer
30	Rocker switch
40	Fuse
60	PCB
90	Front panel
120	Collar
130	Standard attachment MS 1.1
140	O- ring
150	Rubber plate T 55x25
155	Rubber plate P 98x12
170	Rubber base
200	Leaf spring
210	Motor compl.
220	Rocker assy.
230	Panel (only MS 1)
235	Cap
250	Tension roller
260	Operating knob
270	Receptacle
290	Compression spring
300	Snab knob (only MS 1)
310	Leaf spring (only MS 1)
1902	1.32 Test tube insert (only MS 1)
1903	1.3 Test tube attachment
1904	1.4 Microtiter plate attachment (only MS 1)
1905	1.21 One-hand insert

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