DILUMAT 3mk2

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DILUMAT 3mk2

EC COMPLIANCE CERTIFICATE

AES Laboratoire declares that
the DILUMAT 3 mk2 (ref: AESAP1055)
complies with the following European directives:

- **89/336/CEE Directive** (Electromagnetic Compatibility)
Updated by 92/31/CEE 93/68/CEE directives
3 MK2
Reference Standard:
 - Emission
Generic Standard: NF EN 50 081 – 1 (June 1992)
Basic Standard: NF EN 55011 (conducted and radiated)
 - Immunity
Generic Standard: NF EN 50 082-1 (June 1992)
Basic Standard: NF EN 61 000-4-2 (Electromagnetic Burst)
NF EN 61 000-4-3 (Immunity to radiant fields)
NF EN 61 000-4-4 (Fast transient bursts)
- **73/23/CEE Directive** (Electrical Equipment Safety)
Updated by the 93/68 CEE directive
Reference Standard:
 - NF EN 61010-1 (08/93 + A2 11/95)
- **98/37/CEE Directive** (Machines Safety)

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You have just bought a Dilumat 3mk2.

This diluter was designed to combine high performances and simplicity in use.

The Dilumat 3mk2 was made to fulfil any of your expectations. Complete tests have been carried out to ensure the quality and reliability of the unit.

The functioning and technical features comply with the EC regulations (a certificate is enclosed in this manual).

We made great efforts to design a machine that is convenient to use.

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The quality and reproductibility of microbiological analysis results depend on the sample treatment (sampling, dilution, and homogenisation).


The Dilumat 3mk2 aims at simplifying these operations.

We have made great efforts to produce a machine that is both outstanding in performance and user friendly:

- **The LCD screen** allows a simple and intuitive approach.
- It is easy to scroll up and down the main pages.
- The default **peristaltic pump** ensures a fast and accurate dispensing
- **A second pump** unit can easily be added as an option.
- **A printer** (standard or compact) can be connected to ensure traceability.(option)
- **A bar code reader** can be connected so that the operator's task becomes even easier (Labmonitor option).
- Dispensing can be activated via **a foot switch** (option).
- **The load cell** ensures complete accuracy for samples dilution.

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4.1 Installation of the Dilumat 3mk2

The Dilumat 3mk2 should be placed on a flat and horizontal surface.

The room power supply should be 230 VAC or 115 VAC with fuses of F1A L 250V. (limits: 110/120V +/-10% 0.5A, or 220-240V +/-10% 0.3A).

Caution: the bag holder is sensitive to ultraviolet.

- ✓ Screw the arm with the screws that are provided.
- ✓ Place the bag holder stand.

The bag holder may be placed either parallel or perpendicular to the operator, as shown on the photos below:

CAUTION :

To move the bag holder from one position to the other one, it is essential to lift the bag holder up before turning it into the other position.

- ✓ Place some sticky tacs on the bag holder.

Before plugging the unit in, check that the voltage switch position (at the rear) fits the mains supply voltage of the room.

- ✓ Plug the Dilumat 3mk2 in, using the provided cable.
- ✓ Switch the unit on (main switch at the back of the unit). The screen lights up and the machine initialises.

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4.2 Installing the tube on the pump

To ensure good operation and accurate dispensing, **it is essential to place correctly the sterile dispensing tubing set.**

Follow the 3 steps procedure that is schematically represented below:

1. Turn the tubing locking clip anticlockwise to get the plastic strip out.

2. Slide the tubing into the holding strip.

3. Place the tubing around the pump rotor.

Caution : Fitting should be perfect.

The triangular shapes of the strip should fit into the pump head.

It is important to install tubing according to the liquid circulation (as shown with arrows on the photo beside).

Pump rotation:

Both pumps, either the default one or the optional one, turn anticlockwise.

4. Finally, turn the locking clip clockwise.

5. Once fitted, tubing should be neither bent nor pinched. Install the dispensing nozzle onto the arm.

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↵

5. OPERATION

When switching the diluter on, an initialisation screen appears :

**AES LABORATOIRE
DILUMAT 3 V 5.20**

Then the main screen is displayed :

**Imp. ON Fac. --
Dilution <Y> / <N> ↘**

This screen informs the operator :

- About the printer (if it is disabled or not : ON / OFF).
- About the dilution factor that has been used for the previous operation.
- About the pump that is currently selected (↘ front or ↗ back), if the diluter is equipped with the second pump option.

5.1 Getting used to the keyboard and moving from one menu to the other.

Getting used with the keyboard :

During a dilution, numerical keys (0,1, 2, 3) are used as shortcut keys :

- **0** to go out of a menu
- **1** to allow for the tare
- **2** to change pump
- **3** to modify the dilution factor

Press « Yes » to enter a menu or carry on with an operation.

Press « No » to go out of a menu or switch to another one.


When dialing, you may use this key to move the cursor backward and modify figures.

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
Moving from one menu to another:

The various menus are easy of access.

The flow chart below shows how you may go from one menu to another.

By pressing "NO", you may go from menu to the next one.

Please note:


In the "**Parameters**" menu you may enter dilution or weighing settings.

- dilution factor, sample minimum and maximum weight (in dilution mode)
- amount to be dispensed (in dispensing mode)

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
Dilution settings are determined within the « **Parameters** » menu.

Please follow this procedure:

The sample is already in the bag when the latter is placed on the stand.
The tare is not equal to 0 but it is "0 - the weight the registered bag".
In this case the weight of the bag is identified prior any operation.
This is valid if the bags that are used are identical.

BAG WEIGHT
-- , --g <ENTER>

Enter the actual mean weight of the empty bag.
This weight will then be taken into account.
Ex: P= 10g, the screen displays P= -10g once you have pressed the 1 key to tare.

The Dilumat is in the normal functioning mode.

For each dilution an empty bag is hung on the stand, then a sample is placed into the bag.

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5.2 Starting a dilution

- ⇒ Once the machine is switched on, wait for the load cell to settle down (between 15 and 30 minutes). Place the sterile tubing onto the pump rotor (cf. page 9).
- ⇒ Place the dispensing nozzle onto the diluter arm.
- ⇒ Then insert the priming nozzle into the sterile diluent container, or the 9l/19l/ 60l Masterclave.
- ⇒ Place the bag on the bag holder stand. It is essential that the upper part of the bag is level with the stand upper sides. The dispensing nozzle should not touch the bag, as this could affect dilution or dispensing.
- ⇒ Pinch both sides of the stand against the bag, and let it go. The bag opens up and is maintained by the sticky tacs.
- ⇒ Select the Dilution menu (scrolling down the menus by answering « No »).
- ⇒ Then check that the dilution settings are suitable.

In the example below, the screen indicates that:

- A printer is connected and set as switched on (option).
- The dilution factor is 10.
- The diluter is equipped with a second optional pump, and the front pump is selected (arrow pointing downwards).

If the displayed settings are not appropriate to the dilution that you wish to perform, please refer to the « Parameters settings » chapter, in order to modify them. You also may modify them during dilution, via the numerical shortcut keys.

Print. ON	Fac.10,0
Dilution	<Y> / <N> ⬇

- ⇒ Press the <Yes> key to enter the dilution menu.
- ⇒ If the diluter is connected to a printer, the operator code is asked (from 01 to 20). This is used to identify the person that makes the dilution. This code appears on the ticket that is printed once the dilution is performed. (if the diluter is equipped with the printer option, and if it is ON).
- ⇒ The following screen appears. On the upper line 2 messages are alternately displayed :

Place Bag	T <1>
Weight = ----,- g	⬇

Put arm down	P <2>
Weight = ----,- g	⬇

- ⇒ Place an empty bag on the stand. Press <1>, to allow for the tare. The message « **Tare in progress** » appears.

The displayed weight is 000,0 g.

- ⇒ Place enough sample for analysis into the bag, or replace the empty bag by another one containing a sample. In the second case, the « Automatic tare » mode should be set previously in the « Parameters » menu. (see previous page).

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- ⇒ When the sample is on the diluter, its actual weight is displayed on the second line of the screen :

Place Bag T <1> Weight = 10,10 g	Put the arm down P <2> Weight = 10,10g
-------------------------------------	---

- ⇒ Lower the diluter arm to perform dilution.

When the dilution is over, the operator is warned by an audible signal. On screen this appears:

Dilution : 0101,0g
Weight = 0010,10g

On the upper line a message, the message alternately appears : « Lift the Arm ».

- ⇒ If a connected and switched on, a dilution report is automatically printed out, as you lift the diluter arm up. On this report you can read :

- the date
- The dilution factor that was used
- The date of the last calibration
- A dilution number, and the time at which the dilution was performed.
- The sample weight and the total weight after dilution.

(see example on page 23).

- ⇒ Remove the bag, close it with closing clips, and place it on a bag holder before blending.

- ⇒ You can read alternately on the screen:

Place Bag T <1> Weight = 10,10 g	Put the arm down P <2> Weight = 10,10g
-------------------------------------	---

When there is no bag on the stand, a negative weight is displayed (in the case that a tare of 0g was allowed, with an empty bag on the stand).

Then you may place a new bag on the stand and carry on with a new dilution without allowing the tare.

You simply have to insert the sample into the empty bag, to lower the diluter arm to start the dilution. Once the dilution is over, lift the arm up. You may perform as many dilutions as you wish, as long as dilution parameters do not change (dilution factor, pump, type of bag...).

If you wish, you may take time to make a tare before each dilution, by pressing <1>.

A QUICK REMINDER !

During a dilution, the 0,1,2,3 keys become shortcut keys :

- **0** is used to go out of a menu / stop temporarily the dilution (ex: to change bottle)
- **1** is used to make a tare
- **2** is used to change pump (if the optional pump is installed on the diluter)
- **3** is used to modify the dilution factor.

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5.3 Dispensing mode

Dilumat 3mk2 can also be used for dispensing doses.

- ⇒ Once the diluter has been switched on, wait for the load cell to settle down, and place the sterile tubing on the pump rotor. (cf. page 9).
- ⇒ Place the dispensing nozzle on the Dilumat 3 mk2 arm.
- ⇒ Then insert the priming tube into the sterile diluent container.
- ⇒ Place the bottle (or bowl) that will be filled in on the holding plate. It is advised to use a bowl (AESDI0043) stand with a spacer (AESDI0098 or AESDI0099) if required, so that bowl or bottles are properly held.
- ⇒ Select the « **Parameters** » menu to go and enter the weight (volume) to be dispensed. (cf. page 14).
- ⇒ Select the dispensing menu (scrolling down the menus by pressing « No »).

Print. ON
Dispense <Y> / <N> ➡

Note :

If the diluter is equipped with a second pump, the pump that was selected for this dispensing is indicated on screen. In the example below, the front pump (➡) was selected.

- ⇒ Press the "Yes" key to enter the dispensing menu.
- ⇒ If the diluter is connected to a printer, the operator code is asked (from 01 to 20). This is used to identify the person that makes the dispensing. This code appears on the ticket that is printed once the dilution is performed. (if the diluter is equipped with the printer option, and if it is ON).
- ⇒ This screen appears. Two messages are alternately displayed on the top line:

Place Bag T <1>
Dose = 0,0 g ➡

Lower arm P <2>
Dose = 0,0 g ➡

P <2> : If the diluter is equipped with a second pump, you can select the pump that will be used for this dispensing.

- ⇒ Place the container to be filled in onto the plate. Its weight is displayed on the second line.
- ⇒ Press the <1> key, to make a tare. The message « **Tare in progress** » appears. The displayed weight is 0,0 gr.
- ⇒ Then you may start dispensing by lowering the diluter arm. An audible signal informs you that dispensing is over.

Lift arm
Dose = 16,0 g ➡

- ⇒ Lift the diluter arm up and remove the container.
- ⇒ The next dispensing can be performed, according to the same procedure: container placing, tare, 15 arm lowering.

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6. DILUMAT 3MK2 SETTINGS

When the diluter is equipped with options (printer, 2nd pump, foot switch), some additional menus appear to give you the possibility to change these accessories settings.

Some parameters such as operation language, date or calibration may be set on the Dilumat 3 mk2.

6.1 Printer Menu

If the Dilumat 3 mk2 is equipped with a printer (AESDI0042A Epson printer for A4 printouts, or AESDI0042A Star Printer for small printouts), it is in this menu that you can activate it.

Within this menu you may determine the type of printer that you are using, and the number of tickets you wish to print. You may also take samples numbers into account.

- ⇒ Go to the « Printer » menu (cf. chapter 5.1. «Moving from one menu to another ».)

PRINTER MENU
<Y>/<N>

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PRINTER USED
<Y>/<N>

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PRINTER SELECTION
CITI. <0> STAR <1>

- ⇒ Dial the number 0 or 1 according to the type of printer that is connected.
<0> : CITIZEN 120 D+ printer
or EPSON LX300 printer.

<1> : STAR small printer.

TICKET PRINTING
<Y>/<N>

YES

NO

TICKETS NUMBER
-- <ENTER>

Dial the number of tickets (1 to 9)
That you need, then confirm with
ENTER.

Reports are
printed on
paper rolls.

SAMPLE REFERENCE
<Y>/<N>

If you answer YES, for each future dilution, you will have to enter a numeric code to identify samples.

If you answer NO, the actual sample reference is not taken into account.

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6.2 Pump menu

In this menu you may assign a volume of diluant that is available for each pump.

This is quite helpful when the Dilumat is connected to an autopreparator (such as S8000, Masterclave 19 or 60), as the remaining volume cannot be checked visually.

- ⇒ Display the "Pump" menu (cf. Chapter 5.1 "moving from one menu to another".)

Diluant --- ml	<Y>/<N> ↘
Pompes ?	

- ⇒ The volume that remains in the Autopreparator or container is displayed.
The arrow indicates which pump is selected.
Press YES to modify something.

Choix pompe	↗ <0> ↘ <1>
-------------	------------------

- ⇒ Dial the pump number that you want to set.
Ex: <0>.

VOLUME DILUANT	<ENTER> ↗
---- ml	

- ⇒ You may enter the diluent volume that is available in the autopreparator or container. Press ENTER to confirm.

Diluant --- ml	POMPES ?
<Y> / <N> ↗	

YES

NO

CHOIX POMPE	↗ <0> ↘ <1>
-------------	------------------

Impr. ON	Fac. --
Dilution	<Y>/<N> ↘

If you wish to enter the second pump settings

- ⇒ Select the second pump, to enter its settings.

Access to the main screen.

NOTE:

If a volume of diluant is taken into account, the diluter takes out each dispensed volume, and calculates automatically the remaining volume.

You are informed when the container is empty. You have the possibility to stop the operation in order to replace the bottle of diluent:

DILUENT EMPTY!	<Y>/<N> ↘
Dilution	

If you answer "YES":

The operation carries on.

Volume management is stopped.

If you answer "No":

The operation is stopped.

The pump menu appears to give you the possibility to enter the new volume of diluent.

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PEDALE ?
<Y>/<N>

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6.3 Foot switch menu

If the Dilumat 3 mk2 is equipped with a foot switch (option : AESPP0019), it is in this menu that you can activate it.

With a foot switch it is possible to leave the arm in the low position and start a dilution by pressing the foot switch.

This might be useful when the Dilumat is operating under a class II Safety Cabinet, whose protection window does not make it easy to lower the diluter arm.

- ⇒ Go to the **“Foot switch “ menu** (cf. chapter 5.1 going from one menu to another one)

FOOT SWITCH ?
<Y>/<N>

- ⇒ If you want to start a dilution or a dispensing via the foot switch, answer “Yes” to enter the menu.

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6.4 Language setting

To change the operation language, you have to :

- ⇒ Switch the machine off (button at the back of the diluter).
- ⇒ Then switch the machine back on (ON/OFF button), pressing simultaneously the “Enter” key ↵.
- ⇒ This screen appears:

CALI. <0> LANG. <1>
DATE <2> TEST <3>

- ⇒ Press the <1> key to have access to the **“Language”** menu.

ENGLISH ?
<Y>/<N> ↵

- ⇒ Use the “No” and “Yes”, to select you operation language.

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6.5 Date and time adjustment

In order to change the displayed date on the diluter, you have to:

- ⇒ Switch the diluter off (switch at the back of the machine).
- ⇒ Then switch it back on, pressing simultaneously the "Enter" key ↵.
- ⇒ This screen appears:

CALI. <0>	LANG. <1>
DATE <2>	TEST<3>

- ⇒ Press <2> to have access to the "**DATE / TIME**" menu.

DAY/ MONTH / YEAR
27/09/02

- ⇒ Use the keyboard to modify the date, if required. Confirm your entry by pressing "Enter".

- ⇒ The date is displayed.

HOUR : MINUTE
10 : 29 : 33

- ⇒ You may change it if you need to. Press "Enter" to confirm. Switch the unit off and on to go back to the operation mode.

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6.6 Calibration

To access the « **calibration** » menu, you have to :

- ⇒ Switch the diluter off (switch at the back of the machine).
- ⇒ Then switch it back on, pressing simultaneously the "Enter" key ↵.

⇒ This screen appears:

CALI. <0>	LANG. <1>
DATE <2>	TEST <3>

⇒ Press <0> to go to the « **CALIBRATION** » menu.

EMPTY STAND <Y> / <N>

⇒ Press "Yes" to confirm that the bag holding stand is actually empty.
The message "Tare in progress" appears, followed by this screen:

TARE WEIGHT ---,-- g	<ENTER>
-------------------------	---------

⇒ Place the calibration weight on the bag holding stand, and dial its certified value.
⇒ Confirm your entry by pressing "Enter" ↵.

⇒ The message "Tare in progress" is displayed again.
Then the diluter requires the calibration weight a second time, for checking.

Ex: 500g

PLACE 500, 00g 500 , 06 OK? <Y>/<N>
--

⇒ If displayed values are sufficiently close, press "Enter" to confirm.

RESET A 0 g 1,0g	<ENTER>
---------------------	---------

Caution:

For usual operations (weights from 10 to 25gr), we advise the default value that is 1gr.
For very light weighings (0,5 to 1gr), this value can be dropped down to 0,3 (especially when using powder sample).

⇒ Once calibration is over, you are offered to print a calibration report (if your diluter is equipped with the optional printer).
An example of calibration report is given in page 23.
Switch the unit off and then back on to go back to the "Operation" mode.

Calibration should be carried out every year.

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7. TECHNICAL FEATURES

Dimensions without package: 50 x42 x 39 cm
Dimensions with package: 35 x 50 x 50 cm

Weight without package: approximately 11 kg
Weight with package: approximately 13 kg

Environment conditions to ensure a good working order:

- Needs to be used indoors.
- Up to 2000 meters height.
- Temperature: from 5 to 40°C
- 80% Maximum humidity for temperatures up to 31°C, with linear decrease up to 50% of relative humidity at 40°C.
- Power supply: 110/220V or 220/240V +/-10% A, 50-60Hz; 0.5 / 0.3A .
- Class 1 equipment that can be exclusively connected to TT or NT, and connected to earth.
- Category II installation according to IEC664 directive.
- Pollution: degree 2

Dilution factor: from 2,0 to 150,0 (1/2 to 1/150).

Minimal sample weight: from 1 to 25 grams.

Usual operational range for dilution of solid samples: from 100 to 400g (diluent + sample), provided that samples weight does not exceed 40 grams.

Dispensing speed: 800g/min

Dilution accuracy:

Dilution: from 0 to 40g: +/- 0.3g
 from 40 to 100g: +/- 0.4g
 from 100g to 1 Kg: +/- 0.5g

Dispensing: from 0 to 40g: +/- 0.3g
 from 40 to 100g: +/- 0.4g
 from 100g to 1 Kg: +/- 0.5g

Maximal weight allowed on the holding plate: 1500gr.

Minimal weight:

- in automatic mode: 5 gr
- in manual mode: 1 gr (with automatic tare before each dilution).

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8. OPTIONS AND ACCESSORIES

The Dilumat 3 mk2 may be equipped with the following options:

SECOND PUMP OPTION


- ❑ 2nd peristaltic pump option.
Ref: AESDI0040
Factory assembled.

or

- ❑ 2nd peristaltic pump kit
Ref: AESDI0040A
Delivered and installed later on.

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
With this printer, you can print traceability reports about dispensing, dilutions, calibrations or any tests that are performed.

Any information mentioned on these reports depend on the operations that are performed.

```

FICHE DE CALIBRATION
DILUMAT3 V5.26 NUMERO DE SERIE: 23331091

DATE: 19/12/02 CODE OPERATEUR: 01
PROCEDURE DE TARAGE 0g-1000,00g: OK
LECTURE POIDS VERIFIE: 999,96g
CALIBRAGE VALIDE
  
```

```

# CALIBRATION REPORT
DILUMAT3 V3.2 SERIAL NUMBER: 10400061

DATE: 20/03/01 OPERATOR IDENT.: 01
TARE PROCEDURE 0g-1000,00g: OK
DISPLAY CHECKING WEIGHT: 1000,03g
CALIBRATION FINISHED
  
```

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```

SAMPLE REFERENCE: 000002
OPERATOR IDENT.: 01
DATE: 20/03/01 10:50
LAST CALIBRATION: 20/03/01
WEIGHT SAMPLE: 19,94g
DILUTION FACTOR: 10
WEIGHT DILUTION: 199,5g
  
```

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□ **PC connection**
(ref : AESDI0312)

The Dilumat 3mk2 may be connected to a pc in order to record traceability reports with Hyperterminal .

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□ **Installation in Labmaster network mode**
(ref: AESDL0107)

Connected to a Labmonitor and a bar code reader are installed, the Dilumat 3mk2 can operate in a LABMASTER network mode.

The Dilumat 3mk2 communicates with the Labmaster application via the Lonworks network.

- When the operator scans the identification bar code on the Labmonitor, the program's parameters are transmitted from the Labmonitor to the Dilumat 3mk2 via the Labmaster central PC.
- Then the operation results are sent to the database that manages them.
- Monitoring and automatic traceability are achieved.

To get more information, please contact AES Laboratoire.

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FIRST OF ALL, for your own safety it is essential , to switch the Dilumat 3mk2 off and unplug it . Disconnect any machine that may be connected to it.

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The Dilumat 3mk2 does not require a constraining daily cleaning.

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- ✓ In order to work in sterile conditions, the dispensing tubes needs to be autoclaved.
- ✓ The bag holder and bowl holder stands should be removed from the machine and cleaned once a day with a wet cloth and hot water. Dry them afterwards. Do not use abrasive cloth on the bag holder.

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CAUTION :

To remove the bag holder, it is essential to lift the bag holder up beforehand.
If not the bag holder base may break.

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- ✓ The cover needs also to be cleaned daily.
- ✓ A weekly cleaning with 70° alcohol is necessary.

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In case of spillage, make sure that no liquid runs to the upper part of the sensor. The bag holder and deflector special shapes avoid liquid infiltration to the unit sensor.

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If you need to remove the bag holder, you can protect the sensor by placing a plastic stopper (eg. GL45)- ref: BIB36382.

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- ✓ The precut tubing (AESDI0071 / AESDI0071A) that is installed on the pump, should be changed after 100 hours of use.

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ERROR MESSAGES DISPLAYED ON SCREEN :

« Reset »

Switch the machine off (switch at the back of the diluter). Then switch it back on.

« 0001 »

There might be a problem with the 24Volts power supply. Please refer to chapter « Main parts dismantling », in order to check the FU2 fuse.(page 28)

If the FU2 fuse is in good working order, check the transformer .

If no cause can be identified, a problem with the PCB may generate this message. In this case, you should contact the After Sales Service.

THE SCREEN IS COMPLETELY BLANK

Make sure that the power supply cable is correctly connected.

Switch the diluter off, then switch it back on.

If the problem occurs, please contact the After Sales Service.

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TROUBLESHOOTING

? The arm cannot be lifted up.

The arm shaft retaining spring may be broken.

Or

If the message "Lower arm" remains on screen, the optical sensor might be damaged.

In both cases, refer to the "Main parts dismantling" chapter, in order to get to these elements.

? Dilutions keep on being outside tolerances.


First of all, make sure that the top of the bag is not above the upper part of the holding stand. If so, the dispensing nozzle may touch the bag and this may false weighing.

If despite this precaution, dilutions are still outside tolerances, perform a calibration as explained in page 20.

? The pump motor is blocked.

Remove the pump holding strip, and try to turn manually the rotor. If rotation is difficult, please call the After Sales Service for detailed explanation.

Check the pump rotor. The rotor may be worn (A). If so, the dispensing tube may not fit properly and this stops the pump motor.
Please refer to the "Main parts dismantling" section.

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Figures display is unstable

or

Calibration is outside tolerances

or

the diluter does not take into account a weight placed on the plate. Displayed weight is 0.

Weight =00 ,0g
Dilution <Y> / <N>

Perform a calibration as explained on page 20.

If the problem recurs, check that the load cell is in good working order. Cf. « Main parts dismantling » .

If the problem does not come from the load cell, you may have to change the AD7710AN component. You should check with the After Sales Service before any intervention.

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Turn the locking clip anticlockwise to remove the pump holding strip.


Remove both holding screws, to release the pump. Then remove the circlip that holds the rotor.

Separate the axle from the rotor.

Rotor reference: AESDI0061

An exploded view of the pump is given at the end of this manual as an appendix.

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
Remove the bag holding stand.

CAUTION :

The bag holder may be placed either parallel or perpendicular to the operator. To move from one position to the other one, it is essential to lift the bag holder up before turning it into the other position. Otherwise the device base will be damaged.

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Remove the 3 holding screws at the back of the diluter.

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
Lift the back of the cover, and tilt it forwards.

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

You may lay the cover on its left side. The inside parts are now accessible (see next page).

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


Load cell dismantling and adjustment

- ➡ Remove the Dilumat 3 mk2 cover as shown in page 29.

- ➡ Tilt the Dilumat3 mk2 backwards to be able to unscrew the 2 screws under the machine.

- ➡ Remove the PCB armour-plate so that you can unplug the load cell connector.

- ➡ You may replace the load cell if required.

- ➡ **CAUTION:** in order to place properly the load cell, it is **IMPORTANT** that is centred on the hole made in the cover.


Load cell adjustment :

- ➡ Before placing the cover back on, switch the diluter on and perform a calibration (see procedure in page 20).
- ➡ Place a 2000 g weight on the load cell. Then screw or unscrew until the diluter displays 1800 +/- 200g).
- ➡ Then you may reassemble the diluter.

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Electronic board access

Once the dilumat 3mk2 cover has been taken off, you can have access to the PCB by removing the armour-plate.

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
AD7710AN component

To get to this component, unscrew the two screws that hold the printer board. Then remove the printer board.

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This component should be removed only with an appropriate tool.

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11. WARRANTY CONDITIONS

AES Laboratoire certifies that the delivered machine corresponds to the features described in this manual. Our equipment is guaranteed against manufacturing or components problems for twelve months from the delivery date.

AES Laboratoire will repair broken parts free of charge in its premises during the warranty period. Any instrument with a working problem should be, with the After Sales Department's approval, sent back carriage paid by AES agent or customer when necessary.

The warranty will not be valid for damages due to a transport problem, an accident, a wrong use or internal or internal interference caused by a faulty electricity network or non-authorized users.

12. LISTE DE PIECES DETACHEES**12. SPARE PARTS LIST**

REFERENCE	Français	English
AESDI0002	SUPPORT SAC EN ALTUGLASS	BAG HOLDER (PERSPEX)
AESDI0008	DEFLECTEUR DILUMAT 3	DILUMAT 3 DEFLECTOR
AESDI0011	EMBOUT DE DISTRIBUTION	FILLING NOZZLE
AESDI0011A	EMBOUT DISTRIBUTION COUDE 18°	FILLING NOZZLE (F. BLENDER CUP
AESDI0013	SILENTBLOC MOTEUR	MOTOR RUBBER STAND
AESDI0016	POMPE COMPLETE SR25	PUMP (COMPLETE)
AESDI0017	RACCORDS INOX	CONNECTORS (STA.STEEL) 2U
AESDI0021	PONT DE JAUGE (CAPTEUR BALANCE	LOAD CELL
AESDI0022	AFFICHEUR	DISPLAY SCREEN
AESDI0023	CARTE ELECTRONIQUE (DILUMAT 3)	ELECTRONIC P.C BOARD (DIL 3)
AESDI0024A	CARTE IMPRIMANTE SERIE	INTERFACE BOARD SERIAL PRINTER
AESDI0025	MOTEUR PAS A PAS	STEPPING MOTOR
AESDI0026A	CORDON AFFICHEUR DILUMAT3	DILUMAT 3 DISPLAY CABLE
AESDI0027	FOURCHE OPTIQUE (POSITION BRAS	OPTICAL SENSOR (ARM POSITION)
AESDI0040	OPTION DEUXIEME POMPE DILUMAT3	SECOND PUMP KIT (COMPLETE)
AESDI0040A	KIT DEUXIEME POMPE DILUMAT3	SECOND PUMP KIT (COMPLETE)
AESDI0041	JEU DE TUYAU COMPLET	TUBING WITH NOZZLES
AESDI0041A	JEU TUYAU SANS EMBOUT PRELEVE.	SILICONE TUBE W/O NOZZLES DIL3
AESDI0042	IMPRIMANTE EPSONLX300 DILUMAT3	PRINTER (SERIAL MODELL)
AESDI0042A	IMPRIMANTE DILUMAT3 PETIT MOD	DILUMAT3 PRINTER SMALL SIZE
AESDI0042B	PAPIER LIST. PETITE TAILLE 5R	LISTING PAPER 114mm SMALL SIZE
AESDI0042C	ETIQUETTE PETITE TAILLE 5RLX	5x880 STICKER PAPER SMALL SIZE
AESDI0042F	RUBAN ENCREUR IMPRIMANTE STAR	STAR PRINTER INKING RIBBON
AESDI0043	PLATEAU SUPPORT BOL	BOWL STAND (COMPLETE)
AESDI0044	4 PIEDS CAOUTCHOUC	FEET (RUBBER) 4 U
AESDI0045	PLATEAU SUPPORT SAC COMPLET	BAG HOLDER (COMP.)W.FITTINGS
AESDI0045A	SUPPORT SAC COMPLET PETIT SAC	BAG HOLDER (COMP.)F.SMALL BAG
AESDI0046	TUBE DE PRELEVEMENT SEUL	PRIMING TUBE
AESDI0047	TUBE FIN ASPIRATION	AIR VENTING STA.STEEL TUBE
AESDI0048	FILTRE PURADISC 0.2 U (50 UNIT	FILTER (0.2um) 50 U
AESDI0049A	RUBAN ENCREUR IMPRIMANTE LX300	LX300 PRINTER INKING RIBBON
AESDI0052C	CABLE LIAISON IMPRIMANTE	CONNECTION CABLE FOR PRINTER
AESDI0053	TRANSFORMATEUR TORIQUE DILUMAT	TRANSFORMER
AESDI0054	FILTRE SECTEUR AVEC BOUTON M/A	ON/OFF SWITCH WITH FILTER
AESDI0055	BOUCHON SILICONE	SILICONE CAP
AESDI0057	RALLONGE CLAVIER DILUMAT	EXTENSION FOR DILUMAT KEYBOARD
AESDI0058	TUYAU SILICONE 3/16(5-8MM) 2m	TUBING (SILIC)3/16(5-8MM) 2M
AESDI0058A	TUYAU SILICONE 5-8MM (25 M)	TUBING (SILIC)3/16 5-8MM 25MET
AESDI0059	BANDEAU MAINTIEN TUYAU POUR PO	TUBING HOLDER (PLASTIC)F.PUMP
AESDI0060	CLIP DE VERROUILLAGE BANDEAU	CLOSING CLIPS F.PUMP TUB.HOLD
AESDI0061	ROTOR DE POMPE POUR DILUMAT 3	PUMP ROTOR
AESDI0062	EMBOUT DE PRELEVEMENT	PRIMING NOZZLE
AESDI0063	PATE ADHESIVE SACHET (78u)	STICKERS (YELLOW)T.HOLD BAGS

AESDI0064	LISTING ETIQUETTE 250FEUILLES	STICKERS PAPER FOR PRINTER
AESDI0070	COLLIER INOX	TIRE-RAPS (STA. STEEL) 2 U
AESDI0071	TUYAUX PREDECOUPE POMP DIL3	TUBING (PRECUT)F.PUMP 5u
AESDI0071A	TUYAU ROTOR POMPE DILUM. (1M)	TUBING ROTOR PUMP FOR DILUMAT3
AESPP0019	INTERRUPTEUR A PIED	FOOT SWITCH
AESDI0085	RESSORT SUPPORT EMBOUT DIST	SPRING CLIP FOR DISPENS NOZZLE
AESDI0086	AXE BRAS DILUMAT3 MK2	DILUMAT 3MK2 ARM SHAFT
AESDI0088	DRAPEAU AXE BRAS DILUMAT3 MK2	FLAG DILUMAT 3MK2 ARM SHAFT
AESDI0090	BOITIER INFERIEUR DILUMAT3 MK2	LOWER CASING DILUMAT 3MK2
AESDI0091	PORTE EMBOUT DE DISTRIBUTION	DISPENSING NOZZLE DILUMAT3
AESDI0092	RESSORT DE RAPPEL AXE DE BRAS	ARM SHAFT RETAINING SPRING
AESDI0098	CALE HT 70mm Pr PLATEAU BOL	SPACER FOR BOWL STAND
AESDI0099	CALE HT 125mm Pr PLATEAU BOL	SPACER FOR BOWL STAND
AESDI0103	CLIPS SUPPORTS TUYAUX DILUMAT3	TUBING SUPPORTING CLIPS DIL3
AESDI0224	CLE SIX PANS DE 2,5	HEXAGONAL KEY (2.5)
AESDI0226	ROULEMENT INT POMPE DILUMAT 3	INTERNAL BALL BEARING DIL3PUMP
AESDI0227	ROULEMENT EXT POMPE DILUMAT 3	EXT BALLBEARING DIL 3 PUMP
AESDI0229	ENJOLIVEUR BRAS DILUMAT3 MK2	DILUMAT 3MK2 ARM CAP
AESDI0300	CHAPPE ENJOLIVEUR DILUMAT3 MK2	DILUMAT 3MK2 CAP
AESDI0301	TOLE SUPPORT FOURCHE OPTIQUE	SHEET METAL HOLDG OPTIC SENSOR
AESDI0303	CLAVIER DILUMAT 3 MK2	DILUMAT 3 MK2 KEYBOARD
AESDI0305	CAPOT SUP DILUMAT3 MK2 COMPLET	UPPER COVER DILUMAT 3 MK2
AESDI0306	CAPOT SUPERIEUR AVEC CLAVIER	UPPER COVER WITH KEYBOARD
AESDI0307	BRAS COMPLET DILUMAT 3 MK2	DILUMAT 3 MK2 COMPLETE ARM
AESDI0308	BOUCHON SILICONE PERCE D36/44	SILICONE CAP W/ HOLE D36/44
AESDI0312	KIT COMMUNICATION DILUMAT3/PC	DIL 3 COMMUNICATION KIT W/ PC

DILUMAT 3MK2

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Which tubing set ?

AESAP1055

AESDI0041

AESDI0308

AEB611499L
(Tryptone Sel)
AEB110308M
(Eau Peptonée)
AEB610419L
(Fraser Demi)

AESAP1055

AESDI0041

SAM2521000
SAM2522000
SAM2525000

AESAP1055

AESDI0041A

AESAT00117AA

AESAP1047

AESAP1055

AESDI0041A

AESAP0422

AESAP1060

AESDI0041

AESDI0041A

+

AESDI0062

AESDI0041A

AESDI0071

+

AESDI0017

+

AESDI0011

+



AESDI0058 (2m)
AESDI0058A

ANNEXES

APPENDICES


POMPE PERISTALTIQUE PERISTALTIC PUMP


IT.5.5.2 DU 05/06/2000

COPIE ORIGINALE

FICHE D'IMPLANTATION		AES	
 <p>The diagram shows the layout of the RS232 'DILUMAT 3' card. At the top is a large rectangle labeled K2. Below it are six vertical rectangles labeled F6, F3, F5, F4, F2, and F1 from left to right. To the left of these are three circles labeled C3, C4, and C1, each with a '+' sign above it. Below C1 is a circle labeled R1 with a vertical line through its center. To the right of the F rectangles is a horizontal rectangle. Below that is a circle labeled C5 with a '+' sign, and below it is a circle labeled C2 with a '+' sign. At the bottom is a horizontal rectangle labeled K1.</p>			
CREÉ LE : 20/12/00 PAR : RENAULT E		CARTE RS232 "DILUMAT 3"	
MODIFIÉ LE : PAR :		VERIFIÉ LE : PAR :	
		INDICE :	PAGE 1/1

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







COPIE ORIGINAL F

EMBALLAGE**PACKING**

ACCESSOIRES INCLUS / ENCLOSED ACCESSORIES

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If you have any remark or suggestion about this manual or its content, please feel free to contact us by fax, mentioning the manual and equipment references.

.....

.....

.....

.....