


# Microfuge® 18 Centrifuge

**Instruction Manual** 



Symbol Simbolo Symbol 記号 Symbole 符号 Símbolo	Title / Titel / Titre / Titulo / Titolo / 名称 / 名称
4	Dangerous voltage Gefährliche elektrische Spannung Courant haute tension Voltaje peligroso Pericolo: alta tensione 危険電圧 危险电压
<u></u>	Attention, consult accompanying documents Achtung! Begleitpapiere beachten! Attention, consulter les documents joints Atención, consulte los documentos adjuntos Attenzione: consultare le informazioni allegate 注意、添付資料を参照のこと 注意,请参阅附帯的文件
	On (power) Ein (Netzverbindung) Marche (mise sous tension) Encendido Acceso (sotto tensione) 入(電源) 开 (电源)
	Off (power) Aus (Netzverbindung) Arrêt (mise hors tension) Apagado Spento (fuori tensione) 切(電源) 关 (电源)
	Protective earth (ground) Schutzleiteranschluß Liaison à la terre Puesta a tierra de protección Collegamento di protezione a terra 保護アース(接地)
<u></u>	Earth (ground) Erde Terre Tierra Scarica a terra アース(接地)



This safety notice summarizes information basic to the safe operation of the equipment described in this manual. The international symbol displayed above is a reminder that all safety instructions should be read and understood before installation, operation, maintenance, or repair of this centrifuge. When you see the symbol on other pages, pay special attention to the safety information presented. Observance of safety precautions will also help to avoid actions that could damage or adversely affect the performance of the centrifuge.

#### Safety During Installation and/or Maintenance

Any servicing of this equipment that requires removal of any covers can expose parts which involve the risk of electric shock or personal injury. Make sure that the power switch is off and the centrifuge is disconnected from the main power source, and refer such servicing to qualified personnel.

Do not replace any centrifuge components with parts not specified for use on this instrument.

#### **Electrical Safety**

To reduce the risk of electrical shock, this equipment uses a three-wire electrical cord and plug to connect the centrifuge to earth-ground. To preserve this safety feature:

- Make sure that the matching wall outlet receptacle is properly wired and earth-grounded. Check that the line voltage agrees with the voltage listed on the name-rating plate affixed to the centrifuge.
- Never use a three-to-two wire plug adapter.
- Never use a two-wire extension cord or a two-wire non-grounding type of multipleoutlet receptacle strip.

Do not place containers holding liquid on or near the chamber door. If they spill, liquid may get into the centrifuge and damage electrical or mechanical components.

#### Safety Against Risk of Fire

Certain electrical circuits within this equipment are protected by fuses against overcurrent conditions. For continued protection against the risk of fire, replace only with the same type and rating specified.

This centrifuge is not designed for use with materials capable of developing flammable or explosive vapors. Do not centrifuge such materials (such as chloroform or ethyl alcohol) in this centrifuge nor handle or store them within the 30-cm (1-ft) area surrounding the centrifuge.

#### **Mechanical Safety**

For safe operation of the equipment, observe the following:

- Use only the rotors and accessories designed for use in this centrifuge.
- Before starting the centrifuge, make sure that the rotor tie-down screw is securely fastened.
- Do not exceed the maximum rated speed of the rotor in use.
- NEVER attempt to slow or stop the rotor by hand.
- Do not lift or move the centrifuge while the rotor is spinning.
- NEVER attempt to override the door interlock system while the rotor is spinning.
- Maintain a 7.6-cm (3-in.) clearance envelope around the centrifuge while it is running. During operation come within the envelope only to adjust instrument controls, if necessary. Never lean on the centrifuge or place items on the centrifuge while it is operating.

#### **Chemical and Biological Safety**

Normal operation may involve the use of solutions and test samples that are pathogenic, toxic, or radioactive. Such materials should not be used in this centrifuge, however, unless all necessary safety precautions are taken.


- Observe all cautionary information printed on the original solution containers prior to their use.
- Handle body fluids with care because they can transmit disease. No known test offers complete assurance that they are free of micro-organisms. Some of the most virulent—Hepatitis (B and C) and HIV (I–V) viruses, atypical mycobacteria, and certain systemic fungi—further emphasize the need for aerosol protection. Handle other infectious samples according to good laboratory procedures and methods to prevent spread of disease. Because spills may generate aerosols, observe proper safety precautions for aerosol containment. Do not run toxic, pathogenic, or radioactive materials in this centrifuge without taking appropriate safety precautions. Biosafe containment should be used when Risk Group II materials (as identified in the World Health Organization Laboratory Biosafety Manual) are handled; materials of a higher group require more than one level of protection.
- Dispose of all waste solutions according to appropriate environmental health and safety guidelines.

It is your responsibility to decontaminate the centrifuge and accessories before requesting service by Beckman Coulter.



# Microfuge® 18 Centrifuge

**Instruction Manual** 



# **Contents**

	Page
INTRODUCTI	ON
	CertificationNScope of ManualNConventionsVCFC-Free CentrifugationVRadio InterferenceVRecycling LabelV
SECTION 1:	DESCRIPTION
	Centrifuge Function, Specifications, and Safety Features1-1Centrifuge Function1-1Specifications1-2Available Rotors1-3Safety Features1-3Ground Wire Check1-4Name Rating Plate1-4Chassis1-4Housing1-4Door1-4Rotor Chamber1-5Controls and Indicators1-5Power Switch1-5Control Panel1-5
SECTION 2:	INSTALLATION Installing the Centrifuge. 2-1 Electrical Requirements 2-2 Test Run 2-3
SECTION 3:	OPERATION  Run Procedure

	Page	
SECTION 4:	TROUBLESHOOTING Troubleshooting	
	Error Messages 4-2	
	Accessing the Rotor in Case of Power Failure	
SECTION 5:	CARE AND MAINTENANCE	
	Maintenance5-2Preventive Maintenance5-2Replacing Fuses5-2Cleaning5-3	
	Decontamination	
	Sterilization and Disinfection 5-4	
	Storage and Transport	
	Storage	
	Returning a Centrifuge	
	Replacement Parts	
	Supplies	
Warranty	Illustrations and Tab	ole:
F: 4.4		
Figure 1-1. Figure 1-2.	The Microfuge 18 Centrifuge	
Figure 2-1.	Dimensions of the Microfuge 18 Centrifuge	
Table 4-1. Table 4-2.	Troubleshooting Chart	

## Introduction

#### **CERTIFICATION**

To ensure full system quality, Beckman Coulter Microfuge<sup>®</sup> 18 series centrifuges have been manufactured in an ISO 9001 or 13485 facility. They have been designed and tested to be compliant (when used with Beckman Coulter rotors) with the laboratory equipment requirements of applicable regulatory agencies. Declarations of conformity and certificates of compliance are available at www.beckmancoulter.com.

#### **SCOPE OF MANUAL**

This manual is designed to familiarize you with the Beckman Coulter Microfuge 18 centrifuge, its functions, specifications, operation, and routine operator care and maintenance. We recommend that you read this entire manual, especially the SAFETY NOTICE and all safety-related information, before operating the centrifuge or performing instrument maintenance.

- Section 1 contains system specifications and a brief physical and functional description of the centrifuge, including the operating controls and indicators.
- Section 2 provides instructions for installing and connecting the centrifuge.
- Section 3 contains centrifuge operating procedures.
- Section 4 lists possible malfunctions, together with probable causes and suggested corrective actions.
- Section 5 contains procedures for routine operator care and maintenance, as well as a brief list of supplies and replacement parts.



If the centrifuge is used in a manner other than specified in this manual, the safety and performance of this equipment could be impaired. Further, the use of any equipment other than that recommended by Beckman Coulter has not been evaluated for safety. Use of any equipment not specifically recommended in this manual and/or the appropriate rotor manual is the sole responsibility of the user.

#### **CONVENTIONS**

Certain symbols are used in this manual to call out safety-related and other important information. These international symbols may also be displayed on the centrifuge and are reproduced and described below and on the inside of the front cover.

#### NOTES, CAUTIONS, AND WARNINGS



Used to call attention to important information that should be followed during installation, use, or servicing of this equipment



Used to indicate a potentially hazardous situation which, if not avoided, may result in minor or moderate injury and/or mechanical damage. It is also used to alert against unsafe practices.



#### WARNING

Used whenever an action or condition may potentially cause personal injury or loss of life. Mechanical damage may also result.



#### WARNING

Indicates high voltage or risk of electric shock. Refer servicing of all areas displaying either symbol to service personnel.

#### **CFC-FREE CENTRIFUGATION**



To ensure minimal environmental impact, no CFCs are used in the manufacture or operation of Microfuge 18 centrifuges.

#### RADIO INTERFERENCE

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with this instruction manual, may cause interference to radio communications. Operation of this equipment in a residential area is likely to cause interference, in which case the user will be required to correct the interference at his own expense.

#### **CANADIAN REGULATIONS**

This equipment does not exceed the Class A limits for radio noise emissions from digital apparatus as set out in the radio interference regulations of the Canadian Department of Communications.

Le présent appareil numérique n'émet pas de bruits radioélectriques dépassant les limites applicables aux appareils numériques de Classe A prescrites dans le reglement sur le brouillage radioelectrique édicté par le Ministère des Communications du Canada.

#### RECYCLING LABEL



Note: On the instrument, the triangle background is yellow rather than gray.

This symbol is required in accordance with the Waste Electrical and Electronic Equipment (WEEE) Directive of the European Union. The presence of this marking on the product indicates:

- 1) the device was put on the European market after August 13, 2005 and
- 2) the device is not to be disposed via the municipal waste collection system of any member state of the European Union.

It is very important that customers understand and follow all laws regarding the proper decontamination and safe disposal of electrical equipment. For Beckman Coulter products bearing this label please contact your dealer or local Beckman Coulter office for details on the take back program that will facilitate the proper collection, treatment, recovery, recycling and safe disposal of the device.



# **Description**

### CENTRIFUGE FUNCTION, SPECIFICATIONS, AND SAFETY FEATURES

#### **CENTRIFUGE FUNCTION**

The Beckman Coulter Microfuge 18 (Figure 1-1) is a microprocessor-controlled compact benchtop centrifuge that generates centrifugal forces required for a wide variety of applications. Together with the Beckman Coulter F241.5P rotor, designed specifically for use in this centrifuge, applications include:

- Nucleic acid plasmids and bacteriophages isolation.
- Routine processing such as sample preparations, pelleting, extractions, purifications, concentrations, phase separations, and receptor binding.
- · Virus isolation.
- Rapid sedimentation of protein precipitates, large particles, and cell debris.
- Preparation of subcellular organelles such as mitochondria, granules, and crude microsomes.
- · Cell isolation.



Figure 1-1. The Microfuge 18 Centrifuge

#### **SPECIFICATIONS**

Only values with tolerances or limits are guaranteed data. Values without tolerances are informative data, without guarantee.

Set speed
Width
Depth
Height, door closed
Height, door open
Weight
Ventilation clearances (sides and rear) 7.6 cm (3.0 in.)
Electrical requirements
100–120 VAC, 4 A, 50/60 Hz
Electrical supply
Maximum heat dissipation into room
under steady-state conditions
Noise level 0.91 m (3 ft) in front of centrifuge
Installation (overvoltage) category
Pollution degree

 $<sup>^{\</sup>ast}$  The Microfuge 18 centrifuge can be used in a humidity controlled cold room.

 $<sup>^\</sup>dagger$  Normally only nonconductive pollution occurs; occasionally, however, a temporary conductivity caused by condensation must be expected.

#### AVAILABLE ROTORS

The F241.5P, rated for 14 000 rpm, is a fixed-angle rotor with a tube angle of 45 degrees from the axis of rotation. The rotor can centrifuge up to twenty-four 1.5 to 2.0-mL reaction-vials, or 250 to 750- $\mu$ L vials in adapters. Refer to the rotor manual (MM-TB-001) for detailed information about the rotor.



Use only the F241.5P rotor in the Microfuge 18.

#### **F241.5P Rotor Specifications**



Maximum speed
Density rating at maximum speed
Relative Centrifugal Field* at maximum speed
· · · · · · · · · · · · · · · · · · ·
at $r_{\text{max}}$ (80 mm)
Conditions requiring speed reductions see RUN SPEEDS
(in rotor manual)
Number of tube cavities
Nominal tube dimensions (1.5 to 2.0-mL vials) $10 \times 31 \text{ mm}$
Nominal tube capacity
Nominal rotor capacity
Approximate acceleration time to maximum speed
(fully loaded)
Approximate deceleration time from maximum speed
(fully loaded)
Weight of fully loaded rotor 0.57 kg (1.16 lb)
Rotor material polypropylene
Lid material polycarbonate

#### **SAFETY FEATURES**

The Microfuge 18 centrifuge has been designed and tested to operate safely indoors at altitudes up to 2 000 m (6 562 ft). Safety features include the following.

Safety features include an electromechanical door-locking mechanism that prevents operator contact with spinning rotors. When the door is closed it locks automatically. It can be unlocked only by pressing the DOOR key, and opened only when the power is on and the rotor is at rest.

#### **GROUND WIRE CHECK**

There is a ground screw on the rear panel. A ground wire check can be performed using an appropriate measuring instrument.

#### NAME RATING PLATE

The name rating plate is affixed to the rear of the centrifuge. Check that the line voltage agrees with the voltage listed on this name rating plate before connecting the centrifuge. Always mention the serial number and the model number shown when corresponding with Beckman Coulter regarding your centrifuge.

#### **CHASSIS**

#### HOUSING

The housing is a sheet-steel casing. The control panel is covered by a protective overlay made of coated polycarbonate.

#### **DOOR**

The sheet-steel door has a strobe port for speed verification. A gasket system around the chamber opening assures sealing. (Centrifuge gaskets have not been designed as bioseals for aerosol containment.) An electromechanical door lock system keeps the door locked when a run is in progress and can be opened only when the rotor is stopped. In the event of a power failure, the door lock can be manually tripped for sample recovery (see Section 4, Troubleshooting).

#### ROTOR CHAMBER

The rotor chamber is made of high-impact thermoplastic for safety and easy cleaning.

#### DRIVE

The asynchronous direct-drive motor is brushless for low maintenance. A tie-down screw secures the rotor to the drive shaft. The resilient suspension ensures that loads will not be disturbed by vibration, and prevents damage to the drive shaft if an imbalance occurs during centrifugation.

#### **CONTROLS AND INDICATORS**

#### **POWER SWITCH**

A two-position rocker switch (I, on; O, off), located on the right side panel, controls electrical power to the centrifuge.

#### **CONTROL PANEL**

The control panel (Figure 1-2) is mounted at an angle on the centrifuge front for easy visibility and access. Run time and speed are selected using front-mounted knobs and displayed on the panel.






Figure 1-2. The Control Panel

#### **SPEED**


The RPM knob is used to select the run speed, up to 14 000 rpm, in 500-rpm increments. The entered speed is indicated by an LCD light bar on the display. The upper scale of the display shows the speed in rpm (revolutions per minute); the lower scale shows the equivalent rcf (relative centrifugal force).



Actual speed is indicated during the run by a sequence of lighted bars.




If the speed is set below the existing run speed while a run is in progress, the new set speed is indicated by a flashing bar and the drive decelerates to the new value.




TIME

Using the TIME knob, operation in any of three modes can be selected:


• A timed run can be set for any run time up to 30 minutes in 1-minute increments. Lighted bars count down time remaining and "min" is displayed. Runs can be stopped at any time by pressing the STOP key.



• A continuous (hold) run can be selected by turning the timer to ∞. Lighted bars count elapsed run time in 1-minute increments, and "min" is displayed. After 30 minutes the timer will stop showing elapsed time, but the run will continue until the STOP key is pressed.



• A pulse run can be selected by setting the timer to the short run position. The centrifuge will accelerate to the maximum speed when START is pressed and held, decelerating when the key is released. Lighted bars count elapsed run time in 1-second increments, and "sec" is displayed.



• An alternate pulse-style run can be selected with the timer set to the short run position. When the START key is pressed and released, the centrifuge will accelerate to the selected speed and run for 30 seconds, then decelerate to stop. Lighted bars count elapsed run time in 1-second increments, and "sec" is displayed.

**DOOR** 



When pressed, this key releases the latch mechanism and allows the door to be opened (the rotor must be stopped, the power must be on, and the key light on).

**START** 



• With the TIME control set for a timed or continuous run, pressing and releasing the switch starts the centrifuge run; the rotor accelerates to the set speed and continues until the time runs down (or the STOP switch is pressed). A deceleration in process can be stopped by pressing the START key; this will restart the centrifuge.

 With the TIME control set for a short run (pulse), pressing and holding the switch causes the installed rotor to accelerate to maximum speed and spin as long as the key is pressed. Deceleration begins when the key is released.


**STOP** 




Pressing this switch will cause the centrifuge to decelerate to a complete stop.

Selectable slow acceleration and deceleration, or slow deceleration only, can be used to double the acceleration and deceleration times to preserve sample-to-gradient interfaces.

- At the time of run set-up:
  - Pressing the STOP key (before pressing START) will select slow acceleration *and* slow deceleration profiles, and the word SOFT will be displayed.



• Pressing the STOP key a second time (before pressing START) will select slow deceleration only. The words SOFT STOP will show on the display, and deceleration time will be doubled.



- Pressing the STOP key a third time will deselect the slower start *and* stop.
- During a run, when the STOP key is lit, you can press it to terminate the run. During deceleration, you can press it again to select the slow acceleration and slow deceleration feature. The word SOFT will show on the display. You can press the STOP key again to select the slow deceleration feature. The words SOFT STOP will show on the display. Pressing the STOP key again will deselect the slower stop profile, and SOFT STOP will no longer be displayed.



## **Installation**

#### INSTALLING THE CENTRIFUGE



#### WARNING

Do not place the centrifuge near areas containing flammable reagents or combustible fluids. Vapors from these materials could enter the centrifuge's air system and be ignited by the motor. No hazardous materials should be handled or stored within the 30-cm (1-ft) boundary surrounding the centrifuge. Maintain a 7.6-cm (3-in.) clearance envelope around the centrifuge while it is running. No persons should be within this clearance envelope while the centrifuge is operating, except to adjust centrifuge controls, if necessary.

Carefully remove the centrifuge and accessories from the shipping container. Save the container and packing materials for possible future relocation or storage.

- Select a location away from heat-producing laboratory equipment, with sufficient ventilation to allow for heat dissipation.
- Position the centrifuge on a level surface, such as a sturdy table or laboratory bench that can to support the weight of the centrifuge (13 kg/28.7 lb) and resist vibration.
- In addition to space for the centrifuge itself (see Figure 2-1 for dimensions), allow 7.6-cm (3-in.) clearances at the sides and back of the centrifuge to ensure sufficient air circulation. The centrifuge must have adequate air ventilation to ensure compliance to local requirements for vapors produced during operation.




Figure 2-1. Dimensions of the Microfuge 18 Centrifuge

• Ambient temperatures during operation should not be lower than 4°C (39°F) or higher than 35°C (95°F). Relative humidity should not exceed 80% (noncondensing).



During transport between areas with varying temperatures, condensation may occur inside the centrifuge. Allow sufficient drying time before running the centrifuge.

### **ELECTRICAL REQUIREMENTS**

100-V centrifuge	90–110 VAC, 4 A, 50/60 Hz
120-V centrifuge	108–132 VAC, 4 A, 60 Hz
230-V centrifuge	207–253 VAC, 2 A, 50 Hz
Power cord	a 2.5-m (8-ft) power cord with grounded
	plug is supplied with the centrifuge

To reduce the risk of electrical shock, this centrifuge uses a three-wire electrical cord and plug to connect the centrifuge to earth-ground. To preserve this safety feature:

- Make sure that the matching wall outlet receptacle is properly wired and earth-grounded. Check that the line voltage agrees with the voltage listed on the name rating plate affixed to the centrifuge. Then plug in both ends of the centrifuge power cord.
- Never use a three-to-two wire plug adapter.
- Never use a two-wire extension cord or a two-wire non-grounding type of multiple-outlet receptacle strip.
- If there is any question about voltage, have a qualified service person measure it under load while the drive is operating.

To ensure safety the centrifuge should be wired to a remote emergency switch (preferably outside the room where the centrifuge is housed, or adjacent to the exit from that room), in order to disconnect the centrifuge from the main power source in case of a malfunction.

#### **TEST RUN**

We recommend that you make a test run to ensure that the centrifuge is in proper operating condition following shipment. See Section 3 for instructions on operating the centrifuge.

After completing the test run, return the preaddressed warranty card included with this literature. This will validate the centrifuge warranty and ensure your receipt of further information regarding new accessories and/or modifications as they become available.



# **Operation**



#### WARNING

Handle infectious samples according to good laboratory procedures and methods to prevent spread of disease. Because spills, operator error, or tube failure may generate aerosols, observe proper safety precautions for aerosol containment.

Do not run toxic, pathogenic, or radioactive materials in this centrifuge without taking appropriate safety precautions. Biosafe containment should be used when Risk Group II materials (as identified in the World Health Organization *Laboratory Biosafety Manual*) are handled; materials of a higher group require more than one level of protection.



#### WARNING

The centrifuge must not be used in the vicinity of flammable liquids or vapors, and such materials should not be run in the centrifuge. Never bring any flammable substances within the 30-cm (1-ft) boundary surrounding the centrifuge. Do not lean on the centrifuge or place items on the centrifuge while it is operating. During operation you should not come within the 7.6-cm (3-in.) clearance envelope except to adjust centrifuge controls, if necessary.

#### **RUN PROCEDURE**



The Microfuge 18 centrifuge and F241.5P rotor were developed, manufactured, and tested for safety and reliability as part of a Beckman Coulter centrifuge/rotor system. Safety or reliability of the centrifuge cannot be assured if used with any other rotor. Use of any equipment not specifically recommended in this manual is the sole responsibility of the user.

#### PREPARATION AND LOADING

Refer to the rotor manual for instructions on preparing the rotor for centrifugation.

#### Action

#### Result

- 1. Press the power switch to on (I).
- Power is applied to the centrifuge; all LEDs on the control panel light up momentarily.
- 2. Press **DOOR** and lift the door up.

The door will remain open.




#### CAUTION

Do not drop the rotor onto the drive shaft. The shaft can be bent if the rotor is forced sideways or dropped onto it. Install the rotor by centering it over the shaft and carefully lowering it straight down.

#### Action

#### Result

3. Install the rotor according to instructions in the rotor manual. *Always run the rotor with a balanced load*.



4. Secure the rotor to the drive shaft with the tie-down screw by holding the rotor with one hand while turning the T-handle rotor wrench to the right (clockwise).



### IIII ■ NOTE

If the rotor is left in the centrifuge between runs, make sure it is seated on the drive shaft and the tie-down screw is tight before each run. Approximately every 20 runs, or once a day, loosen the tie-down screw and retighten it to ensure proper connection between the rotor and the shaft.


- 5. Install the rotor lid and secure it by turning the knob clockwise.
- 6. Close the centrifuge door and push firmly down on it until you hear the latch engage.

#### STARTING A TIMED OR HOLD RUN


Action

Result

1. Turn the **RPM** knob to the selected speed (up to 14 000 rpm).




Turn the **TIME** knob to the required run time (up to 30 minutes) for a timed run or to ∞ for a hold run.




3. Check that the door is shut and securely latched.

Selectable slow acceleration and deceleration, or slow deceleration only, can be used to double the acceleration and deceleration times to preserve sample-to-gradient interfaces.

- At the time of run set-up:
  - Pressing the STOP key (before pressing START) will select slow acceleration *and* slow deceleration profiles, and the word SOFT will be displayed.



• Pressing the STOP key a second time (before pressing START) will select slow deceleration only. The words SOFT STOP will show on the display, and deceleration time will be doubled.



• Pressing the STOP key a third time will deselect the slower start *and* stop.


Action

Result

• During a run, when the STOP key is lit, you can press it to terminate the run. During deceleration, you can press it again to select the slow acceleration and slow deceleration feature. The word SOFT will show on the display. You can press the STOP key again to select the slow deceleration feature. The words SOFT STOP will show on the display. Pressing the STOP key again will deselect the slower stop profile, and SOFT STOP will no longer be displayed.

4. Press **START**.


• Timed run—Lighted bars count down time remaining, in 1-minute increments, and "min" is displayed. Lighted bars on the speed scale indicate rotor speed (rpm × 1000 on upper scale, rcf × 1000 on lower scale).





Do not lift or move the centrifuge while the rotor is spinning.

 Continuous (hold) run—Lighted bars count up elapsed run time in 1-minute increments and "min" is displayed. (After 30 minutes the timer will stop showing elapsed time, but the run will continue until the STOP key is pressed.) Lighted bars on the speed scale indicate rotor speed (rpm × 1000 on upper scale, rcf × 1000 on lower scale).



Action	Result		
	WARNING		
	Do not attempt to override the door interlock system while the rotor is spinning.		
	A timed run ends automatically when the set time counts down to zero. You can end a hold or timed run in progress for any reason b pressing the <b>STOP</b> key.		
5. After the rotor stops spinning press <b>DOOR</b> .	The door latch is released. Open the door and unload the rotor.		
	IIII NOTE		
	Depending upon the duration of the event, transient power-line interruptions may cause this equipment to decelerate or reset with the possible loss of in-process operation. If you suspect this has occurred, simply repeat the runs.		
STARTING A PULSE R	UN		
STARTING A PULSE R	UN  Two pulse-run modes are available.		
STARTING A PULSE R			

and securely latched.

2. Turn the **TIME** knob to

the short run position.

#### Action

#### Result

3. Press and hold the **START** key.

The centrifuge accelerates to the maximum speed and continues as long as the key is pressed, decelerating when the key is released. Lighted bars count elapsed run time in 1-second increments, and "sec" is displayed.



#### CAUTION

Do not lift or move the centrifuge while the rotor is spinning.

4. After the rotor stops spinning press **DOOR**.

The door latch is released. Open the door and unload the rotor.



WARNING

Do not attempt to override the door interlock system while the rotor is spinning.

#### Pulse Mode 2


This pulse-run mode will accelerate to the set speed when START is pressed and run at that speed for 30 seconds, then decelerate to stop.

#### Action

#### Result

- 1. Turn the **RPM** knob to the selected speed (up to 14 000 rpm).
- 2. Turn the **TIME** knob to the short run position.

A lighted bar is at the "short run" position, and "sec" is displayed.




Action

Result

- 3. Check that the door is shut and securely latched.
- 4. Press and release the **START** key.

The centrifuge accelerates to the selected speed and continues for 30 seconds, then decelerates to stop. Lighted bars count elapsed run time in 1-second increments, and "sec" is displayed.





CAUTION

Do not lift or move the centrifuge while the rotor is spinning.

5. After the rotor stops spinning press **DOOR**.

The door latch is released. Open the door and unload the rotor.



WARNING

Do not attempt to override the door interlock system while the rotor is spinning.

#### **UNLOADING**

After completing a run, unload the rotor according to instructions in the rotor manual.



CAUTION

If disassembly reveals evidence of leakage, you should assume that some fluid escaped the rotor. Apply appropriate decontamination procedures to the centrifuge and accessories.



## **Troubleshooting**

IIII NOTI	Ε
-----------	---

It is your responsibility to decontaminate the centrifuge, as well as any rotors and accessories, before requesting service by Beckman Coulter representatives.

#### TROUBLESHOOTING

Malfunctions that may occur are described in Table 4-1, along with probable causes and corrective actions required. Possible causes for each problem are listed in the probable order of occurrence. Perform the recommended corrective action in sequence, as listed. If you are unable to correct the problem, call your local Beckman Coulter representative.


To help diagnose and correct the problem, try to gather as much information about the situation as you can:

- Note the operating situation when the error occurred (such as speed or load type).
- Note any unusual environmental and/or operating conditions (such as ambient temperature or voltage fluctuations).
- Add any other information that may be helpful.

Table 4-1. Troubleshooting Chart

Problem	Probable Cause	Recommended Action
No indication on displays	Power not on	Plug in power cord; turn power on (I).
	Fuse blown	Change both fuses as described in Section 5, Care and Maintenance.
	Mechanical failure	Call Beckman Coulter Field Service.
Centrifuge cannot be started	Power not on	Plug in power cord; turn power on (I).
Starteu	Fuse blown	Change both fuses as described in Section 5, Care and Maintenance.
	Mechanical failure	Call Beckman Coulter Field Service.
DOOR key light blinks and centrifuge will not start	Door not closed properly	Press the DOOR key to open the door, then close the door again, pressing down firmly to ensure both latches engage.
Rotor cannot achieve set	Line voltage below rating	Measure line voltage while centrifuge is operating.
speed	Electrical failure	Check connections; call Beckman Coulter Field Service.
	Mechanical failure	Call Beckman Coulter Field Service.
Door will not open	Rotor spinning	Wait until rotor stops.
	Power not on	Plug in power cord; turn power on (I). If power cannot be restored, see ACCESSING THE ROTOR IN CASE OF POWER FAILURE, below, to retrieve your sample.

#### **ERROR MESSAGES**



This would indicate Error Number 3.

If a problem occurs during operation, the three icon lights (START, DOOR, and STOP) will flash and a lighted bar on the TIME display will point to the error number. Refer to Table 4-2 to determine the nature of the problem and recommended actions. If you are unable to correct the problem, call your Beckman Coulter Field Service representative.

#### ACCESSING THE ROTOR IN CASE OF POWER FAILURE

If the facility power fails you will have to restart the run when the power is restored. In the event of an extended power failure, it may be necessary to trip the door-locking mechanism manually to remove the rotor and retrieve your sample.

Table 4-2. Error Message Chart.

If the recommended action does not correct the problem, call Beckman Coulter Field Service.

Error Number	Description	Recommended Action
1	Bad tachometer reading	Turn the power off (O) and wait for the rotor to come to a complete stop, then turn the power back on (I) to reset.
2	No tachometer signal	complete stop, then turn the power back on (i) to reset.
3	Drive is not enabled	
4	Door does not open after door relay is active	Turn the power off (O). Verify that the rotor is completely stopped, then open the door using the procedures under ACCESSING THE ROTOR IN CASE OF POWER FAILURE, below. Close the door firmly, and turn the power back on (I).
5	Door is open and rotor is spinning.	Press STOP and wait for the rotor to come to a complete stop. Turn the power off (O), open and close the door, then turn the power back on (I) to reset.
6	MFR signal failure	Turn the power off (O) and wait for the rotor to come to a complete stop, then turn the power back on (I) to reset.
7	Overspeed condition	Press STOP and wait for the rotor to come to a complete stop.
8	CPU failure	Turn the power off (O), wait 3 minutes, then turn the power back on (I) to reset.
9	EPROM bad checksum	
10	Software problem	
11	Wrong jumper setting	



#### WARNING

The following procedure may expose the operator to the possibility of contact with a spinning rotor. Turn the power off (O) and disconnect the centrifuge from the main power source before proceeding. Never attempt to override the door interlock system while the rotor is spinning.

## Action

## Result

Turn the power switch to off
 (O) and unplug the power
 cord from the power source.

No indicators are lit.

#### Action


- 2. Make sure that the rotor is not spinning.
- 3. Tilt the centrifuge back and locate the emergency access device at the bottom front of the front panel.
- 4. Turn the screw to the right (clockwise) until you can grip the release device with your fingers.
- 5. Pull down on the device to release the latch. After the latch is released, reinsert the device.

6. Open the door and remove

the rotor.

## Result

No sound or vibration comes from the centrifuge.



The latch releases and the door can be opened.

If the rotor is still spinning, close the door and wait until it stops before attempting to remove it.





# **Care and Maintenance**

For maintenance not covered in this manual, contact Beckman Coulter Field Service (1-800-742-2345 in the U.S.A.; customers outside the United States should contact their local Beckman Coulter representative). Refer to the rotor manual and Chemical Resistances (publication IN-175) for instructions on the care of rotors and rotor accessories.



It is your responsibility to decontaminate the centrifuge, as well as any rotors and accessories, before requesting service by Beckman Coulter representatives.



## WARNING

Any maintenance procedure or servicing of this equipment that requires removal of any covers can expose parts which involve the risk of electric shock or personal injury. Make sure that the power switch is off (O) and the centrifuge is disconnected from the main power source, and refer such servicing to qualified service personnel.

Do not use alcohol or other flammable substances in or near operating centrifuges.

## **MAINTENANCE**

#### PREVENTIVE MAINTENANCE

Perform the following procedures regularly to ensure satisfactory performance and long service life of the centrifuge.

- Regularly inspect the interior of the rotor chamber for accumulations of sample or dust. Clean as required (see CLEANING, below), as these accumulations can result in rotor vibrations.
- Regularly check the air intake and exhaust vents for obstructions. Keep vents clear and clean.



Before using any cleaning or decontamination methods except those recommended by the manufacturer, users should check with the manufacturer that the proposed method will not damage the equipment.

#### REPLACING FUSES

The fuse holder is located on the back of the centrifuge, below the power cord receptacle.



#### CALITION

For continued protection against the risk of fire, replace fuses only with the listed fuses (T2A, 365639; T4A, 367200).

#### Action

#### Result

Turn the power switch to off
 (O) and unplug the power
 cord from the power source.

No indicators are lit.

#### Action

#### Result


2. Use your fingernails or a small screwdriver to press in on the fuse holder latches.

Latches release and the fuse holder can be pulled out.

3. Pull the fuse holder straight out.

Fuses are accessible.

4. Pull the fuse straight out of the fuse holder.



- 5. Carefully insert a new fuse. Press it straight in until it seats in the fuse holder.
- 6. Repeat steps 4 and 5 for the other fuse.
- 7. Position the fuse holder in the opening and press it in until it seats.

### **CLEANING**



Frequent cleaning will ensure proper operation and prolong the life of the centrifuge. Always clean up spills when they occur to prevent corrosives or contaminants from drying on component surfaces.\*

- To prevent accumulations of sample and/or dust, keep the interior of the rotor chamber clean and dry by frequent wiping with a cloth or paper towel.
- Clean the drive shaft, shaft cavity, threads, and the tie-down screw
  at least once a week using a mild detergent such as Beckman
  Solution 555<sup>TM</sup> (339555) and a soft brush. Dilute the detergent 10
  to 1 with water. Rinse thoroughly and dry completely. Lubricate
  the drive shaft with Spinkote (306812) after cleaning.
- Wash the bowl using a mild detergent such as Solution 555. Rinse thoroughly and dry completely.
- Clean the centrifuge case and door by wiping with a cloth dampened with Solution 555. Do not use acetone or other solvents.

#### **DECONTAMINATION**





If the centrifuge and/or accessories are contaminated with radioactive or pathogenic solutions, perform appropriate decontamination procedures. Refer to *Chemical Resistances* to be sure the decontamination method will not damage any part of the centrifuge.

## STERILIZATION AND DISINFECTION

Ethanol (70%)<sup>†</sup> may be used on the centrifuge surface. See *Chemical Resistances* for more information regarding chemical resistance of centrifuge and accessory materials.

While Beckman Coulter has tested these methods and found that they do not damage the centrifuge, no guarantee of sterility or disinfection is expressed or implied. When sterilization or disinfection is a concern, consult your laboratory safety officer regarding proper methods to use.

<sup>\*</sup> Before using any cleaning methods except those recommended by the manufacturer, users should check with the manufacturer that the proposed method will not damage the equipment.

<sup>†</sup> Flammability hazard. Do not use in or near operating centrifuges.

## STORAGE AND TRANSPORT

#### **STORAGE**

Before storing a centrifuge for an extended period, return it to the original shipping container to protect it from dust and dirt. Temperature and humidity conditions for storage should meet the environmental requirements described under SPECIFICATIONS, in Section 1.

#### RETURNING A CENTRIFUGE



Before returning a centrifuge or accessory for any reason, prior permission (a Returned Goods Authorization form) must be obtained from Beckman Coulter, Inc. Contact your local Beckman Coulter office to obtain the RGA form and for packaging and shipping instructions.

To protect our personnel, it is the customer's responsibility to ensure that all parts are free from pathogens and/or radioactivity. Sterilization and decontamination must be done before returning the parts.

All parts must be accompanied by a signed note, plainly visible on the outside of the box, stating that they are safe to handle and that they are not contaminated with pathogens or radioactivity. Failure to attach this notification will result in return or disposal of the items without review of the reported problem.

## **SUPPLY LIST**



Publications referenced in this manual can be obtained by calling 1-800-742-2345 in the United States, or by contacting your local Beckman Coulter office.

Contact beckman Coulter Sales (1-800-742-2345 in the United States; worldwide offices are listed on the inside back cover of this manual) for information about ordering parts and supplies. See the Beckman Coulter *Benchtop Rotors, Tubes & Accessories* catalog (BR-9742, available at www.beckmancoulter.com) for detailed information on

ordering rotors, tubes, and accessories. For your convenience, a partial list is given below. Refer to the rotor manual for materials and supplies needed for rotors.

## REPLACEMENT PARTS

F241.5P rotor
Rotor tie-down screw
T-handle wrench
Fuse, T2A, 230 V
Fuse, T4A, 100–120 V

## **SUPPLIES**

Spinkote lubricant (2 oz)	306812
Silicone vacuum grease (1 oz)	335148
Beckman Solution 555 (1 qt)	339555

## MICROFUGE® SERIES CENTRIFUGE WARRANTY

Subject to the exceptions and upon the conditions specified below and the warranty clause of the Beckman Coulter terms and conditions in effect at the time of sale, Beckman Coulter agrees to correct either by repair or, at its election, by replacement, any defects of material or workmanship which develop within one (1) year after delivery of a Microfuge® series centrifuge (the product), to the original buyer by Beckman Coulter or by an authorized representative, provided that investigation and factory inspection by Beckman Coulter discloses that such defect developed under normal and proper use.

Some components and accessories by their nature are not intended to and will not function for as long as one (1) year. A complete list of such components or accessories is maintained at the factory and at each Beckman Coulter District Sales Office. The lists applicable to the products sold hereunder shall be deemed to be part of this warranty. If any such component or accessory fails to give reasonable service for a reasonable period of time, Beckman Coulter will repair or, at its election, replace such component or accessory. What constitutes either reasonable service and a reasonable period of time shall be determined solely by Beckman Coulter.

### REPLACEMENT

Any product claimed to be defective must, if requested by Beckman Coulter, be returned to the factory, transportation

charges prepaid, and will be returned to Buyer with the transportation charges collect unless the product is found to be defective, in which case Beckman Coulter will pay all transportation charges.

#### **CONDITIONS**

Beckman Coulter shall be released from all obligations under all warranties, either expressed or implied, if the product(s) covered hereby are repaired or modified by persons other than its own authorized service personnel, unless such repair in the sole opinion of Beckman Coulter is minor, or unless such modification is merely the installation of a new Beckman Coulter plug-in component for such product(s).

#### **DISCLAIMER**

IT IS EXPRESSLY AGREED THAT THE ABOVE WARRANTY SHALL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND OF THE WARRANTY OF MERCHANTABILITY AND THAT NEITHER BECKMAN COULTER, INC. NOR ITS SUPPLIERS SHALL HAVE ANY LIABILITY FOR SPECIAL OR CONSEQUENTIAL DAMAGES OF ANY KIND WHATSOEVER ARISING OUT OF THE MANUFACTURE, USE, SALE, HANDLING, REPAIR, MAINTENANCE, OR REPLACEMENT OF THE PRODUCT.



Beckman Coulter, Inc. • 250 S. Kraemer Blvd. • Brea, California 92821 Sales and Service: 1-800-742-2345 • Internet: www.beckmancoulter.com

