

TGGE MAXI System

(order number 024-200)

Manual

Version 3.02

October 2000



!! Warning !!
Please read this manual carefully before
using the apparatus



Biometra
biomedizinische Analytik GmbH
Rudolf-Wissell-Str. 30, D-37079 Goettingen
Tel.: 0551/50 686-0; Fax: 0551/50 686-66
email: info@biometra.co.uk

<http://www.biometra.com>

1. WHAT'S NEW IN THIS MANUAL	2
1.1 SAFETY INSTRUCTIONS / GENERAL REMARKS.....	4
2. INTRODUCTION.....	5
2.1 PRINCIPLE OF THE METHOD.....	5
2.2 SPECIAL FEATURES OF THE BIOMETRA TGGE SYSTEM	6
3. BEFORE YOU START: THE TGGE FLOW CHART OF SUCCESS.....	7
4. THE TGGE MAXI SYSTEM	9
4.1 SYSTEM OVERVIEW.....	9
4.2 INSTALLATION	9
4.3 ADAPTATION OF PROTOCOLS FROM THE TGGE "MINI" SYSTEM.....	9
5. CASTING OF GELS	10
5.1 ASSEMBLY OF THE GEL CUVETTE	10
5.2 PREPARING GEL SOLUTION.....	12
5.2.1 <i>Preparing gel solution for TBE buffer system</i>	12
5.2.2 <i>Preparing gel solution for MOPS / EDTA buffer system</i>	12
5.3 POURING GELS	13
6. ELECTROPHORESIS.....	14
6.1 ELECTROPHORESIS CONDITIONS.....	14
6.2 PRE-RUN FOR SAMPLE LOADING AND TEMPERATURE EQUILIBRATION.....	14
6.3 SETUP ELECTROPHORESIS UNIT.....	15
6.3.1 <i>Prepare prior to assembly of the electrophoresis unit</i>	15
6.3.2 <i>Gel setup for electrophoresis</i>	17
6.4 PERPENDICULAR TGGE.....	18
6.5 PARALLEL TGGE	19
6.6 HOW TO IDENTIFY THE OPTIMUM TEMPERATURE RANGE FROM A PERPENDICULAR GEL.....	19
7. PROGRAMMING THE TGGE CONTROLLER.....	21
7.1 CREATE / EDIT PROGRAM.....	21
7.2 SELECT PROGRAM	21
7.3 NAME PROGRAM	22
7.4 ENTER TEMPERATURES FOR THE GRADIENT BLOCK.....	22
7.5 ENTER ELECTROPHORESIS PARAMETERS	23
7.6 START ELECTROPHORESIS	24
7.7 STOP/PAUSE ELECTROPHORESIS	25
7.8 VIEW TEMPERATURES OF THE GRADIENT	25
7.9 SPECIAL FUNCTIONS.....	26
7.9.1 <i>Print programs</i>	26
7.9.2 <i>Select / de-select signal</i>	26
7.9.3 <i>Select language</i>	26
7.9.4 <i>Set block type</i>	26
8. STAINING.....	28
8.1 SILVER STAINING	28
8.1.1 <i>Silver staining protocol</i>	28
8.2 ETHIDIUM BROMIDE-STAINING.....	29
8.3 AUTORADIOGRAPHY	29
8.4 ELUTION OF DNA FROM THE TGGE GEL.....	30
9. SAMPLE PREPARATION.....	31
9.1.1 <i>Purity of samples</i>	31
9.1.2 <i>Sample preparation for direct DNA analysis</i>	31
9.1.3 <i>Denaturation / Renaturation for heteroduplex analysis of DNA</i>	31
10. OPTIMIZATION OF TGGE.....	32

10.1	DESIGN OF DNA FRAGMENT FOR TGGE	32
10.1.1	<i>Poland analysis</i>	32
10.1.2	<i>GC clamps</i>	35
10.1.3	<i>Chemical clamp with Psoralen (Furo[3,2-g]coumarin, C₁₁H₆O₃)</i>	35
10.1.4	<i>Use of SSCP primers</i>	35
10.2	FIND CORRECT TEMPERATURE GRADIENT.....	35
10.3	PARALLEL ANALYSIS OF MULTIPLE SAMPLES	36
10.4	OPTIMIZATION OF PARALLEL TGGE	37
10.5	WHAT TO DO, IF DIFFERENT SAMPLES HAVE VERY SIMILAR MELTING POINTS: HETERODUPLEX ANALYSIS WITH TGGE	38
10.5.1	<i>Principle of heteroduplex analysis</i>	38
10.5.2	<i>Evaluation of a heteroduplex analysis</i>	40
11.	THE TGGE TEST KIT (ORDER NUMBER 024-050)	41
11.1	PERPENDICULAR TGGE USING THE BIOMETRA TGGE TEST KIT	41
11.2	PARALLEL TGGE USING THE BIOMETRA TGGE TEST KIT	42
12.	TECHNICAL SPECIFICATION.....	43
13.	ORDERING INFORMATION:.....	44
14.	TROUBLE-SHOOTING	45
15.	REFERENCES.....	47
16.	APPENDIX.....	50
16.1	BUFFERS	50
17.	INSTRUCTIONS FOR RETURN SHIPMENT.....	52
18.	EQUIPMENT DECONTAMINATION CERTIFICATE	53
19.	WARRANTY	54

1. What's new in this manual

The current version of the TGGE maxi manual reflects the strong progress that has been made with the TGGE method. New insights concerning handling and separation conditions have been integrated in the proven step by step protocol.

The major news are:

- 1) The volume of thermal coupling solution is a critical factor for an even migration front and should be held as small as possible (section 6.3.2)
- 2) Overlaying the gel with buffer should be omitted (section 6.3.2)
- 3) The gel cover film should be applied right after putting the gel on the block (section 6.3.2)
- 4) Identical salt concentration in all samples is important for an even migration front.
- 5) Instead of buffer wicks, household sponges may be used (section 6.3.2)
- 6) Extended Silver staining protocol for optimum gels

The **TGGE** method is covered by patents issued to Diagen (now QIAGEN GmbH).

The polymerase chain reaction (**PCR**) process is covered by patents issued to Hoffman-La Roche.

Acryl-Glide is a trademark of Amresco Inc.

Biometra is a trademark of Biometra GmbH.

Whatman is a trademark of Whatman International Ltd.

The **POLAND** software service established by Gerhard Steger, Department of Biophysics, University of Duesseldorf, is available by internet

www.biophys.uni-duesseldorf.de/POLAND/poland.html.

Please read the TGGE MAXI manual carefully before starting operation



1.1 Safety instructions / general remarks

- Do not fill buffer chambers above marking for maximum level
- If buffer has been spilled on the electrophoresis unit, clean it carefully before start of electrophoresis
- Never run device without gel cover plate
- The thermoblock is covered with a Teflon film. Avoid damaging this film.
- For cleaning of the thermoblock do not use aggressive chemicals or strong detergents.
- Do not use paraffin oil on the thermoblock.
- In case of strong condensation under the safety lid stop run, dry instrument and re-start
- Switch off power before removing the safety lid
- Do not move instrument during operation
- Do not lift the electrophoresis unit by holding it on the white frame. Instead, lift the lower part of the device (red corpus).

2. Introduction

Temperature Gradient Gel Electrophoresis is a powerful technique for the separation of nucleic acids or proteins. The TGGE method, which is covered by patents, uses the temperature dependent changes of conformation for separating molecules (for review see Reference 1).

Since the introduction of the first commercial available TGGE apparatus in 1989, temperature gradient gel electrophoresis has gained high interest in scientific and clinical research laboratories due to the unprecedented resolution capability and easiness of analysis. The range of scientific publications using the TGGE method is broad and covers all disciplines which use molecular biology methods: e.g. Oncology²⁻⁴, Virology^{5,6}, Immunology^{7,8}, RNA Viroid Research⁹⁻¹², Prion Research¹³, Population Analysis¹⁴⁻¹⁵. The TGGE method has also been used for quantitative analysis in industry¹⁶⁻¹⁷ and for conformational analysis of proteins¹⁸⁻¹⁹.

2.1 Principle of the method

Conventional protein or nucleic acid electrophoresis separates molecules according to their size or charge. **TGGE** adds a new parameter for separation, namely the **melting behavior of a molecule**. The melting behavior is determined by primary sequence and secondary and tertiary structure of the molecule and can be changed by external influences like temperature, salt concentration, pH etc.

During electrophoresis the sample migrates along a temperature gradient. As the temperature rises the molecules start to denature. Working with PCR fragments for example electrophoresis starts with double stranded molecules. At a certain temperature the DNA starts to melt, resulting in a fork-like structure (partial single strand, see Figure 1). In this conformation the migration is slowed down compared to a completely double stranded DNA fragment (of same size). Since the melting temperature strongly depends on the base sequence, DNA fragments of same size but different sequence can be separated. This is used in mutation detection where PCR fragments of identical size but different sequence are separated. Thus TGGE not only separates molecules but gives additional information about melting behavior and stability.




Figure 1: Different conformations of DNA during temperature gradient gel electrophoresis.

2.2 Special features of the Biometra TGGE System

The most powerful characteristic of the Biometra TGGE system is the high reproducibility of the temperature gradient. In contrast to conventional systems using chemical gradients (DGGE) the temperature gradient of the Biometra TGGE system establishes the same denaturing gradient time after time after time. The microprocessor driven gradient block of the TGGE System allows strictly defined linear temperature gradients with high resolution and reproducibility.

Because of the small amount of material used for separation, DNA or RNA fragments appear as fine bands which can be clearly distinguished from each other. Even complex band patterns can be analyzed due to the high resolution capability of the gradient block. Comparing the TGGE method with other screening methods like SSCP the superior performance of the TGGE method becomes evident²⁰⁻²².

The Biometra TGGE system is available in two formats. The standard TGGE “mini” system (024-000) operates small gels and is therefore ideally suited for fast, serial experiments. The TGGE maxi system (024-200) provides a large separation distance and allows high parallel sample throughput.

Using the Biometra TGGE system it is very easy to separate samples either **parallel** or **perpendicular** to a temperature gradient. All that has to be changed is the position of the buffer tanks. Whereas perpendicular TGGE is mainly used for the optimization of separation conditions, parallel TGGE allows fast analysis of multiple samples.

PERPENDICULAR TGGE	<p>Temperature gradient is perpendicular to the electrophoretic migration:</p> <p>One sample is separated over a broad temperature range to determine the optimum temperature gradient or to analyze temperature dependent changes in conformation</p>
PARALLEL TGGE	<p>Temperature gradient is parallel to electrophoretic migration:</p> <p>multiple samples are separated in parallel</p>




Figure 2: Typical results after perpendicular TGGE (A: temperature gradient from left to right) and parallel TGGE (B: temperature gradient from top to bottom).

3. Before you start: The TGGE flow chart of success

TGGE is a powerful technique to separate molecules of same size, but different sequence. Nevertheless, every DNA fragment has its own characteristics and three steps have to be taken before successful analysis of multiple samples in parallel TGGE can begin. Each of the following steps is described in detail in section 10, Optimization of TGGE analysis.

Step 1 Check your DNA sequence in Poland analysis. The Poland computer program (<http://www.biophys.uni-duesseldorf.de/POLAND/poland.html>) calculates the melting behavior of dsDNA molecules. Poland analysis can predict, whether a fragment is suited for TGGE or not.

Melting profile is ok
Poland analysis shows a satisfying profile.

Proceed with step 2

Melting profile is not ok

If Poland analysis shows that the fragment in its current state is not suited for TGGE, optimize your primer design.

Never try to separate samples in TGGE if the calculated melting profile is not ok.

Step 2 If the Poland analysis shows a suitable melting profile you should test separation conditions in a perpendicular TGGE. In perpendicular TGGE, a large aliquot of the sample runs over a broad temperature range. The result of parallel TGGE allows identification of the temperature gradient for parallel analysis.

Perpendicular gel is ok
Perpendicular TGGE shows a nice melting curve.

Proceed with step 3

Perpendicular gel is not ok

If perpendicular analysis does not show the expected melting profile, check sequence again in Poland analysis. Also check purity of chemicals and electrophoretic conditions.

Do not try samples in parallel TGGE, if the perpendicular gel does not show a defined melting curve.

Step 3 Set up a parallel gel with the temperature gradient derived from the perpendicular gel. Separation can be optimized by varying the temperature gradient and voltage.

4. The TGGE maxi system

4.1 System overview

The TGGE MAXI system consists of three components:

- 1) **electrophoresis unit**
including thermoblock, buffer chambers, safety lid
- 2) **controller**
control of electrophoresis parameters (Voltage) and temperature gradient
- 3) **power supply**
power supply for electrophoresis unit and controller

4.2 Installation

Connect electrophoresis unit and controller.

Connect controller and power supply

Connect power cables from the safety lid to the controller

The buffer chambers can be placed in two orientations, depending on the direction of the temperature gradient. Be sure that the orientation is correct. The markings on the gradient block (L0 to L10) indicate the direction of the temperature gradient. The direction of electrophoresis (minus to plus for nucleic acids) is indicated on the safety lid.

4.3 Adaptation of protocols from the TGGE “mini” system

In contrast to the „mini“ system, the TGGE maxi system works with gels of 1mm thickness. Therefore electrophoretic parameters have to be adapted. This can be done by running a perpendicular gel or by performing a time chase experiment. Both techniques are described in detail in section 10.2 and 10.3.

With the TGGE maxi system the settings for T1 and T2 have been omitted. The temperature gradient is set by defining a temperature for L0 and a temperature for L10. The corresponding lines are marked on the thermoblock. This distance between L0 and L10 is effectively used for the separation. For the evaluation of stained gels we have included a plastic film where lines L0 to L10 are indicated. The gel is placed on the film and the position of the bands can be identified in correlation to the temperature lines (see section 6.6).


5. Casting of gels

5.1 Assembly of the gel cuvette

There are two kind of gels for parallel TGGE analysis. The standard gel contains 30 slots for 5µl of sample each. This system is described in the bellow text.

Alternatively, a gel without slots can be casted (use glass plate without slot formers, 024-227). In this case, samples are loaded with an applicator strip (024-223). This silicone applicator strip contains 32 holes and is simply placed on top of the gel. Samples are loaded in the holes. After turning on voltage samples diffuse into the gel and migrate along the gel matrix.

The standard gel cuvette consists of one glass plate with spacers and slot formers (024-228 or 024-229), cover plate without spacers (024-221) and a silicone sealing (024-230). The sandwich is fixed with 9 metal clamps.



Note: The gel is poured on a support film to optimize temperature transition between block and gel. The gel is bound covalently to this hydrophilic film (024-234).

Note: The gels sticks to the support film throughout the whole procedure, including staining.

- 1) Clean glass plates with 70% ethanol and a soft tissue.

Note : Be careful not to damage the slot formers when cleaning the glass plates. Never use strong detergents in the area of the slot formers. Avoid strong mechanical contact. Self adhesive slot forming units for replacement (024-222) are included in the starter kit.

- 2) Pre-treatment of glass plate with spacers and slot formers:
The glass plates with the slot formers are treated with a solution that makes the surface hydrophobic and therefore facilitates the removal of the gel (the gel remains on the support film)
* Apply approx. 2ml AcrylGlide (211-319) on the glass plate and spread it with a soft tissue, especially between the slot formers (alternative hydrophobic solutions may be used)
* wait for 2 minutes
* polish glass plate with a soft tissue

Note: Never apply Acryl Glide onto the spacers, because this will lead to leakage of the gel cuvette.

- 3) Place polybond film (024-234) on the glass plate **without** spacers.
- 4) Attach polybond film by carefully wiping it with a soft tissue. The film should attach uniformly to the glass plate. To improve contact between glass plate and film a drop of water may be applied to the glass plate.

Note: The support film may be fixed along the upper side of the cover glass plate with ordinary adhesive tape. This way no gel solution can accidentally get behind the support film.

- 5) Place the silicone sealing around spacers.
- 6) Assemble gel sandwich and fix it with 3 metal clamps on each side.

Note: the clamps should be placed directly on the spacers.

- 7) Set gel sandwich upright on the 3 bottom clamps (see Figure 3)




Figure 3: Final setup of the gel cuvette.

5.2 Preparing gel solution

The choice of the buffer system has a strong impact on TGGE analysis. Concentration of salt and denaturing agents (urea or formamide) strongly affects the melting temperature of DNA and proteins. In general, urea is used for the separation of nucleic acids in a concentration between 7 and 8M. Urea reduces the melting temperature and thus enables a separation at lower temperatures (which is favorable, because at higher temperatures the gel tends to dry out). To further reduce the melting temperature (deionized) formamide may be used in concentrations of up to 20%. The most popular buffer systems for TGGE are TBE, TAE and MOPS. In the following, two standard protocols for TBE and MOPS buffer systems are listed. Please note, that the buffer system should be adapted to each special kind of application.

5.2.1 Preparing gel solution for TBE buffer system

For one TGGE maxi gel prepare 50 ml gel solution

gel composition final concentrations	STOCK SOLUTION	FOR 50 ML	FOR 100 ML
Acrylamid [8%]	40 % (37,5 :1)	10 ml	20ml
urea [7M]	solid	21 g	42 g
TBE [0.1x]	1 x	5 ml	10 ml
Glycerol [2%]	40%	2,5 ml	5 ml
adjust with aqua bidest		to 50 ml	To 100ml

- Stir solution at 50°C until urea is completely dissolved.
- Carefully degas gel solution
- let cool down to room temperature and start polymerization with

APS	10 %	80 µl	160 µl
TEMED	100%	110 µl	220 µl

- Load gel solution in a syringe and attach a 0,4µm or 0,25µm sterile filter
- Pour gel through sterile filter into the glas sandwich.

5.2.2 Preparing gel solution for MOPS / EDTA buffer system

Preparation of buffers and stock solutions is described in section 16.1.

For one TGGE maxi gel prepare 50 ml gel solution

gel composition final concentrations	STOCK SOLUTION	FOR 50 ML	FOR 100 ML
Acrylamid [8%]	40 % (37,5 :1)	10 ml	20ml
urea [7M]	solid	21 g	42 g
MOPS [1 x]	50 x	1 ml	2 ml
Glycerol [2%]	40%	2,5 ml	5 ml
fill up with aqua bidest		to 50 ml	to 100ml

- Stir solution at 50°C until urea is completely dissolved.
- Carefully degas gel solution
- let cool down to room temperature and start polymerization with

APS	10 %	80 µl	160 µl
TEMED	100%	110 µl	220 µl

- Load gel solution in a syringe and attach a 0,4µm or 0,25µm sterile filter
- Pour gel through sterile filter into the glas sandwich.

5.3 Pouring gels

- 1) Pour gel solution slowly into the sandwich. Avoid bubbles!
- 2) Let polymerize for approx. 3 h at room temperature

Note: Polymerized gels may be stored for 3 days or even longer at **room temperature**. Remove clamps and wrap gel sandwich including glass plates in wet paper towels. Store in a tight plastic bag.

Note: Prepare electrophoresis unit prior to disassembling the gel cuvette. The gel should not be exposed to the air for extended periods since this may lead to drying of the gel.

- 3) After polymerization remove clamps. Remove the glass plate without spacers by sliding the glass plate away from the rest of the sandwich (if you have fixed the support film with adhesive tape, remove or cut tape first). The gel must stick on the support film.
- 4) Remove gel together with support film carefully from the glass plate with spacers. Be careful not to damage the slots.

6. Electrophoresis

6.1 Electrophoresis conditions

The electrophoresis unit of the TGGE System has been designed to accommodate TGGE and all related applications like CTGE, TTGE and SSCP without cumbersome changes. The buffer tanks can be positioned in two orientations, allowing a temperature gradient parallel or perpendicular to the electrophoresis direction (see section 6.4 and 6.5).

Electrophoresis conditions in general depend on

- the kind of sample, e.g. protein, nucleic acid, fragment size
- the kind of application, e.g. parallel or perpendicular TGGE
- the sample preparation, e.g. high salt or low salt preparation,
- the buffer system.

Any recommendations should be regarded as guidelines to start with. Further improvement of the analysis should be done by adjusting the run conditions to individual needs.

Voltage:	100 V - 400 V	start with 300 V
(Current:	5 mA - 25 mA	Approx. 10-20 mA)*
Run Time:	30 min - 4 h	Approx. 3h

***Note 1:** The controller is designed to control voltage rather than amperage (set [mA] and [Vh] to maximum values).

Note 3: To minimize run time in a parallel gel, start with a temperature right below the melting temperature of the fragment.

6.2 Pre-run for sample loading and temperature equilibration

Let samples migrate into the gel at homogenous temperature. For this purpose a multistep program can be set up in the controller:

Step 1:	10 minutes, 300Volt, 20°C (L0 and L10)	sample migrate from slots into gel
Step 2:	10 minutes, 0Volt, temperature gradient	temperature gradient equilibrates
Step 3:	main run	samples are separated

For creating a multi step program see section 7.5.

6.3 Setup electrophoresis unit

The Biometra TGGE system is a horizontal electrophoresis system. The buffer bridges to the gel are established by layering one side of each buffer wick (024-215) on the gel and submerging the other in the buffer inside the tank. To protect the gel from drying, it is covered with a gel cover film (024-232). The complete setup, consisting of gel with cover film and buffer wicks, is covered with the gel cover plate. The gel cover plate has two sealings and fits tightly onto the thermoblock. It holds the buffer wicks in place and helps to build a humidity chamber around the gel. This is important to prevent evaporation during the run.



Important: Never run a gel without gel cover plate. (This could lead to massive condensation under the safety lid. Danger of electric shock.)




Figure 4: Set up of the gel for electrophoresis

6.3.1 Prepare prior to assembly of the electrophoresis unit

Note: Be sure to have everything on hand, to avoid extended handling times. Don't let the disassembled gel dry, during setup of the electrophoresis system.

Prepare:

- ➡ parallel or perpendicular gel (see section 5)
- ➡ samples in loading buffer (for sample preparation see chapter 9)
- ➡ 1000 ml running buffer (for recipes see chapter 16.1)
- ➡ 4 buffer wicks (024-215)

- ➡ 1 cover film (024-232)
- ➡ thermal coupling solution (0.1% Triton or 0.1% Tween 20, degas carefully)
- ➡ gel cover plate

6.3.2 Gel setup for electrophoresis

- 1) Adjust electrophoresis chamber with the 4 leveling feet.
- 2) Fill 500 ml **running buffer** in each buffer tank (check orientation of the buffer tanks: parallel or perpendicular TGGE! see section 6.4 and 6.5).

Note: Wipe off any spilled buffer from the electrophoresis unit. Never run device if buffer has been spilled.

- 3) Soak 4 **buffer wicks** with running buffer
- 4) Disassemble **gel sandwich**. Clean backside of the gel support film with a soft tissue.
- 5) Apply not more than 2 ml of **thermal coupling solution** (0.1% Triton or 0.1% Tween 20) on the thermoblock

Note: The volume of coupling solution should be as small as possible. Excess coupling solution leads to an irreproducible temperature distribution under the gel. The result is a wavelike migration front and poor separation of fragments.

- 6) Place the **gel on the thermoblock**. The thermal coupling solution should spread over the whole block. Avoid formation of bubbles. Wipe off any residual coupling solution along the edges of the gel support film.

Note: The thermal coupling solution is essential for efficient heat transfer from block to gel. If bubbles are entrapped under the gel support film, remove support film with gel from the block and place it back again on the block.

- 7) Cover the gel with a cover film. The cover film should be placed just beneath the slots.
- 8) Attach pre-soaked **buffer wicks** on top and bottom of the gel (fold 2 sandwiches of two wicks each in the middle; the folded side should face the gel)




Figure 5: Setup of gel , cover film and buffer wicks.

- 9) Load **samples** (approx. 5µl each for parallel gel with 32 slots, approx. 200µl for perpendicular gel with)

Note: Be careful not to touch the samples with the buffer wick! Otherwise the samples will diffuse into the wick.

- 10) Attach **gel cover plate** (the cover plate should have contact to the wicks, but must not squeeze gel or wicks).

- 11) Close **safety lid** and start **run**.

Note: For parallel TGGE let temperature gradient equilibrate for approx. 10 minutes, then start main run. This step may be omitted for a perpendicular gel.

6.4 Perpendicular TGGE

In perpendicular TGGE one sample is separated over a broad temperature range. This application is mainly used to check the melting behavior of a sample (see section 10.2). For casting of the gel use a glass plate with one large slot former (024-228). The temperature gradient must be orientated perpendicular to the migration of the sample. The buffer tanks must be positioned as described in Figure 6.

The migration of DNA / RNA molecules is indicated by the arrow on the safety lid of the electrophoresis unit.




Figure 6: Positioning of buffer tanks for perpendicular TGGE. Be sure that the direction of electrophoresis is perpendicular to the temperature gradient. The temperature gradient is indicated by the lines on the edge of the block.

6.5 Parallel TGGE

In parallel TGGE multiple samples are separated along the temperature gradient. For casting of the gel use glass plate with 32 slot formers (024-229). The buffer tanks for parallel TGGE must be positioned as depicted in Figure 7.





Figure 7: Positioning of buffer tanks for parallel TGGE. Be sure that the direction of electrophoresis is parallel to the temperature gradient. The slots of the gel should be at the same side as the markings on the block.

Note: For a parallel TGGE a 10 minutes equilibration of the temperature gradient may be included after the pre-run.

6.6 How to identify the optimum temperature range from a perpendicular gel

The theoretical background for the separation of DNA fragments in a perpendicular gel is described in section 10.2.

Place the stained perpendicular gel on the plastic film with the printed lines (L0 to L10). Identify the line where the double strand starts to melt (T1) and the line where the double strands separates into the single strands (T2).



The calculation of the corresponding temperatures is simple, since there is a linear temperature gradient between L0 and L10 (i.e. the temperature increment from one line to next line is always the same).

Calculation: Divide range of gradient by ten, this is the temperature increment from one line to the adjacent lane.

Example: calculation of temperature at line 6 (L6) in a temperature gradient from 40°C (L0) to 60°C (L10)

- **subtract** temperature at L0 from temperature L10 (range of gradient: $60-40^{\circ}\text{C} = 20^{\circ}\text{C}$)
- **divide** temperature by 10 (increment from line to lane: $20^{\circ}\text{C}/10 = 2^{\circ}\text{C}$)
- **multiply** increment by 6 (6 increments from L0 to L6: 12°C)
- **add** this value to the temperature at L0 ($40^{\circ}\text{C} + 12^{\circ}\text{C}$)

result: temperature at L6 is 52°C

7. Programming the TGGE controller

All parameters of the run are controlled by the TGGE system controller. This includes electrophoretic parameters (voltage, amperage, time) as well as control of the temperature gradient.

Note: Each program can consist of several steps. Thus you can program pre-run, equilibration and main run in the same program (see also section 6.2).

7.1 Create / edit program

Main screen

```
L0: 20.0°C   L10: 20.0°C
maxi block off
A?  BElpho  Cprograms  D+
```

Press **C** [programs] to enter the programming mode

```
program no:

Alist Bdel Cquit Denter
```

7.2 Select program

Enter a program directly by number or press **A** [list] to view a list of the existing programs.

```
0 test 1 maxi
1 parallel mini
2 empty
A↑ B↓ C quit D enter
```

Scroll through the list with **A**↑ **B**↓ and accept highlighted program with **D** [enter].

Note: Each program depends on the type of TGGE system (mini or maxi). The type (mini or maxi) for which the program was written, is displayed behind the program name. Programs can only be written /edited in the respective mode. Programs that are created in the maxi mode are automatically saved as maxi programs and can only be run with the maxi system (and vice versa).

To set type of system (mini or maxi) see section 7.9.4.

7.3 Name program

Each program is specified by name and a program number. To facilitate retrieval of a program, you can enter a name for each program existing of letters, numbers and symbols.

```
name:>          <
ABCDEF GHIJ KLMNOPQRST
UVWXYZ -( )αδ/,<>&+.%
A→ BABC C quit D enter
```

Press **B** [ABC] to enter the mode for the selection of letters.

```
name:>          <
ABCDEF GHIJ KLMNOPQRST
UVWXYZ -( )αδ/,<>&+.%
A← B→ C quit D enter
```

Move to the desired letter with **A** [←] and **B** [→]. Accept highlighted letter with **D** [enter].

```
name:>test      <
ABCDEF GHIJ KLMNOPQRST
UVWXYZ -( )αδ/,<>&+.%
A← B→ C quit D enter
```

If the program name is complete, confirm name with **D** [enter]. In the following screen you can set the temperatures for the gradient block.

7.4 Enter temperatures for the gradient block

```
1:L0:          L10:
A? B          C quit D →
```

Enter temperature for L0 and accept with **D** [enter].

```
1:L0:30.0      L10:
A? Bdelete Cquit Denter
```

Enter temperature for L10 and accept with **D** [enter].

```

1:L0: 30.0      L10: 70.0
ok?
A  B no  C quit D yes

```

Confirm settings with **D** [yes]. If settings are not correct press **B** [no] and repeat entry of temperature settings.

After you have confirmed the temperature settings, the following screen is displayed. Here you can enter all parameters for electrophoresis.

7.5 Enter electrophoresis parameters

```

1:L0: 30.0      L10: 70.0
time:
El:  0V  500mA  50W
A?  B V*h  C quit D →

```

Enter time for electrophoresis and accept with **D** [enter].

Note: There is a convention on how time settings are entered in all BIOMETRA instruments:

hours • minutes • seconds

If you enter a number without “dot” this value will be interpreted as seconds (“300” => 5 minutes). To program minutes enter a “•” after the number of minutes. To enter hours enter • after the number. You can also enter any combination of hours, minutes and seconds. Example: for 1 hour, 30 minutes, 20 seconds enter 1•30•20.

The time values will be displayed in the following format: 00 m 00s

Accept time setting with **D** [enter]

```

1:L0: 30.0      L10: 70.0
time: 10m 0s
El:  0V  500mA  50W
A?  Bdelete C quit D enter

```

Enter Voltage and accept with **D** [enter]

```

1:L0: 30.0      L10: 70.0
time: 10m 0s
El:  300V  500mA  50W
A?  Bdelete C quit D enter

```

Note: The values for amperage [mA] and wattage [W] are set to maximum level as default. If you enter lower values, these parameters may become limiting during electrophoresis.

Accept default settings for voltage and wattage with **D** [enter] or enter different values.

Note: We recommend to control electrophoresis by constant voltage rather than by constant amperage (set [mA] to maximum value, respectively accept default value).

Note: Each program can consist of several steps. Thus it is possible to program complex protocols including a pre-run, a pause for handling of the gel and the main run..

In the following screen you can program a second step for your protocol.

```

2:L0:                L10:

A?   B           C quit D →

```

If you do not want to program another step, accept program with **C** [quit].

The program name, number of steps and the total run time is displayed.

```

program no. 8
pgm end: 1 step(s)
run time: 0h10m 0s
A?   B           C quit D →

```

7.6 Start electrophoresis

Main screen:

```

L0: 20.0°C   L10: 20.0°C
maxi block off

A?   BElpho   Cprograms D+

```

To start a program press **B** [Elpho]

```

start program:

Alist Bdel Cquit Denter

```

Enter program number or choose a program from the list with **A** [list].

```

start program: 8

Alist Bdel Cquit Denter

```

Confirm program number with **D** [enter].

The program starts and parameters of gradient block and electrophoresis are displayed. During temperature equilibration of the gradient block the elapsed time is displayed.

Note: Electrophoresis starts as soon as the set temperature in the block is achieved.

```
L0: 30.0°C   L10: 70.0°C
hold: 1   2m12s  11.4Vh
El: 300V   8mA   20.3W
A?  BElpho  Cprograms  D+
```

7.7 Stop/pause electrophoresis

```
L0: 30.0°C   L10: 70.0°C
hold: 1   2m12s  11.4Vh
El: 300V   8mA   20.3W
A?  BElpho  Cprograms  D+
```

To stop/pause the active program press **B** [Elpho]

```
program 8 test
pause ?
stop ?
A?  Bpause  Cquit  Dstop
```

Press **B** [pause] to pause program

Press **D** [stop] to stop program

Press **C** [quit] to return to the active program.

7.8 View temperatures of the gradient

```
L0: 30.0°C   L10: 70.0°C
hold: 1   2m12s  11.4Vh
El: 300V   8mA   20.3W
A?  BElpho  Cprograms  D+
```

To display the temperatures in the block during a run press **A** [?]:

```
L0 : 30.1
L1 : 34.2

L2 : 38.3
A↑  B↓  C quit  D enter
```

You can scroll through the different lines with **A**↑ **B**↓.

```

L6 : 54.8
L7 : 58.9
L8 : 62.9
A↑ B↓ C quit D enter

```

Note: There is a difference in the number of lines depending on the type of TGGE system (mini or maxi) that is installed. In the maxi mode 10 different temperatures are displayed (L0 to L10) in the mini mode 5 different temperatures are displayed (L1 to L6).

7.9 Special functions

Main screen

```

L0: 20.0°C    L10: 20.0°C
maxi block off
A? BElpho Cprograms D

```

Press **D** [+] to enter the menu for special functions

```

1 print programs
2 signal
3 language
A↑ B↓ C quit D enter

```

Scroll through the list with **A**↑ **B**↓.

7.9.1 Print programs

Connect controller to a dot matrix printer. Select option 1 in the above menu and confirm with **D** [enter].

7.9.2 Select / de-select signal

Select option 2 in the special functions. Press **A** [on] to activate the signal, press **B** [off] to inactivate signal.

7.9.3 Select language

Select option 3 in the special function screen. Choose between German and English

7.9.4 Set block type

Select option 6 in the special function screen. Choose between mini and maxi system.

Note: The selection of the TGGE system type (mini or maxi) is saved together with each individual program. Programs that have been written in the MAXI mode can only be edited

and run in the maxi mode. Programs that have been written in the mini mode can only be edited and run in the mini mode.

8. Staining

8.1 Silver staining

Aside from autoradiography silver staining is the most sensitive method for detecting small amounts of DNA, RNA or proteins in polyacrylamid gels. Other staining protocols may be used, but generally exhibit less sensitivity. This must be considered in relation to the amount of DNA loaded on the gel.

All incubation steps are done in small plastic containers which are agitated on a rocking platform (e.g. order number 042-400 or 042-500).

Wear non-powdered protective gloves during all steps of the silver staining protocol to avoid staining artifacts due to the high sensitivity of the staining protocol.

The quality of chemicals is essential in silver staining. Prepare solutions freshly, use only chemicals of high quality (p.a.) and fresh double distilled water.

- **Important:** Remove the protective plastic sheets from the gel.
- Carefully remove any residual thermal coupling solution from the back of the gel (gel support film) prior to staining.
- Put the polyacrylamid gel with the gel side upwards into the staining tray. Avoid air bubbles during all staining steps.
- It's recommended to prepare at least 400 ml solution for each incubation step.
- Prepare stopping solution prior to developing.

8.1.1 Silver staining protocol

Step	Time	Solutions*
Fixation	30 min	300-400 ml 10% glacial acid, 30% EtOH
Sensitization	2 x 30 min	300-400ml 30% EtOH
Washing	Rins gel 30 seconds under running water, then wash 5 x 10 min	Fresh aqua dest
Silver Binding	30 min	400ml 0,1 % AgNO ₃ , prepare freshly add 400µl Formaldehyde (37%) prior to use

Washing	Rins 30 seconds then wash 1 min Rins again 30 seconds	Fresh aqua dest.
Developing	Until bands become visible, can take several minutes, don't let gel unattended!	Solution 1: dissolve 2g Sodium thiosulfate ($\text{Na}_2\text{S}_2\text{O}_3$) in 10ml bidest, Solution 2: dissolve 10g Sodium Carbonate (Na_2CO_3) in 400ml bidest Add 400µl solution 1 to solution 2 Add 400µl Formaldehyde (37%)
Stopping	30 min	Dissolve 5,84g EDTA and 8g Glycine in 400ml bidest
Storage	Up to several days	10% Glycerol

8.2 Ethidium bromide-staining

Incubate the gel in staining solution (0.5 µg/ml ethidium bromide in 1 X TBE) for 30 - 45 min. Analyze under UV radiation (27).

8.3 Autoradiography

TGGE gels can also be directly exposed to x-ray films if radiolabeled samples are analyzed.

Direct exposure:

Incubate the TGGE gel for 15 min. in Fixation solution (see 6.5 Silver staining). Optional: Silver stain the gel.

Remove residual buffer from the gel. Expose to an x-ray film at room temperature.

Exposure of dried TGGE gels:

Incubate the TGGE gel for 15 min. in Fixation solution (see 6.5 Silver staining). Optional: Silver stain the gel.

Incubate the gel in 2-5% glycerol for 10 minutes to prevent the gel from cracking. Incubate an appropriate sheet of cellophane (no Saran wrap!!!!) in 2 - 5% glycerol. Layer the cellophane on the gel. Air dry at room temperature for one day or use a gel dryer at 50°C for at least 3h. Exposure to an X-ray film.

8.4 Elution of DNA from the TGGE gel

DNA fragments which have been separated on TGGE, for example, different alleles of one gene, can be eluted from silver-stained TGGE gel and re-amplified by PCR.

Using a Pasteur pipette, puncture the gel and extract a μl piece containing the particular DNA duplex. Incubate in 20 μl TE buffer overnight. Use a 1 μl aliquot for re-amplification.

9. Sample preparation

9.1.1 Purity of samples

Due to the high sensitivity of the staining procedure after TGGE it is recommended to use purified DNA, RNA or protein samples. Any impurities might be misinterpreted after TGGE, thereby making the analysis of gels difficult. Nevertheless it is possible to use even crude mixtures for TGGE analysis.

PCR-amplified DNA fragments usually can be analyzed without further purification. Please note, that the presence of high amounts of nonspecific, secondary PCR products may result in difficulties with interpretation of band pattern, melting profile, etc. For example, in parallel TGGE, nonspecific bands with a higher molecular weight than the specific PCR product may be misinterpreted as heteroduplices, or analogs with lower thermal stabilities. Therefore, prior to TGGE check the PCR product in a conventional agarose gel. If necessary, purify your specific PCR product, e.g., by agarose gel electrophoresis and subsequent gel extraction.

9.1.2 Sample preparation for direct DNA analysis

1 volume of DNA/RNA samples is mixed with 1 volume of TBE or Na-TAE loading buffer or with 0.1 volume of the total loading volume ME loading buffer (see Appendix). The resulting mixture is loaded directly on to the polyacrylamid gels. Be sure that the slots are filled up to maximum (if necessary, add 1x loading buffer to fill up the slots to maximum).

In case of low-concentration samples we recommend to prepare 5x conc. loading buffer. 0.2 volume of this concentrated loading buffer is mixed with 0.8 volumes of the sample and loaded onto the gel.

9.1.3 Denaturation / Renaturation for heteroduplex analysis of DNA

Mix sample with equal amount of standard DNA and heat to **95°C for 5 minutes (denaturation)**. Then let slowly cool down to **50°C (renaturation)**. This can be done by programming a thermocycler to 94°C for 5 minutes and then 50 °C for 15 minutes with a ramping rate of -0.1°C/second. The sample is then loaded directly to the gel. In order to achieve the recommended loading volumes for diagonal or perpendicular TGGE , the sample volume should be adjusted with running buffer.

10. Optimization of TGGE

There are 3 steps in the setup of a new TGGE experiment:

- 1) Design of the PCR fragment
- 2) Identification of the correct temperature gradient
- 3) Parallel analysis of multiple samples

10.1 Design of DNA fragment for TGGE

The design of the DNA fragment is an important step for successful TGGE. Starting with the gene fragment of interest PCR primers should be designed with a conventional computer program. The melting behavior of the resulting fragment should then be checked with the Poland software. It is essential that the DNA fragment shows different melting domains. If there is only one single melting domain, an artificial higher melting domain (called GC clamp) must be added during PCR.

10.1.1 Poland analysis

The melting profile of a DNA fragment can be analyzed with a computer program. The Poland software calculates the melting behavior of a DNA fragment according to its base sequence. This software is free accessible via the internet.

(<http://www.biophys.uni-duesseldorf.de/POLAND/poland.html>).

How to perform a Poland analysis

- Open start page (URL see above)
- 1) enter a name for the query
- 2) copy / paste DNA sequence in the sequence window
- 3) choose the T_m plot (de-activate all other plots)
- 4) submit query
- retrieve T_m plot (melting curve)

Poland service request form

The Poland server will calculate the thermal denaturation profile of double-stranded RNA, DNA or RNA/DNA-hybrids based on sequence input and parameter settings in this form.

NEW: Thermodynamic parameters set for dsDNA in 75 mM NaCl (Blake & Delcourt) added.

Calculation is based on **Poland**'s algorithm in the implementation described by **Steger**. Graphics results are directly sent to your WWW client.


1)

Sequence title line:

2)

Sequence:

(plain format;
no numbers;
max. 1000 nts;
min. 5 nts)



Mismatched positions: (comma-separated numbers)

Thermodynamic parameters:

Dissociation constant β : Oligonucleotide (β is function of seq.length) Long double strand (default: $\beta=1.0E-3/M$)

☐ ☐

Strand concentration:
(default: 1.0E-6 M)

Temperature range: Low temperature limit: (default: 40.0 °C) High temperature limit: (default: 110.0 °C) Temperature step size: (default: 2.0 °C)

3) **Which graphics do you want:**

$T_m(p=50\%)$ plot ☒ 3d plot ☐ Mobility plot ☐ Melting curve ☐ Diff. melting curve ☐

Graphics size:
(GIF format)

4) Click here to , or Click here to the form to defaults.

The T_m plot (second order, red color) shows the melting profile of the DNA fragment according to the base sequence. The ideal fragment shows at least two distinct melting domains. Note that mutations can be detected in all but the highest melting domain. This means that in a DNA fragment with two melting domains, mutations can only be detected in the lower melting domain.




Figure 8: Tm plot of a 140bp DNA fragment resulting from Poland analysis. The second order curve (red color in the original) shows two different melting domains.

If the fragment consists of a single melting domain only, or if you want to scan the entire fragment for mutations, add a so called GC clamp to one end of the PCR fragment.

10.1.2 GC clamps

A GC clamp is an artificial, high melting domain which is attached to one end of the fragment during PCR. The name “GC clamp” implies that this short stretch will hold the DNA fragment together, preventing a dissociation into the single strands at higher temperatures. The optimum location for the GC clamp at the PCR fragment (5′ or 3′) can be easily checked with the Poland software. Copy / paste the GC sequence to either side of your sequence and repeat Poland analysis. In the following box you will find different examples for a GC clamp.

short GC-clamp (23 bp):	cccgc cgcgc cccgc cgcgc gcc
long GC-clamp (40 bp) ⁴⁴ :	cgccc gccgc gccgc gccgc cgccc cgccc cccc gccgc
long GC-clamp (39 bp) ⁴⁵ :	ccccg cccc gccgc cccc ccgcg cccgg cgccc ccgc

To integrate a GC clamp into a PCR fragment, one of the two primers has to be modified. The non-specific GC sequence is added to the 5′-end of the primer. Thus the GC sequence is incorporated in the fragment during PCR.

10.1.3 Chemical clamp with Psoralen (Furo[3,2-g]coumarin, C₁₁H₆O₃)

In addition to “clamping” a fragment with an artificial high melting domain it is as well possible to covalently fix the end of a PCR fragment. To achieve this, one of the primers carries a Psoralen molecule. Psoralen is a high reactive group when exposed to UV radiation. Thus it is possible to covalently close one end of the PCR fragment. The optimal primer sequence may be 5′(Pso)pTaPnpnp.....3′, given the preference of Psoralen for binding between TpA and ApT pairs^{13,46,47}. Crosslinking of the PCR product is done e.g. in a flat-bottom microtiter plate using a 365 nm UV source. Working with small volumes it may be necessary to minimize evaporation by cross-linking at 4 – 10°C. The yield is not affected by temperature. The distance of the sample from the UV source affects the yield. 15 min at 0.5 cm distance of the sample from an 8 W UV lamp is sufficient.

10.1.4 Use of SSCP primers

In many cases primer from SSCP may be used for TGGE analysis. Nevertheless, the resulting DNA fragments should be checked in the Poland analysis. If there is only one melting domain add a GC clamp to one of the primers (see section 10.1.2)

10.2 Find correct temperature gradient

Poland analysis gives the first indication, which temperature gradient should be applied for

parallel analysis of multiple samples. Under experimental conditions separation is performed in the presence of high concentrations of urea. Urea lowers the melting temperature of the DNA. This is important because gel electrophoresis at very high temperatures may lead to partial drying of the gel, resulting in a disturbed separation pattern. Therefore it is necessary to identify the optimum temperature gradient under experimental conditions.

To identify the optimum temperature gradient the DNA fragment is separated in a perpendicular TGGE. This means the temperature gradient is perpendicular to the migration of samples (see section 6.4). Thus the migration of a fragment can be checked simultaneously at different temperatures in a single run. If the PCR fragment has been designed properly the separation in a perpendicular temperature gradient leads to a distinct melting curve (see Figure 9)




Figure 9: Identification of the optimum temperature gradient in a perpendicular TGGE. At low temperature (below T1) DNA migrates as a double strand (left side). At intermediate temperature (between T1 and T2) the DNA opens at one side, the partial double strand is increasingly slowed down. Above T2 the DNA separates into the single strands.

At T1 the double strand starts to melt and forms a branched structure. At T2 the partial double strand separates irreversibly into the single strands. Analysis of samples in parallel TGGE should be performed precisely in this temperature range between T1 and T2.

How to identify the optimum temperature range from a perpendicular gel:

Place the stained gel on the plastic film with the printed lines (I0 to L10). Identify the line where the double strand starts to melt (T1) and the line where the double strands separate into the single strands (T2). For the calculation of temperature at the corresponding lines see section 6.6.

10.3 Parallel analysis of multiple samples

After identification of T1 and T2 in a perpendicular TGGE this temperature gradient is spread over the whole block for parallel analysis.




Figure 10: Application of T1 and T2 in a parallel gel.

Note: The DNA fragments are separated by their melting behavior. They can be distinguished as soon as the fragments begin to melt, i.e. they form a fork like structure (temperature higher than T1). During electrophoresis the fragments should not separate into single strands. This is an irreversible transition resulting in diffuse bands.

Note: If there are only small differences in the migration of different samples, perform a heteroduplex analysis (see chapter 10.5)

10.4 Optimization of parallel TGGE

To improve separation in parallel TGGE the gradient should start directly at the temperature where the fragments start to melt (see perpendicular gel) and should be rather flat. This means there should be only a moderate temperature increase over the whole gel. Different fragments in one sample separate as soon as the first fragment starts to melt. At a certain (higher) temperature the next fragment starts to melt. In a moderate gradient, the temperature increase per centimeter is smaller than in a steeper gradient. This means, the distance between two temperatures (i.e. locations in the gel) is bigger than in a steeper temperature gradient. This results in a wider separation of fragments that melt at different temperatures (see Figure 11).




Figure 11: Parallel TGGE using a steep (A) or a flat (B) temperature gradient. With a smaller temperature gradient (30 to 40 °C, B) the separation of samples is much wider. (Note: the temperatures in this figure are only for demonstration.)

10.5 What to do, if different samples have very similar melting points: Heteroduplex analysis with TGGE

10.5.1 Principle of heteroduplex analysis

If the difference in melting temperature between wildtype and mutant is very small, heteroduplex analysis is a rewarding approach. Heteroduplex analysis makes it very easy to distinguish between the wildtype and mutant form of a DNA fragment. The basic principle is to mix each sample with an external standard. In most cases this standard is a PCR fragment without mutations, for example amplified from the wild type. After mixing the standard DNA fragment with the PCR fragment from the sample the mixture is heated and subsequently slowly cooled down (for protocol see section 9.1.3).




Figure 12: Principle of heteroduplex analysis.

The re-annealing of sample and standard results in 4 different DNA fragments. 1) The wildtype homoduplex (AA), 2) the mutant homoduplex (aa), 3) and 4) two different heteroduplexes (Aa and aA). These heteroduplexes carry at least one mismatch (disturbed base pairing) and have a significant lower melting temperature than the homoduplexes.

This procedure results in a complete denaturing of both double stranded PCR fragments and a subsequent re-annealing. If the sample is different from the standard, re-annealing leads to 4 different double stranded DNA fragments (see Figure 12): 1) the homoduplex of the standard (wildtype AA), 2) the homoduplex of the sample (mutant aa), 3) and 4) two heteroduplexes between standard and sample (Aa and aA). Due to the differences between sample and standard these heteroduplexes display mismatches in their base pairing in least one position. Such mismatches have a strong impact on the melting behavior because the number of base pairs between the two strands is reduced. Therefore the heteroduplexes can be easily separated from the homoduplexes using TGGE.

The identification of the optimum temperature gradient for the separation of a heteroduplex analysis is absolutely the same as for a single fragment. The separation of a heteroduplex sample in a perpendicular TGGE results in 4 different melting curves. The 2 heteroduplexes have a lower melting temperature and denature at a lower temperature compared to the homoduplexes.




Figure 13: Separation of a heteroduplex sample in **perpendicular TGGE**.

The temperature gradient can then be adapted in the same way as for a conventional sample (see chapter 10.2). In parallel TGGE, the samples melt as they migrate along the temperature gradient. The heteroduplexes (with mismatch) melt at a lower temperature than the homoduplexes. Thus they open earlier in the partial single strand and are slowed down in the

gel matrix. The homoduplexes migrate a longer distance as complete double strands and start to melt at a higher temperature (i.e. later in respect to the temperature gradient). Therefore the lower bands in parallel TGGE are the homoduplexes, whereas the higher bands are the heteroduplexes.




Figure 14: Schematic drawing of a screening multiple samples in a **parallel TGGE**. Both homoduplexes (AA, aa) have a higher melting temperature and migrate further in the gel. The heteroduplexes melt at a lower temperature resulting in a slower migration.

10.5.2 Evaluation of a heteroduplex analysis

There are two possible states in heteroduplex analysis: 1) the sample is identical to the standard (wildtype) 2) the sample is different from the wildtype. In the former case, the denaturation / renaturation procedure results in one (the same) homoduplex. The subsequent separation in parallel TGGE shows only a single band. In the latter case, denaturation / renaturation leads to the four different populations depicted in Figure 12. Separation in parallel TGGE results in up to four different bands (see

Figure 13). If the temperature gradient has not been correctly optimized, or if separation time was too short, there may as well only be two or three bands.

This makes heteroduplex analysis very easy to evaluate:

number of bands	result
one	sample is identical to the standard no mutation
more than one (up to 4 bands)	sample is different from the standard mutation

11. The TGGE test kit (order number 024-050)

The TGGE test kit was developed to get familiar with the TGGE system . It consists of 3 different DNA samples:

- a wild type sample (DNA fragment without mutation)
- a mutant sample (DNA fragment that differs in on position from the wild typ)
- the heteroduplex sample (sample has been prepared as described in section 10.5.1)

The samples are separated in a **8% PAA** gel with **8M Urea** and a **1 x TAE** buffer system (for preparation of gel solution and buffer see section 5.2.1).

11.1 Perpendicular TGGE using the Biometra TGGE test kit

1) sample preparation

mix **100 µl** heteroduplex sample with
100 µl loading buffer TAE (see section 16.1)

- 2) Load heteroduplex sample to the broad slot of a perpendicular gel
- 3) Let sample migrate into the gel with 350V for approx. 20 minutes
- 4) Cover gel with cover film, assemble buffer wicks, cover plate and safety lid
- 5) start run

Temperature gradient	30 to 70°C
Voltage:	300V
Run time:	4h

6) silver stain gel

11.2 Parallel TGGE using the Biometra TGGE test kit

- 1) sample preparation:
mix **5µl sample** (wildtype or mutant or heteroduplex) with
10µl running buffer (TAE) and add
15µl loading buffer TAE (see section 16.1)
- 2) Assemble electrophoresis unit, cover gel with cover film (beneath the slots), assemble cover plate and safety lid
- 3) Load 5 µl of wildtype, mutant and heteroduplex samples
- 6) Let sample migrate into the gel with 400V for approx.. 5 minutes
- 4) After pre-run, cover gel like described in section 6.3.2
- 5) start main run

Temperature gradient	30 to 60 °C
Voltage:	400V
Run time:	3 h

- 4) silver stain gel

12. Technical specification

TGGE MAXI System	Electrophoresis unit with temperature gradient block and two removable buffer chambers, Controller, Power supply, Starter kit.
TGGE MAXI Electrophoresis unit	
Temperature gradient formation	High performance Peltier technology
Block size	20 x 20 cm
Temperature range	5 – 80 °C
Linear temperature gradient	maximum 45 °C
Temperature accuracy	± 0.3 °C
Temperature uniformity	± 0.5 °C
Glass plate size	23.5 x 23.5 cm
Gel size	approx. 20 x 20 cm
Separation distance	parallel: 16 cm perpendicular: 19 cm
Dimensions (L x W x H)	42.3 x 42.3 x 33.3 cm
Weight	22 kg
TGGE MAXI System Controller	Microprocessor driven control of temperature gradient and electrophoretic parameters
Current	maximum 500 mA
Voltage	maximum 400V
Wattage	maximum 50W
Program stores	up to 100 programs can be stored
Program modes	constant Voltage V/h integration Programs can contain different steps (pre-run, pause, run)
Display	LCD display, 4 lines, English / German
Interfaces	parallel (Centronics), serial (RS 232)
Dimensions	31 x 22 x 11.5 cm
Weight	3.5 kg
Voltage	110 / 230 V
TGGE MAXI Power Supply	
Dimensions	29.5 x 22 x 8 cm
Weight	6.5 kg
Transformer	450 VA

13. Ordering information:

TGGE MAXI System: Electrophoresis unit with gradient block, controller, power supply, manual, starter kit	024-200
TGGE MAXI starter kit 1 glass plate perpendicular, 1 glass plate parallel, 2 glass plates without Spacer, 12 clamps, 2 silicone sealings, polybond film (25), gel cover film (10), buffer wicks (100), AcrylGlide (100ml)	024-204
TGGE test kit Test DNA for perpendicular and parallel test runs: wild type DNA (control), mutant DNA, heteroduplex DNA, sample buffer.	024-050

Self adhesive slot forming units (8 strips with 28 units á 8µl each for parallel gels, 9 strips with one broad slot á 200µl each for perpendicular gels)	024-121
TGGE MAXI buffer wicks, 100/pkg (18 x 20 cm)	024-215
TGGE MAXI glass plate without Spacer 23.5 x 23.5 cm	024-221
Self adhesive slot forming units for replacement, 28 pcs. á 8µl	024-222
Applicator strips 240mm (3 pcs.), 43 slots á 8µl each	024-223
TGGE Maxi glass plate with spacers and no slot formers for use with applicator strips (024-223)	024-227
TGGE MAXI glass plate perpendicular, with Spacer (1mm) and slot former (1 slot, 75 µl)	024-228
TGGE MAXI glass plate parallel, with Spacer (1mm) and slot former (32 slots, 5 µl)	024-229
TGGE MAXI silicon sealing for casting of gels, 1mm	024-230
TGGE MAXI gel cover film 25/pkg	024-232
TGGE MAXI polybond film 25/pkg	024-234
TGGE MAXI polybond film 100/pkg	024-235

14. Trouble-shooting

The following trouble-shooting guide may be helpful in solving any problem that you may encounter. If you need further assistance, please do not hesitate to contact your local Biometra distributor or Biometra.

Problem	Cause	Solution
Electrophoresis		
Leakage of gel cuvette	Inaccurate positioning of sealing Dust on spacers Acryl glide on spacers	Check positioning of silicone sealing, clean Spacer, do not apply Acryl Glide onto the spacers
Slotformer fall off	After heavy use this may happen from time to time	Self adhesive slot forming units (024-221) are included in the TGGE maxi system
Acrylamide solution gets behind the support film during pouring the gel	Support film is not properly attached to the glass plate	Fix polybond film with adhesive tape along the upper edge of the glass plate
Teflon film peels away from the thermoblock	Block has been cleaned with strong detergents or aggressive chemicals	Contact Biometra
No current	Amperage and Wattage have been set to "0"	Set Amperage and Wattage to maximum values (Electrophoresis should be controlled by Voltage)
Current oscillates	Coating of the thermoblock is damaged, the safety shutoff is activated	Contact Biometra
Gel interpretation		
Wavelike migration front	Temperature inhomogeneity under the gel due to excess thermal coupling solution	Use as little as possible thermal coupling solution (not more than 2ml)
No sigmoid melting curve (perpendicular TGGE)	Fragment melts completely	Perform Poland analysis Optimize primer design
No separation of hetero duplex samples (parallel TGGE)	Inappropriate fragment Wrong temperature gradient Acrylamide of poor quality	Perform Poland analysis Perform perpendicular TGGE Use only high quality chemicals (p.a.)
Irreproducible gels	Erratic temperature distribution over and under the gel	Use only minimum volume of thermal coupling solution under the gel

		Do not overlay gel with buffer
Silver staining		
Bad silver stain	Chemicals of poor quality Stale water Too much silver nitrate	Use only high quality chemicals Use only freshly prepared aqua bidest Refer to the staining protocol
Strong background	Insufficient washing after incubation in staining solution	Extend wash step, change water frequently
Weak staining of bands	Excessive washing after binding of staining solution	Reduce wash step after staining

15. References

1. Riesner, D., Henco, K. and Steger, G. (1990): Temperature-Gradient Gel Electrophoresis: A method for the analysis of conformational transitions and mutations in nucleic acids and protein. Page 169-250 In Chrambach, A., Dunn, M.J., Radola, B.J.: *Advances in Electrophoresis*, Vol. 4, VCH Verlagsgesellschaft Weinheim
2. Kappes, S. et al.(1995): p53 mutations in ovarian tumors, detected by temperature-gradient gel electrophoresis, direct sequencing and immunohistochemistry. *Int. J. Cancer* 64: 52-59
3. Milde-Langosch, K. et al. (1995): Presence and persistence of HPV and p53 mutation in cancer of the cervix uteri and the vulva. *Int. J. Cancer* 63: 639-645
4. Horn, D. et al.(1996): Three novel mutations of the NF1 gene detected by temperature gradient gel electrophoresis of exons 5 and 8. *Electrophoresis* 17: 1559-1563
5. Wieland, U. et al.(1996): Quantification of HIV-1 proviral DNA and analysis of genomic diversity by polymerase chain reaction and temperature gradient gel electrophoresis. *J. Virology Methods* 57: 127-139
6. Kuhn, J.E. et al. (1995): Quantitation of human cytomegalovirus genomes in the brain of AIDS patients. *Journal of Medical Virology* 47: 70-82
7. Linke, B. et. al. (1995): Identification and structural analysis of rearranged immunoglobulin heavy chain genes in lymphomas and leukemia. *Leukemia* 9: 840-847.
8. Menke M.A. et al. (1995): Temperature gradient gel electrophoresis for analysis of a polymerase chain reaction-based diagnostic clonality assay in the early stages of cutaneous T-cell lymphomas.
9. Hecker, R. et al. (1988): Analysis of RNA structure by temperature-gradient gel electrophoresis: viroid replication and processing. *Gene* 72: 59-74
10. Baumstark, T. and Riesner, D. (1995): Only one of four possible secondary structures of the central conserved region of potato spindle tuber viroid is a substrate for processing in a potato nuclear extract. *Nucleic Acids Research* 23: 4246-4254
11. Loss, P., Schmitz, M., Steger, G. and Riesner, D. (1991): Formation of a thermodynamically metastable structure containing hairpin II is critical for the potato spindle tuber viroid. *EMBO Journal* 10: 719-728
12. Riesner, D. (1998): Nucleic acid structures. In: *Antisense Technology. Practical Approach Series*. Oxford University Press. p1-24 (in press)
13. Wiese, U. et al. (1995): Scanning for mutations in the human prion protein open reading frame by temporal temperature gradient gel electrophoresis. *Electrophoresis* 16: 1851-1860
14. Nubel, U. et al. (1996): Sequence heterogeneities of genes encoding 16S rRNAs in *Paenibacillus polymyxa* detected by temperature gradient gel electrophoresis.
15. Lessa, E.P. and Applebaum, G. (1993): Screening techniques for detecting allelic variation in DNA sequences. *Molecular Ecology* 2: 119-129
16. Richter, A., Plobner, L., Schumacher, J. 1997: Quantitatives PCR-Verfahren zur Bestimmung der Plasmidkopienzahl in rekombinanten Expressionssystemen. *BIOforum* 20: 545-547
17. Henco, K. and Heibey, M. (1990): Quantitative PCR – the determination of template copy numbers by temperature gradient gel electrophoresis. *Nucleic Acids Research* 18: 6733-6734
18. Birmes, A. et al. (1990): Analysis of the conformational transition of proteins by temperature-gradient gel electrophoresis. *Electrophoresis* 11: 795-801
19. Arakawa, T. et al. (1993): Analysis of the heat-induced denaturation of proteins using temperature gradient gel electrophoresis. *Analytical Biochemistry* 208: 255-259
20. Chen, X. et al. (1995): High resolution SSCP by optimization of the temperature by transverse TGGE. *Nucleic Acids Research* 23: 4524-4525
21. Scholz, R.B. et al. (1993): Rapid screening for Tp53 mutations by temperature gradient gel electrophoresis: a comparison with SSCP analysis. *Human Molecular Genetics* 2: 2155-2158

22. Elphinstone and Baverstock, P.R. (1997): Detecting mitochondrial genotypes by temperature gradient gel electrophoresis and heteroduplex analysis. *BioTechniques* 23: 982-986
23. Poland, D. (1974): Recursion relation generation of probability profiles for sequence-specific macromolecules with long-range correlations. *Biopolymers* 13:1859-1871
24. Lerman, L.S. and Silverstein, K. (1987): Computational simulation of DNA-melting and its application to denaturing gradient-gel electrophoresis. *Meth. Enzymol.* 155: 482-501
25. Steger, G. (1994): Thermal denaturation of double-stranded nucleic acids: prediction of temperatures critical for gradient electrophoresis and polymerase chain reaction.
26. Schumacher, J, Randels, J.W. and Riesner, D. (1983): A two dimensional electrophoretic technique for detection of circular viroids and virusoids. *Anal. Biochem.* 135, 288 - 295
27. Sambrook, J., Fritsch, E.F. and Maniatis, T. (1989): Molecular cloning, Cold Spring Harbor Laboratory press
28. Steger, G. and Riesner, D. (1992): Temperaturgradienten-Gelelektrophorese: eine Methode zur Analyse von Konformationsübergängen und Mutationen in Nukleinsäuren und Proteinen. In Radola, B.J. (ed) *Handbuch der elektrophorese*, VCH Verlagsgesellschaft, Weinheim
29. Sheffield, V.C., Cox, D.R. and Lerman, R.M. (1989): Attachment of a 40-base-pair G+C-rich sequence (GC-clamp) to genomic DNA fragments by the polymerase chain reaction results in improved detection of single-base changes. *Proc. Natl. Acad. Sci. USA* 86, 232 - 236
30. Hecker R., Wang Z., Steger G. and Riesner D. (1988): Analysis of RNA structure by temperature-gradient gel electrophoresis: viroid replication and processing. *Gene* 72, 59-74
31. Jiang L., Chen W., Tain L.P. and Liu Y. (1991): Temperature-gradient gel electrophoresis of apple scar skin viroid. *Acta Microbiol. Sin.* 30, 278-283
32. Riesner D., Hecker R. and Steger G. (1988): Structure of viroid replication intermediates as studied by thermodynamics and temperature-gradient gel electrophoresis. In Sarma R.H. and Sarma M.H. (eds.) *Structure & Expression, Vol. I: From Proteins to Ribosomes*, Adenine press, 261-285
33. Riesner D., Steger G., Zimmat R., Owens R.A., Wagenhöfer M., Hillen W., Vollbach S. and Henco K. (1989): Temperature-gradient gel electrophoresis of nucleic acids: Analysis of conformational transitions, sequence variations, and protein-nucleic acid interactions. *Electrophoresis* 10, 377-389
34. Rosenbaum V. and Riesner D. (1987): Temperature-gradient gel electrophoresis: thermodynamic analysis of nucleic acids and proteins in purified form and in cellular extract. *Biohys. Chem.* 26, 235-246
35. Schönborn J., Oberstraß J., Breyel E., Tittgen J., Schumacher J., Lukacs N. (1991): Monoclonal antibodies to double-stranded RNA as probes of RNA structure in crude nucleic acid extracts. *Nucleic Acids Res.* 19, 2993-3000
36. Po Tien, Steger G., Rosenbaum V., Kaper J. and Riesner D. (1987): Double-stranded cucumovirus associated RNA5: experimental analysis of necrogenic and non-necrogenic variants by temperature-gradient gel electrophoresis. *Nucleic Acids Res.* 15, 5069-5083
37. Zimmat R., Gruner R., Hecker R., Steger G. and Riesner D. (1991): Analysis of mutations in viroid RNA by non-denaturing and temperature-gradient gel electrophoresis. In R.H. Sarma and M.H. Sarma (eds.) *Structure & Methods, Vol. 3: DNA & RNA*, Adenine Press, 339-357
38. Rosenbaum V., Klahn T., Lundberg Holmgren E., von Gabain A. and Riesner D. (1992): Co-existing structures of an mRNA stability determinant: The 5' region of the Escherichia coli and Serratia marcescens ompA mRNA, *J.Mol.Biol.*, in press
39. Birmes A., Sättler A., Maurer S.O. and Riesner D. (1990): „Analysis of the conformational transitions of proteins by temperature-gradient gel electrophoresis“. *Electrophoresis* 11, 795-801
40. Sättler A., Kanka S., Schrörs W. and Riesner D. (1992): „Random mutagenesis of the weak calcium binding site in Subtilisin Carlsberg and screening for thermal stability by temperature-gradient gel electrophoresis“. Accepted for: 1st International Symposium of Subtilisin Enzymes, EMBL, Hamburg
41. Thatcher D. and Hodson B. (1981): „Denaturation of proteins and nucleic acids by thermal-gradient electrophoresis“. *Biochem. J.* 197, 105-109

42. Wagenhöfer M., Hansen D. and Hillen W. (1988): „Thermal denaturation of engineered tet repressor proteins and their complexes with tet operator and tetracycline studie by temperature-gradient gel electrophoresis“. Analytical Biochem. 175, 422-432
43. Sanguinetti C.J., Neto E.D. and Simpson A.J.G. (1994): BioTechniques 17, 915
44. Kappes, S., Milde-Langosch, K., Kressin, P., Passlack, B., Dockhorn-Dwornczak, B., Röhlke, P. and Löning, T. (1995): "p53 Mutations in ovarian tumors, detected by temperature-gradient gel electrophoresis, direct sequencing and immunohistochemistry". Int. J. Cancer 64, 52 - 59
45. Kluwe, L., MacCollin, M., Tatagiba, M., Thomas, S., Hazim, W., Haase, W. and Mautner, V.-F. (1998): "Phenotypic variability associated with 14 splice-site mutations in the NF2 gene". American Journal of Medical Genetics 77, 228 - 233
46. Lerman, L. S. and Beldjord, C. (1998). "Comprehensive mutation detection with denaturing gradient gel electrophoresis". In R.G.H. Cotton, E. Edkins and S. Forrest (eds.) Mutation Detection, A Practical Approach, Oxford University Press, 35 - 62
47. Gamper, H., Piette, J. and Hearst, J.E. (1984): Photochem. Photobiol. 40, 29 ff

16. Appendix

16.1 Buffers

16.1.1 Running buffers:

TBE Running buffer 0.1 x conc. TBE (up to 1x conc. TBE is possible)

10 x TBE (stock solution)	890 mM Boric Acid
	20 mM EDTA
	890 mM TRIS
	Do not titrate to adjust pH!

TAE Running Buffer 1 x conc. TAE, pH 8.0

50 x TAE (stock solution)	242g Tris base (2M)
pH 8.0	57,1 ml glacial acid
	100ml 0.5M EDTA (pH 8.0)

MOPS-Running Buffer 1 x conc. MOPS

50x MOPS (stock	1M MOPS
	50 mM EDTA
	pH = 8.0

16.1.2 Loading buffers:

Loading buffer TBE	TBE running buffer
	0.1% Triton-X 100
	0.01% Bromophenol Blue dye
	0.01% Xylene Cyanol dye

Loading buffer TAE	TAE running buffer
	0.1% Triton-X 100
	0.01% Bromophenol Blue dye
	0.01% Xylene Cyanol dye
	2 mM EDTA

Loading buffer MOPS	MOPS running buffer
	1 mM EDTA
	0.05% Bromophenol Blue dye
	0.05% Xylene Cyanol dye
	pH = 8.0

16.1.3 Other buffers:

TE buffer	10 mM Tris/HCl 0.1 mM EDTA pH = 8.0
TEMED	Solution of N,N,N',N'tetramethylethylenediamine
APS	10% Ammonium persulfate
Glycerol 40%	40% glycerol in water
Glycerol 50%	50% glycerol in water

17. Instructions for return shipment



If you would like to send the unit back to us, please read the following return instructions.

Should you have any problems with the TGGE System, please contact your local **Biometra** dealer or our service department:

Biometra biomedizinische Analytik GmbH

Service Department

Rudolf-Wissell-Straße 30

D-37079 Göttingen

Phone: ++49 – (0)5 51 / 50 68 6-0

Fax: ++49 – (0)5 51 / 50 68 6-66

Return only defective devices. For technical problems which are not definitively recognisable as device faults please contact the Technical Service Department at **Biometra**.

Use the original box or a similarly sturdy one.

Label the outside of the box with “CAUTION! SENSITIVE INSTRUMENT!”

Please enclose a **precise description of the fault**, which also reveals during which procedures the fault occurred, if possible.



- **Important:** Clean all parts of the instrument from residues, and of biologically dangerous, chemical and radioactive contaminants. Please **include a written confirmation** (use the “**Equipment Decontamination Declaration**” following on the next page) that the device is free of biologically dangerous and chemical or radioactive contaminants in each shipment. If the device is contaminated, it is possible that **Biometra** will be forced to refuse to accept the device.
- The **sender of the repair order will be held liable** for possible losses resulting from insufficient decontamination of the device.
- Please enclose a note which contains the following:
 - a) Sender’s name and address,
 - b) Name of a contact person for further inquiries with telephone number.

18. Equipment Decontamination Certificate

To enable us to comply with german law (i.e. §28 StrlSchV, §17 GefStoffV and §19 ChemG) and to avoid exposure to hazardous materials during handling or repair, will you please complete this form, prior to the equipment leaving your laboratory

COMPANY / INSTITUTE ---

ADDRESS

TEL NO _____

FAX NO _____

E-MAIL _____

EQUIPMENT

Model

Serial No

If on loan / evaluation

Start Date: _____ Finish Date _____

Hazardous materials used with this equipment_

Has the equipment been cleaned and decontaminated? YES / NO (delete)

Method of cleaning / decontamination:

NAME _____

POSITION _____

(HEAD OF DIV./ DEP./ INSTITUTE / COMPANY)

SIGNED _____

DATE _____

PLEASE RETURN THIS FORM TO BIOMETRA GMBH OR YOUR LOCAL BIOMETRA DISTRIBUTOR TOGETHER WITH THE EQUIPMENT.

**PLEASE ATTACH THIS CERTIFICATE OUTSIDE THE PACKAGING.
INSTRUMENTS WITHOUT THIS CERTIFICATE ATTACHED WILL BE
RETURNED TO SENDER.**

19. Warranty

This Biometra instrument has been carefully built, inspected and quality controlled before dispatch. Hereby Biometra warrants that this instrument conforms to the specifications given in this manual. This warranty covers defects in materials or workmanship for 12 month as described under the following conditions:

This warranty is valid for **12 month** from date of shipment to the customer from Biometra or an authorized distributor. This warranty will not be extended to a third party without a written agreement of Biometra.

This warranty covers only the instrument and all original accessories delivered with the instrument. This warranty is valid only if the instrument is operated as described in the manual.

Biometra will repair or replace each part which is returned and found to be defective.

This warranty does not apply to wear from normal use, failure to follow operating instructions, negligence or to parts altered or abused.