



Circulating Baths with Advanced Digital Temperature Controller

INSTRUCTION MANUAL



North American Catalog Number

| | <u>120 V</u> | <u>240 V</u> |
|---------------------------------------|--------------|--------------|
| AD07R-20 Refrigerating / Heating Bath | 89202-970 | 89202-972 |
| AD07R-40 Refrigerating / Heating Bath | 89202-978 | 89202-980 |
| AD7LR-20 Refrigerating / Heating Bath | 89202-962 | 89202-964 |
| AD15R-30 Refrigerating / Heating Bath | 89202-986 | 89202-988 |
| AD15R-40 Refrigerating / Heating Bath | 89202-994 | 89202-996 |
| AD20R-30 Refrigerating / Heating Bath | 89203-002 | 89203-004 |
| AD28R-30 Refrigerating / Heating Bath | 89203-010 | 89203-012 |
| AD45R-20 Refrigerating / Heating Bath | 89203-018 | 89203-020 |
| AD07H200 Heating Only Bath | 89202-926 | 89202-928 |
| AD15H200 Heating Only Bath | 89202-934 | 89202-936 |
| AD20H200 Heating Only Bath | 89202-942 | 89202-944 |
| AD28H200 Heating Only Bath | 89202-950 | 89202-952 |

European Catalog Number

| | <u>240 V</u> |
|---------------------------------------|--------------|
| AD07R-20 Refrigerating / Heating Bath | 462-0226 |
| AD07R-40 Refrigerating / Heating Bath | 462-0228 |
| AD7LR-20 Refrigerating / Heating Bath | 462-0224 |
| AD15R-30 Refrigerating / Heating Bath | 462-0230 |
| AD15R-40 Refrigerating / Heating Bath | 462-0232 |
| AD20R-30 Refrigerating / Heating Bath | 462-0234 |
| AD28R-30 Refrigerating / Heating Bath | 462-0236 |
| AD45R-20 Refrigerating / Heating Bath | 462-0238 |
| AD07H200 Heating Only Bath | 462-0215 |
| AD15H200 Heating Only Bath | 462-0217 |
| AD20H200 Heating Only Bath | 462-0219 |
| AD28H200 Heating Only Bath | 462-0221 |

Version: 1
Issued: 08 March 2013



Legal Address of Manufacturer

United States

VWR International, LLC
Radnor Corporate Center, Building One, Suite 200
P.O. Box 6660, 100 Matsonford Road
Radnor, PA 19087 U.S.A.
800-932-5000
<http://www.vwr.com>

Europe

VWR International bvba
Researchpark Haasrode 2020
Geldenaaksebaan 464
B-3001 Leuven
+ 32 16 385011
<http://be.vwr.com>

Country of Origin

United States of America

Table of Contents

| | |
|--|-----------|
| Introduction | 3 |
| VWR Signature Circulating Baths with the Advanced Digital Temperature Controller | 3 |
| General Safety Information | 4 |
| Safety Recommendations | 5 |
| Regulatory Compliance & Testing | 6 |
| Unpacking Your Circulator | 7 |
| Contents | 7 |
| Controls & Components | 8 |
| Advanced Digital Controller | 8 |
| Refrigerating/Heating Baths | 9 |
| Heating Only Baths | 10 |
| Quick-Start | 11 |
| Installation & Startup | 13 |
| General Site Requirements | 13 |
| Adding Liquid to the Bath Reservoir | 13 |
| Pump Inlet and Outlet Connections | 14 |
| External Closed Loop Circulation | 15 |
| Open Loop Circulation | 15 |
| Refrigeration Control Connections (Refrigerating/Heating Circulators only) | 16 |
| Electrical Power | 16 |
| Refrigerating/Heating Circulators | 16 |
| Heat Only Circulators | 17 |
| Communication | 18 |
| USB Communication | 19 |
| Ethernet | 19 |
| RS232 / RS485 Serial Communication | 19 |
| External (P2) Temperature Probe | 19 |
| Controller Setup | 20 |
| Power | 20 |
| Safety Set Temperature | 21 |
| Normal Operation | 22 |
| Keys and Controls | 22 |
| Turning Your Circulator ON | 22 |
| Main Operational Display | 23 |
| Set-Up Sub-Menus | 25 |
| Adjusting the Temperature Set Point | 26 |
| Selecting the Temperature Unit | 27 |
| Selecting the Pump Speed | 28 |
| Setting the High Limit Temperature | 29 |
| Setting the Low Limit Temperature | 30 |
| Setting the Differential Temperature | 31 |
| Setting the Auto Cool Temperature | 32 |
| Setting the Specific Heat Capacity (SHC) | 33 |
| Selecting Internal or External Control | 34 |
| Calibrating Your Circulator (Offset) | 35 |
| Displaying the Safety Set Temperature | 36 |
| Selecting a Remote Communication and Control Protocol | 37 |
| Using the Timer | 38 |
| Entering a Password | 39 |
| Enabling / Disabling Data Logging | 40 |
| Selecting the Operational Language | 41 |
| Setting Auto Restart | 42 |
| Resetting the Factory Default Values | 42 |
| Changing Your Circulator's Viewing Angle | 43 |
| Inert Gas Purge | 43 |
| Tap Water Cooling | 43 |
| Reservoir Cover Storage | 43 |

| | |
|---|-----------|
| Display Messages and Alarms | 44 |
| Routine Maintenance & Troubleshooting | 45 |
| Maintaining Clear Bath Water..... | 45 |
| Draining the Bath Reservoir..... | 45 |
| Checking the Over-Temperature / Low Liquid Level Safety Systems..... | 46 |
| Over-Temperature Protection..... | 46 |
| Low Liquid-Level Protection | 46 |
| Cleaning Your Circulator..... | 47 |
| Temperature Controller | 47 |
| Bath Reservoir..... | 47 |
| Pump Impeller..... | 47 |
| Condenser, Air Vents, and Reusable Filter (Refrigerating / Heating Circulators only)..... | 47 |
| Temperature Controller Removal and Re-Installation..... | 48 |
| Removal..... | 48 |
| Re-Installation..... | 49 |
| Viewing Component Operating Time..... | 50 |
| Troubleshooting Chart..... | 51 |
| Technical Information | 53 |
| Performance Specifications..... | 53 |
| Reservoir Fluids..... | 54 |
| Application Notes..... | 55 |
| Tubing and Fitting Temperature Ranges..... | 55 |
| RS232/RS485 Configuration..... | 56 |
| RS232/RS485 Communications..... | 58 |
| USB Data Logging..... | 60 |
| USB B Setup, Monitoring, and Control..... | 60 |
| Configuring the Ethernet Connection..... | 61 |
| Direct Computer to Controller Configuration..... | 61 |
| Wired or Wireless Network Configuration..... | 61 |
| Equipment Disposal (WEEE Directive) | 65 |
| Service & Technical Support | 65 |
| Warranty | 66 |
| Local VWR Offices in Europe and Asia Pacific | 67 |

Introduction

Thank you for choosing a VWR Signature Circulating Bath. It is intended for the precise temperature control of suitable Class I non-flammable or Class III flammable liquids (per DIN 12876-1) in a reservoir.

| | |
|---|--|
|  | <p>WARNING: VWR Circulating Baths are not intended for directly controlling the temperature of foods, pharmaceuticals, medicines, or other objects which may be ingested by or injected in humans or animals. Any such objects must be isolated from contact with the bath fluid and bath surfaces.</p> |
|---|--|

It will take you very little time to get your new Circulating Bath installed and running. This Operator's Manual is designed to guide you quickly through the process. We recommend that you read it thoroughly before you begin.

VWR Signature Circulating Baths with the Advanced Digital Temperature Controller

| Model Type | Reservoir Capacity | Temperature Range | |
|---------------------------------------|--------------------|-----------------------|-----------------------|
| | | °C | °F |
| AD07R-20 Refrigerating / Heating Bath | 7 liters | -20° to 200°C | -4° to 392°F |
| AD07R-40 Refrigerating / Heating Bath | 7 liters | -40° to 200°C | -40° to 392°F |
| AD7LR-20 Refrigerating / Heating Bath | 7 liters | -20° to 200°C | -4° to 392°F |
| AD15R-30 Refrigerating / Heating Bath | 15 liters | -30° to 200°C | -22° to 392°F |
| AD15R-40 Refrigerating / Heating Bath | 15 liters | -40° to 200°C | -40° to 392°F |
| AD20R-30 Refrigerating / Heating Bath | 20 liters | -30° to 200°C | -22° to 392°F |
| AD28R-30 Refrigerating / Heating Bath | 28 liters | -30° to 200°C | -22° to 392°F |
| AD45R-20 Refrigerating / Heating Bath | 45 liters | -25° to 135°C | -13° to 275°F |
| AD07H200 Heating Only Bath | 7 liters | Ambient +10° to 200°C | Ambient +20° to 392°F |
| AD15H200 Heating Only Bath | 15 liters | Ambient +10° to 200°C | Ambient +20° to 392°F |
| AD20H200 Heating Only Bath | 20 liters | Ambient +10° to 200°C | Ambient +20° to 392°F |
| AD28H200 Heating Only Bath | 28 liters | Ambient +10° to 200°C | Ambient +20° to 392°F |

General Safety Information

When installed, operated, and maintained according to the directions in this manual and common safety procedures, your Circulating Bath should provide safe and reliable temperature control. Please ensure that all individuals involved in the installation, operation, or maintenance of this Circulating Bath read this manual thoroughly prior to working with the unit.

| | |
|---|--|
|  | This symbol alerts you to a wide range of potential dangers. |
|  | This symbol advises you of danger from electricity or electric shock. |
|  | This symbol indicates that a hot surface may be present. |
|  | This symbol marks information that is particularly important. |
|  | This symbol indicates alternating current. |
|  | These symbols on the Power Switch / Circuit Breaker indicate that they place the main power supply ON / OFF. |
|  | This symbol on the Power Key indicates that it places the unit in a standby mode. It DOES NOT fully disconnect the unit from the power supply. |
|  | This symbol indicates a protective conductor terminal. |

**Read all instructions pertaining to safety, set-up, and operation.
Proper operation and maintenance is the user's responsibility.**

Safety Recommendations

To prevent injury to personnel and/or damage to property, always follow your workplace's safety procedures when operating this equipment. You should also comply with the following safety recommendations:

| | | |
|--|---|--|
|  | <p>WARNING:</p> <ul style="list-style-type: none"> This Circulating Bath is suitable for use with Class III flammable fluids per DIN 12876-1. A fire hazard may be present. Be aware of the chemical hazards that may be associated with the bath fluid used. Observe all safety warnings for the fluids used as well as those contained in the material safety data sheet. Explosive gas mixtures may accumulate if used with insufficient ventilation. Use this Circulating Bath in a well ventilated area or beneath a suitable fume hood only. Use only recommended bath fluids; see <i>Technical Information</i> in the rear of this manual for recommended fluids. Use only non-acid bath fluids. <p>WARNING: When using Class III flammable fluids per DIN 12876-1, the user must attach the following warning labels to the front of the unit so that they are well visible:</p> | |
| | <p>Warning Label W09 Colors: Yellow/black</p> |  |
| <p>Mandatory Label M018 Colors: Blue/white</p> <p>or</p> <p>Semi S1-0701 Table A1-2 #9 Colors: Blue/white</p> |   | <p>Carefully read the user information prior to beginning operation. Scope: EU</p> <p>Carefully read the user information prior to beginning operation. Scope: NAFTA</p> |

| | |
|---|--|
|  | <p>WARNING:</p> <ul style="list-style-type: none"> Always connect the power cord on this Circulator to a grounded (3-prong) power outlet. Make certain that the outlet is the same voltage and frequency as your unit. Never operate the Circulator with a damaged power cord. Always turn the Circulator OFF and disconnect mains power before performing any maintenance or service. |
|---|--|

| | |
|---|---|
|  | <p>WARNING:</p> <ul style="list-style-type: none"> Never operate the Circulator without bath fluid in the reservoir. Periodically check the reservoir to ensure that the liquid depth is within acceptable levels. Always refill the reservoir using the same bath fluid type that is already in the reservoir. Bath oil must not contain any water contaminants and should be preheated to the actual bath temperature before adding as there is an explosion hazard at high temperatures. Always drain all fluid from the reservoir before moving or lifting your Circulator. Be sure to follow your organization's procedures and practices regarding the safe lifting and relocation of heavy objects. |
|---|---|

**WARNING:**

- Always allow the bath fluid to cool to ambient temperature before draining.
- The reservoir cover, top deck, and/or external pump connections may become hot with continuous use. Exercise caution when touching these parts.



WARNING: It is the user's responsibility to properly decontaminate the unit in the event hazardous materials are spilled on exterior or interior surfaces. Consult manufacturer if there is any doubt regarding the compatibility of decontamination or cleaning agents.

Regulatory Compliance and Testing

This equipment is compliant with the European Directive 2002/95/EC and its latest amendments on Restrictions on Hazardous Substances (RoHS) and below the given limits of hazardous substances.

ETL Intertek (60 Hz units)

UL 61010-1 / CSA C22.2 No. 61010-1 — Safety Requirements for Measurement, Control, and Laboratory Use; Part 1: General Requirements

UL 61010A-2-010 / CSA C22.2 No. 61010-2-010:04 — Safety Requirements for Measurement, Control, and Laboratory Use; Part 2-010: Particular Requirements for Laboratory Equipment for the Heating of Materials

UL 61010A-2-051 / CSA C22.2 No. 61010-2-051:04 — Safety Requirements for Measurement, Control, and Laboratory Use; Part 2-051: Particular Requirements for Laboratory Equipment for the Mixing and Stirring

CE (all units)

EC Low Voltage Directive 2006/95/EC

EC Electromagnetic Compatibility Directive 2004/108/EC

IEC 61010-1-2001

IEC 61010-2-2001

IEC 61326:2005 / EN 61326 : 2006

Unpacking Your Circulator

Your Circulator was packed in a special carton or cartons. You should keep the packaging, along with all packing materials, until the unit has been installed and you are certain it is working properly.

| | |
|---|--|
|  | <p>CAUTION: Remove any loose packing material that may have fallen into the heater/pump housing during shipping. Before powering up, check that nothing remains around the heater or Circulator pump.</p> |
|---|--|

We recommend that you begin using your Circulator immediately to confirm proper operation, since beyond one week you may be eligible for warranty repair only (rather than replacement). You'll find complete warranty information in the back of this manual.

In the unlikely event that the unit was damaged or does not operate properly, contact the transportation company, file a damage claim, and contact the company where your Circulator was purchased.

Contents

The items included with your Circulator will vary depending on which model Circulating Bath you purchased.

Models with North American power cord:


| | Refrigerating / Heating Circulators | Heating Only Circulators |
|--|-------------------------------------|--------------------------|
| Resource Disk with Instruction Manual | 1 | 1 |
| Reservoir Lid | 1 (2 on 45L model) | 1 |
| 3-ft / 0.91 m IEC to IEC Power Cord | 1 | N/A |
| 6-ft / 1.82 m IEC to Mains Power Cord | 1 | 1 |
| Refrigeration Control Cable | 1 | N/A |
| 1/4 in. NPT to 3/16 in. barbed adapter | 2 | 2 |
| 1/4 in. NPT to 1/4 in. barbed adapter | 2 | 2 |
| 1/4 in. NPT to 3/8 in. barbed adapter | 2 | 2 |
| Certificate of Compliance | 1 | 1 |
| Quick-Start Guide | 1 | 1 |

Models with European power cords:



| | Refrigerating / Heating Circulators | Heating Only Circulators |
|--|-------------------------------------|--------------------------|
| Resource Disk with Instruction Manual | 1 | 1 |
| Reservoir Lid | 1 (2 on 45L model) | 1 |
| 3-ft / 0.91 m IEC to IEC Power Cord | 1 | N/A |
| 6-ft / 1.82 m IEC to Mains Power Cord – EU plug | 1 | 1 |
| 6-ft / 1.82 m IEC to Mains Power Cord – UK plug | 1 | 1 |
| 6-ft / 1.82 m IEC to Mains Power Cord – Swiss plug | 1 | 1 |
| Refrigeration Control Cable | 1 | N/A |
| 1/4 in. NPT to 3/16 in. barbed adapter | 2 | 2 |
| 1/4 in. NPT to 1/4 in. barbed adapter | 2 | 2 |
| 1/4 in. NPT to 3/8 in. barbed adapter | 2 | 2 |
| 1/4 in. NPT to M16 barbed adapter | 2 | 2 |
| Certificate of Compliance | 1 | 1 |
| Quick-Start Guide | 1 | 1 |

Controls and Components


Advanced Digital Controller



Refrigerating/Heating Baths











Heating Only Baths



Quick-Start

Unless otherwise specified, quick-start instructions apply to all models. See Installation and Startup for additional information.

| | | |
|-----------------|--|---|
| <p>1</p> | <p>Fill reservoir with fluid</p> |  <p>Maximum: 1 in. / 2.54 cm below underside of top deck</p> <p>Minimum: 3.0 in. / 7.6 cm below underside of top deck</p> |
| <p>2</p> | <p>Connect all power cords and control cables</p> | <div style="display: flex; justify-content: space-between;"> <div data-bbox="688 638 997 863">  <p>Heating only models</p> </div> <div data-bbox="688 915 1386 1325">  <p>IEC power cord from Controller to Refrigeration Power Module</p> <p>Refrigerating / Heating models</p> <p>Refrigeration control cable</p> </div> </div> |
| <p>3</p> | <p>Place Power Switch / Circuit Breaker in ON position</p> | <div style="display: flex; justify-content: space-between;"> <div data-bbox="789 1377 1040 1566">  <p>Heating only models</p> </div> <div data-bbox="870 1598 1386 1776">  <p>Refrigerating / Heating models</p> </div> </div> |

| | | |
|----------|------------------------------------|--|
| <p>4</p> | <p>Turn Controller "ON"</p> |  |
| <p>5</p> | <p>Set safety thermostat</p> |  |
| <p>6</p> | <p>Enter temperature set point</p> |  <p>6A. Press SET</p> <p>6B. Touch and hold or slide finger up/down scroll bar</p> |

Installation and Startup

Your Circulating Bath with Advanced Digital Temperature Controller is designed to be simple to set-up and install. The only tools required are a No.1 Phillips-head screwdriver and a container for adding water or other suitable fluid to the bath reservoir.

General Site Requirements

Locate your Circulator on a level surface free from drafts and direct sunlight. Do not place it where there are corrosive fumes, excessive moisture, high room temperatures, or in excessively dusty areas.

Refrigerating / Heating Circulators must be 10.2 cm / 4 inches or more away from walls or vertical surfaces so that airflow is not restricted.


Avoid voltage drops by using properly grounded power outlets wired with 14 gauge or larger diameter wire and if possible, be close to the power distribution panel. The use of extension cords is not recommended; this will reduce the potential for problems caused by low line voltage.

Adding Liquid to the Bath Reservoir

| | | | |
|--|---|---|--|
|  | <p>WARNING: When using Class III flammable fluids per DIN 12876-1, the user must attach the following warning labels to the front of the unit so that they are well visible:</p> | | |
| | <p>Warning Label W09 Colors: Yellow/black</p> |  | <p>Danger Area. Attention! Observe instructions (operating manual, safety data sheet)</p> |
| | <p>Mandatory Label M018 Colors: Blue/white</p> <p>or</p> <p>Semi S1-0701 Table A1-2 #9 Colors: Blue/white</p> |   | <p>Carefully read the user information prior to beginning operation. Scope: EU</p> <p>Carefully read the user information prior to beginning operation. Scope: NAFTA</p> |

| | |
|---|--|
|  | <p>WARNING: Read the safety data sheet for the bath fluid being used carefully before filling reservoir.</p> <p>WARNING: Always use fluids that satisfy safety, health, and equipment compatibility requirements.</p> <p>WARNING: If the proper fluid level is not maintained, the heater coil may become exposed and possibly damaged (fluid level too low) or the bath may overflow (fluid level too high).</p> |
|---|--|

The liquid in the reservoir should be maintained at a depth between 1 inch / 2.54 cm and 3.0 inches / 7.6 cm below the underside of the bath's top deck. Upon start up, it may be necessary to add fluid to the bath to compensate for the fluid required for external circulation. Likewise, be sure to compensate for fluid displacement when placing samples or other materials in the Circulator's reservoir.



WARNING: Always drain all fluid from the reservoir before moving or lifting your Circulator. Be sure to follow your organization's procedures and practices regarding the safe lifting and relocation of heavy objects.



WARNING: To avoid the potential for burns, allow the Circulator to cool completely before cleaning or performing any maintenance.

Pump Inlet and Outlet Connections



WARNING: When connecting tubing to an external application, it is the user's responsibility to make sure that the tubing and fittings connected to the Circulator are suitable for the fluid being used and the temperature range of operation.

CAUTION: The Circulator's bypass tubing is secured to the fluid inlet and outlet connections by high temperature nylon hose clamps, which can be removed by carefully cutting them with diagonal cutters.

CAUTION: Secure the tubing to the inlet and outlet fittings using hose clamps with a minimum ID of 7/8 inch (22 mm). Do not operate the unit without hose clamps.



WARNING: If the Circulating Bath will not be used for external circulation, the inlet and outlet ports should remain connected using the Buna N bypass tubing provided with the unit.

The pump inlet and outlet ports are female ¼ inch NPT connections that permit use of barbed tubing adapters or hard plumbing fittings. ½ inch (13mm) ID tubing may also be slid over these connections and held in place with a hose clamp (minimum 7/8 inch / 22 mm ID).

If the pump inlet and outlet are not used for external circulation, the Bypass Tubing provided with the unit should be left in place in order to optimize fluid mixing within the reservoir.

The nylon barbed tubing adapter fittings supplied with the unit are intended for applications from -40° to 93°C. For applications above 93°C, brass, stainless steel, or Teflon® fittings are recommended. ¼ inch NPT to M16 stainless steel male adapter fittings are provided with all 50Hz models.



NOTE: The use of quick-connect fittings is not recommended as they typically restrict flow rate.

External Closed Loop Circulation

Connect the pump inlet and outlet to the external apparatus. To maintain adequate flow, avoid restrictions in the tubing. When connecting the Circulator to more than two closed loops, the use of a manifold made of “Y” adapters to divide the fluid into multiple banks is recommended. After setting up multiple closed loops, check for adequate flow at the return manifold of each loop and check that the bath fluid is at an adequate level. A booster pump may be added to closed loops without damaging the Circulator’s pump.

The temperature control stability of a closed loop system is better at the external apparatus than in the Circulator reservoir (provided the control point of the apparatus represents a constant load and is well insulated). For example, if you circulate fluid through a viscometer at 50°C, the temperature variation observed in the Circulator reservoir may be $\pm 0.1^\circ\text{C}$ while the temperature variation in the viscometer may be only $\pm 0.05^\circ\text{C}$.

Although temperature stability is generally better at the external apparatus control point, depending on the length of tubing used and the efficiency of the insulation, the actual temperature reading at the external apparatus may be slightly different than the temperature reading at the Circulator reservoir.


Open Loop Circulation

The duplex (pressure/suction) pump permits circulation to and from an external open bath. To prevent siphoning when the Circulating Bath is turned off, position both baths so that the two fluid levels are at the same elevation.


Connect the pump inlet and outlet to the external bath using tubing of the same diameter and length. The same size fittings should also be used on both the inlet (suction) and outlet (pressure). This helps ensure a balanced flow. A restricting valve or pinch clip should be installed in the pressure (outlet) tubing and adjusted to match the return suction (inlet) flow rate. Cut the external end of the suction tube into a “V” shape so that the tube will not seal itself against the wall of the external tank. Both the pressure and suction tubing should be securely fastened to the external tank to prevent movement during use.

When using flexible tubing, the suction tubing must have a wall thickness that will not collapse under vacuum, particularly when going around bends.

Circulating Bath Height Regulation — Position the ends of the pressure and suction tubes at the desired maximum fluid level in the external bath and fill the bath to that level. Fill the Circulating Bath to a height one inch (25mm) below the top of the reservoir. Start the pump and adjust the restricting valve/pinch clip on the pressure tubing until the liquid height in both baths remains constant. Add fluid to the baths as needed to compensate for the fluid in the inlet and outlet lines.



Refrigeration Control Connections (Refrigerating/Heating Circulators only)



Electrical Power

| | |
|---|---|
|  | WARNING: The Circulator's power cord must be connected to a properly grounded electrical receptacle. Make certain that this electrical outlet is the same voltage and frequency as your Circulator. The correct voltage and frequency for your Circulator are indicated on the identification label on the back of the Controller. |
|---|---|

| | |
|---|--|
|  | CAUTION: The use of an extension cord is not recommended. If one is necessary, it must be properly grounded and capable of handling the total wattage of the unit. The extension cord must not cause more than a 10% drop in voltage to the unit. |
|---|--|


Refrigerating / Heating Circulators

Attach the 3-ft / 0.91 m power cord to the IEC electrical connectors on the Temperature Controller (male) and the Refrigeration Power Module (female).

Attach the 6-ft / 1.8 m power cord to the IEC electrical connection on the Refrigeration Power Module and then plug the male connector into the Mains electrical outlet.

Place the Power Switch / Circuit Breaker on the Refrigeration Power Module in the ON position. The LCD on the Controller will light and "Standby" will appear on the display; the Power Key will also light.

| | |
|---|--|
|  | NOTE: To conserve power when not in use, the LCD's backlighting will go out about 5 seconds after "Standby" appears. The Power Key will remain lit to indicate that the Controller is energized and ready to use. |
|---|--|



Heat Only Circulators

Attach the 6-ft / 1.8 m power cord to the IEC electrical connection on the Temperature Controller and then plug the male connector into the Mains electrical outlet.


Place the Power Switch / Circuit Breaker on the Temperature Controller in the ON position. The LCD on the Controller will light and "Standby" will appear on the display; the the Power Key will also light.



NOTE: To conserve power when not in use, the LCD will go black about 5 seconds after "Standby" appears. The Power Key will remain lit to indicate that the Controller is energized and ready to use.

Communication

The Advanced Digital Controller features a variety of connectivity options. Following are some typical ways you can use them to monitor and control the operation of your Circulator.



USB Communication

Two USB ports (A and B) are provided on the rear of the Temperature Controller. USB A is intended for use with a flash drive and allows you to easily log temperature data. USB B allows you to remotely monitor and control your Circulator using a computer. See *Normal Operation, Selecting a Remote Communication and Control Protocol* and the *Technical Information* section of this manual for additional information.

Ethernet

An Ethernet port is provided on the back of the Temperature Controller to enable you to connect your Circulator to a computer network. See *Normal Operation, Selecting a Remote Communication and Control Protocol* and the *Technical Information* section of this manual for additional information,

RS232 / RS485 Serial Communication



CAUTION: Always turn electrical power to the Circulator OFF before making a connection to the serial (DB9) port.


Your Circulator features RS232 / RS485 serial communication for remote data logging and control capability. A DB9 connector is provided on the rear of the Temperature Controller for this purpose. See *Normal Operation, Selecting a Remote Communication and Control Protocol* and the *Technical Information* section of this manual for additional information.

The serial interface should be connected to a serial communication port on a remote PC using an appropriate cable. Information on the RS232 / RS485 command and communication protocol can be found in the *Technical Information* section of this manual.

External (P2) Temperature Probe

Your Circulator is capable of controlling temperature based on either the temperature of the internal bath or that of an external vessel or device. The connection for the optional external temperature probe is on the rear of the Temperature Controller. The Temperature Controller automatically detects the external temperature probe when it is connected. See *Replacement Parts & Accessories* for available lengths and part numbers.

Pin Out Diagrams — External (P2) Temperature Probe Connection



RTD SENSOR: 4 WIRE CIRCUIT, 100 OHMS @ 0 DEGREES C,
MAXIMUM OPERATING TEMPERATURE @ 200 C, CLASS A 0.003850 OHMS/DEGREES C.

Controller Setup

Power

Press . The Circulator will begin running, actual and set point temperatures will be displayed, and the word "SET" will be continuously lit. The pump symbol will also be lit and the heating or refrigerating symbol may be lit or flashing.



Internal (P1) control only – external probe (P2) not connected

Safety Set Temperature


This is a “Do Not Exceed” temperature setting for your Circulator and is the temperature at which the heater will be turned OFF should the liquid level in the bath drop too low or the heater malfunctions. It is normally set about 5° higher than the desired operating temperature. Setting the Safety Set temperature is a simple 3-step procedure.

| | |
|--|--|
| | WARNING: The Safety Thermostat is user-adjustable from approximately 100° to 220°C / 312° to 428°F. Do not force the indicator dial beyond the stops at either end of the dial's range. |
|--|--|

1. Press the  key until SAFETY appears.



2. Using a No.1 Philips head screwdriver, rotate the Safety Thermostat on the rear of the Temperature Controller until the desired Safety Set Temperature is displayed (clockwise to increase; counter-clockwise to decrease).



3. Press  to return to the main operational screen.

Normal Operation

Keys and Controls

| | | |
|------------------|--|---|
| Power |  | Turns the Circulator's Temperature Controller ON. |
| Home |  | Returns the LCD to the Main Operational Display (from any screen). |
| Menu |  | Accesses the Temperature Controller's set-up sub-menus. The items in these sub-menus are used to configure the Controller's general operational parameters (temperature unit, pump speed, upper and low temperature limits, etc. (see Set-up Sub-Menus, below). |
| SET |  | Used in conjunction with the Touch Scroll Bar to change the set point temperature and some operational parameters. |
| Touch Scroll Bar |  | Used to make temperature set point and other operational changes. Slide finger up / down scroll bar or touch upper / lower sections to make minor adjustments; hold your finger on the scroll bar to make large adjustments. |

Turning Your Circulator ON

Press the  key.

When the Circulator begins running, the actual and set point temperatures will be displayed and the pump symbol will be lit. If the optional external temperature probe is connected, a temperature value will also be displayed along the bottom of the screen. This is the bath temperature measured by the non-controlling (monitor) probe.

If the actual bath temperature is lower than the set point temperature, the heat symbol will also be lit.

Refrigerating/Heating Models: If the actual bath temperature is higher than the set point temperature, the refrigerating symbol will be lit. It is normal for both the heating and refrigerating symbols to be lit simultaneously when nearing or maintaining the set point temperature.

Main Operational Display (Home)

This is the Circulators main operational display. You can return to this screen at any time by pressing the  key.



Internal (P1) control only – external probe (P2) not connected



Internal (P1) control – external probe (P2) connected



External (P2) control

Set-Up Sub-Menus

Pressing the  key accesses and scrolls through the Temperature Controller's set-up sub-menus. The Touch Scroll Bar is used to change the current setting / value in the sub-menus.

| Sub-Menu | Description | Selection / Range | Factory Default |
|----------------------|--|--|-----------------|
| UNIT | Unit in which temperature is set and displayed | °C or °F | °C |
| PUMP | Pump speed | Variable – 5 to 100 | 40 |
| HI LIM | High limit temperature | 25° to 202°C / 80° to 395°F | 202°C / 395°F |
| LO LIM | Low limit temperature | -52° to 20°C / -65° to 65°F | -52°C / -65°F |
| SET DIFF | Maximum allowable difference between the measured internal and external temperatures | 2° to 20°C | 10°C |
| AUTOCOOL | Temperature at which refrigeration is activated | 1° to 150°C | 45°C |
| SHC | The specific heat capacity of the bath fluid | 0.10 to 2.00 | 1.00 |
| CONTROL | Internal or external temperature control | P1 (Internal) or P2 (External) | P1 |
| OFFSET / CALIBRATION | Calibration or display offset value Password required | -3.0°C to +3.0°C | 0.0°C |
| SAFETY | Safety Set Temperature | ~40° to 240°C | N/A |
| COM | Remote communication and control | RS232, RS485, USB, Ethernet | RS232 |
| TIMER | Count down timer with audible signal | 1 second to 999 minutes, 59 seconds | 000:00 |
| PASSWORD | Permits access to Calibration and Diagnostic screens | 10 = Diagnostic screen 12 = Access Calibrate | 0 |
| USB LOG | Turns data logging On and Off | No (Off) or Yes (On) | No (Off) |
| ENGLISH | Language in which information will be displayed | English, French, German, Spanish | English |
| pODO | Operational hours on pump | Display only | N/A |
| cODO | Operational hours on compressor | Display only | N/A |
| AUTO-ON | Determines how unit will restart after a disruption in electrical power | ON or OFF | OFF |

To accept a value in a sub-menu, press , , , or allow the LCD to return to the main operational display (approximately 10 seconds).

Adjusting the Temperature Set Point

This is the temperature at which the fluid in your Circulating Bath will be maintained. It may be set to one-hundredth of a degree over a range of -50.00° to +200.00°C / -60.00° to +390°F. The factory default set point is +20.0°C / +68.0°F.



Touch and hold up/down arrows to make large changes

Touch and release arrows or slide finger up/down scroll bar to make small changes

Press SET once to change whole degree values

Press SET twice to make changes of less than a degree (decimal point will flash)

To Change: Press **SET**. The arrow around the word "SET" will begin flashing. To make changes of one degree or more, touch the scroll bar until the desired set point temperature is displayed. To make changes of less than a degree, (e.g., 0.50°C), press **SET** a second time. The decimal point will begin flashing. Touch the scroll bar until the desired value is displayed.

To Accept: Press **Menu**, **Home**, **SET**, or allow the LCD to return to the main operational display (approximately 10 seconds).



NOTE: An audible alarm and the words Low Limit or High Limit flashing on the display indicate that the temperature set point value is outside the Low Limit or High Limit value. The Circulator will continue to heat/cool until the actual bath temperature reaches the Limit value, at which point operation will stop.

Selecting the Temperature Unit

The temperature units sub-menu (°C / °F) allows you to select the temperature unit in which the actual bath temperature and set point temperature are displayed. The factory default is °C.



To Access: Press the  key until UNIT is displayed.

To Change: To select °F, touch the bottom portion of the scroll bar; to select °C, touch the top portion of the scroll bar.

To Accept: Press , , , or allow the LCD to return to the main operational display (approximately 10 seconds).

Selecting the Pump Speed

This sub-menu allows you to set your Circulator's pump speed. The pump speed range is 5 to 100; the factory default is 40.



To Access: Press the  key until PUMP is displayed.

To Change: Touch the scroll bar until the desired pump speed is displayed (5 minimum, 100 maximum).

To Accept: Press , , , or allow the LCD to return to the main operational display (approximately 10 seconds).

Setting the High Limit Temperature

This sub-menu allows you to limit how high the temperature set point may be set. It also serves as a high limit safety, alerting you if bath temperature rises above the high limit temperature setting. The High Limit value may be set from +25° to +202°C / +80° to +395 °F; the factory default is 202°C.

To avoid an unwanted shutdown during regular operation, the High Limit value should be set at least 5° higher than the selected control temperature.



Touch and hold up/down arrows to make large changes

Touch and release to make minor changes

To Access: Press the  key until HI LIMIT is displayed.

To Change: Touch the scroll bar until the desired high limit temperature is displayed.

To Accept: Press , , , or allow the LCD to return to the main operational display (approximately 10 seconds).

Setting the Low Limit Temperature

This sub-menu allows you to limit how low the temperature set point may be set. It also serves as a low limit safety, alerting you if bath temperature falls below the low limit temperature setting. The Low Limit value may be set from -52°to +20°C / -65°to +65°F ; the factory default is -52°C.

To avoid an unwanted alarms or shutdown during regular operation, the Low Limit value should be set at least 5°lower than the selected control temperature.



Touch and hold up/down arrows to make large changes

Touch and release to make minor changes

To Access: Press the  key until LO LIMIT is displayed.

To Change: Touch the scroll bar until the desired low limit temperature is displayed.

To Accept: Press , , , or allow the LCD to return to the main operational display (approximately 10 seconds).

Setting the Differential Temperature

This sub-menu is for use when external temperature control (P2) is being used. It establishes the maximum allowable differential between the external and internal measured temperatures and is intended prevent uncontrolled heating or cooling should the external temperature probe be dislodged or fail. The differential temperature may be set from 2° to 20°C ; the factory default is 10°C.



Touch and hold up/down arrows o make large changes

Touch and release to make minor changes

To Access: Press the  key until SET DIFF is displayed.

To Change: Touch the scroll bar until the desired until the desired P2 - P1 (external – internal) temperature differential displayed.

To Accept: Press , , , or allow the LCD to return to the main operational display (approximately 10 seconds).



CAUTION: The Differential Temperature value is always shown in degrees C, even if degrees F is selected as the temperature unit in which the control and actual bath temperatures are displayed.

Setting the Auto Cool Temperature

This sub-menu is displayed only on Refrigerating / Heating Circulators. It determines the set point temperature at which refrigeration will be activated and permits more precise control when operating at high temperatures as well as more rapid cool downs. For most applications, a set point that is 15°C above room temperature is recommended. The Auto Cool control range is from +1°C to 150°C. The refrigeration system will turn on when the bath temperature set point (150°C maximum) is below the Auto-Cool set point. The factory default is 45°C.

Cool Command™ Refrigeration — -40°C 7 liter Refrigerating/Heating Circulators and 15 liter and larger Refrigerating/Heating Circulators feature the Cool Command™ modulating refrigeration control system. Cool Command allows the refrigeration system to turn on at a fluid temperature up to 150°C when the temperature set point is changed to or below the Auto Cool set point (150°C maximum). As a result, bath fluid cools more quickly.

Conventional Refrigeration — -20°C 7 liter Refrigerating/Heating Circulators use a conventional refrigeration system. The refrigeration system will turn on when the bath fluid temperature and set point are below the Auto Cool set point (70°C maximum).



Touch and hold up/down arrows to make large changes

Touch and release to make minor changes

To Access: Press the  until AUTOCOOL is displayed.

To Change: Touch the touch scroll bar or touch the arrow symbols until the desired auto cool temperature is displayed.

To Accept: Press , , , or allow the LCD to return to the main operational display (approximately 10 seconds).

Setting the Specific Heat Capacity (SHC)

This sub-menu allows you to tune the Circulator's control algorithm for the specific heat capacity (SHC) of the fluid you are using. The factory default is 1.00.



Touch and hold up/down arrows to make large changes

Touch and release to make minor changes

To Access: Press the  key until SHC is displayed.

To Change: Touch the touch scroll bar or touch the arrow symbols until the desired specific heat capacity value is displayed.

To Accept: Press , , , or allow the LCD to return to the main operational display (approximately 10 seconds).

Selecting Internal or External Temperature Control

This sub-menu lets you determine whether temperature control will be based on the fluid temperature within the Circulator's reservoir or the fluid temperature at an external device.



To Access: Press the  key until CONTROL is displayed.

To Change: Touch the scroll bar until the desired temperature probe is displayed. P1 = Internal; P2 = External.

To Accept: Press , , , or allow the LCD to return to the main operational display (approximately 10 seconds).



If P2 is selected but an external temperature probe is not connected to the Temperature Controller, P2 UNPLUGGED will appear on the display.

Calibrating Your Circulator (Offset)

These two sub-menus allow you to match the Circulator's internal and/or external temperature probe to an external reference thermometer. A value from -3.0°C to +3.0°C may be entered; the factory default is 0.0°C.



CAUTION: The Offset Calibration value is always shown in degrees C, even if degrees F is selected as the temperature unit in which the control and actual bath temperatures are displayed. Your Circulator will automatically convert the °C offset calibration value to the correct °F display offset value.



NOTE: If you attempt to calibrate the Circulator when OFFSET is displayed, **Password!** will appear at the bottom of the screen. To proceed, advance to the PASSWORD sub-menu, enter 12, and then return to this sub-menu.



Protected – Password Required



Touch top arrow to increase offset

Touch bottom arrow to decrease offset

Unlocked

To Access: Press the  key until the calibration sub-menu for the temperature probe you wish to calibrate (P1 or P2) is displayed. If this function has been locked, OFFSET will appear and must be unlocked before changes can be made; if this function has been unlocked, CALIBRATE will appear.

To Change: Touch the scroll bar or touch the arrow symbols until the desired calibration temperature displayed.

To Accept: Press , , , or allow the LCD to return to the main operational display (approximately 10 seconds).

Displaying the Safety Set Temperature

This display shows the current Safety Set Temperature.



NOTE: If you attempt to set the Safety Set temperature using the scroll bar, **Rear Adjust** will appear at the bottom of the screen. Use the Safety thermostat on the rear of the Temperature Controller to adjust the Safety Set temperature.



To Access: Press the  key until SAFETY is displayed.

To Change: This is a display only. The Safety Set Temperature is changed using the Safety thermostat on the rear of the temperature controller. The temperature value shown on the display will change as the thermostat is adjusted.

To Exit: Press , , , or allow the LCD to return to the main operational display (approximately 10 seconds).

Selecting a Remote Communication and Control Protocol

This sub-menu allows you to select the protocol you want to use for remote communication and control. The choices are RS232, addressable RS485, USB, and Ethernet; the factory default is RS232.



Touch top arrow for Ethernet

Touch areas between the arrows for USB and RS485

Touch bottom arrow RS232

To Access: Press the  key until COM and the currently selected protocol appears on the display.

To Change: To select RS232, touch the down arrow; to select Ethernet, touch the top arrow. RS485 and USB protocols are selected by touching areas on the scroll bar that lie between the two arrows.

To Accept: Press , , , or allow the LCD to return to the main operational display (approximately 10 seconds).



CAUTION: Although you can accept the displayed protocol by pressing either the Home, Menu, or SET keys, we recommend using the Menu key the first time a protocol is selected as there are additional parameters that must be configured when RS232, RS485, or Ethernet is selected. This will ensure that you are prompted to configure those parameters. See *Technical Information, RS232/RS485 Configuration* and *Technical Information, Ethernet Configuration* for additional information.



NOTE: See *Technical Information, USB Setup, Monitoring, and Control* for information on controlling your Circulator via the USB-B port.

Using the Timer

This sub-menu allows you to set and activate the Temperature Controller's count-down timer. Any time duration from 1 second to 999 minutes, 59 seconds can be entered.



Touch and hold up/down arrows to set long durations

Touch and release set short durations changes

To Access: Press the  key until TIMER is displayed.

To Change: Touch the scroll bar until the desired time duration is displayed.

To Start the Timer: Press , The timer will begin counting down.

To Pause the Timer: Press . The display will revert to the main operational display.

To restart the timer, return to the Timer sub-menu and press .

The timer counts down in one second increments until it reaches 000:00. At this point, the time display flashes and an audible signal is activated. Press  to acknowledge and deactivate.



NOTE: You can access other menu functions while the timer is running without affecting the count down.

Entering a Password

This sub-menu allows access to the Calibration and Diagnostic sub-menus.



To Access: Press the  key until PASSWORD is displayed.

To Change: Touch the scroll bar until the password needed to access the desired sub-menu is displayed. Calibration = 12; Diagnostic = 10)

To Accept: Press , , , or allow the LCD to return to the main operational display (approximately 10 seconds).



IMPORTANT: Once a password has been entered and accepted, the function it enables stays unlocked until either a new password is entered or the Circulator is turned OFF.

Enabling / Disabling Data Logging

This sub-menu allows you to turn data logging On and Off.



IMPORTANT: This Advanced Digital Temperature Controller has both USB A and USB B ports. USB A is used for data logging; USB B is used for remote monitoring and control. Data is stored in a CSV file that can be read in spreadsheet programs such as Microsoft Excel®. To begin data logging, plug a flash drive into the USB A port and then set USB LOG to YES.



Touch top arrow to turn data logging On (Yes)

Touch bottom arrow to turn data logging Off (No)

To Access: Press the  key until USB LOG is displayed.

To Change: Touch the scroll bar until the desired status is displayed.


To Accept: Press , , , or allow the LCD to return to the main operational display (approximately 10 seconds).



NOTE: See the *Technical Information* for complete information regarding data logging file formats and content.

Selecting the Operational Language

This sub-menu allows you select the language in which information will be displayed. The choices are English, Spanish, French, and German; the factory default is English.



To Access: Press the  key until the current language is displayed.

To Change: Touch the scroll bar until the desired language is displayed.

To Accept: Press , , , or allow the LCD to return to the main operational display (approximately 10 seconds).

Setting Auto Restart



WARNING: The unit may start automatically after a disruption in electrical power.

This sub-menu allows you to select how the unit will begin operating after a disruption in electrical power. When Yes is selected, the Circulator will begin running automatically when power is restored. When No is selected, the Circulator will power up in the Standby mode.



Touch top arrow to turn Auto-On ON (YES)

Touch bottom arrow to turn OFF (NO)

To Access: Press the  key until AUTO-ON is displayed.

To Change: Touch the top of the scroll bar to select Yes (ON); touch the bottom of the scroll bar to select No (OFF).

To Accept: Press , , , or allow the LCD to return to the main operational display (approximately 10 seconds).

Resetting the Factory Default Values

To reset your Circulator to its original factory default values, proceed as follows:

1. Press the  key to place the unit in Standby.
2. Place the Power Switch / Circuit Breaker in the Off position.
3. Return the Power Switch / Circuit Breaker to the On position while pressing the  key until "STANDBY" appears on the display.
4. Press the  key.

Changing Your Circulator's Viewing Angle


Your Circulator is equipped with Swivel 180™, an innovative feature which permits viewing of the temperature display from anywhere within a 180° radius.



NOTE: There are positive stops at 45° intervals; however, the viewing angle may be set anywhere within a 180° radius.

To change the viewing angle, slide the release latch to the right and rotate the Temperature Controller to the desired angle. The latch release will automatically return to the locked position at every 45° positive stop.

Swivel 180™ latch release



Inert Gas Purge

A 0.125 in. / 3 mm port on the rear of the Temperature Controller is provided to allow you to blanket the surface of the liquid in the bath reservoir with nitrogen or another inert gas to help prevent condensation and dilution of the bath fluid.

Tap Water Cooling

Tap water cooling allows for more rapid bath cool down from high temperatures and/or more precise operation at temperatures near ambient.

Heat only Circulating Baths feature an integrated cooling coil as standard equipment. The tap water connections are made on the rear of the unit. Two 0.25 inch / 6.4 mm female NPT fittings are provided for these connections.




WARNING: The fluid outlet must be connected and flow to a suitable drain or vessel located at a level below that of the inlet.

Reservoir Cover Storage

Refrigerating/Heating and Heat Only Circulating Baths feature the LidDock® system to eliminate mess when adding fluid or samples to the reservoir. Specially positioned notches in the inner lip of the top deck allow you to stand the reservoir cover up upright over the bath opening, allowing condensate to flow back into the bath.



Display Messages and Alarms

| Message and/or Symbol | Description | Corrective Action |
|--|---|--|
| POWER FAILED | Informational: Indicates that electrical power was lost during operation; appears only when Auto-On is set to Yes. | Press  key to clear the message. |
| WARNING! LO LIMIT! | Warning: The fluid temperature or temperature set point is below the Low Limit value. (Message flashes, audible beep) | Decrease the Low Limit temperature value or increase the temperature set point. If problem is not corrected within about 30 seconds, the Circulator will go into a Low Limit Fault condition and operation will cease. |
| WARNING! HI LIMIT! | Warning: The fluid temperature or temperature set point is above the High Limit value. (Message flashes, audible beep) | Increase the High Limit temperature value or decrease the set point temperature. If problem is not corrected within about 30 seconds, the Circulator will go into a High Limit Fault condition and operation will cease. |
| FAULT! LO LIMIT! | Fault: The bath temperature has fallen below the Low Limit temperature value. Power to the compressor and pump will remain OFF until the problem is corrected. (Message flashes, continuous tone) | Press  to turn power OFF. Restore power and then decrease the Low Limit temperature value or increase the temperature set point. Controller failure; consult factory. |
| FAULT! HI LIMIT! | Fault: The bath temperature risen has above the High Limit temperature value. Power to the heater and pump will remain OFF until the problem is corrected. (Message flashes, continuous tone) | Press  to turn power OFF. Restore power and then increase the High Limit temperature value or decrease the temperature set point. Controller failure; consult factory. |
|   | Fault: The liquid in the bath has dropped too low or the temperature of the bath fluid has exceeded the Safety Set temperature. Power to the heater will remain OFF until the problem is corrected. | Fluid level in reservoir has fallen below minimum level; add fluid as required. Fluid temperature is higher than Safety Set temperature; increase Safety Set temperature setting. Controller failure; consult factory. |
| FAULT! EXT PROBE 2 | Fault: The external temperature probe has been disconnected. Appears only when using External Control. (Message flashes, continuous tone) | Reconnect external probe, turn power OFF and then back ON. If problem persists, replace external probe or operate using Internal Control. |

Informational Messages — Do not disrupt normal operation. Clear by pressing the  key.

Warnings — Circulator operation continues unless left uncorrected for approximately 30 seconds. Press the  key to silence the audible signal. Correct problem and press the  key a second time to clear Warning.

Faults — Circulator operation is halted (heater, pump, and compressor turn OFF). Press the  key to silence the audible signal. Press  to turn power OFF and then restore power and correct the problem.

Routine Maintenance and Troubleshooting



WARNING: Always turn your Circulator OFF and disconnect it from the electrical power outlet before performing any maintenance or service.



WARNING: To avoid the potential for burns, allow the Circulator to cool completely before cleaning or performing any maintenance.

Maintaining Clear Bath Water

Optimum temperature and moisture conditions for algae growth exist when using water as a bath fluid. To prevent algae contamination and minimize the frequency of draining the reservoir, an algaecide should be used.



WARNING: Do not use chlorine bleach.

Draining the Bath Reservoir



WARNING: Bath fluids should be stored and disposed of according to applicable laws and regulations.

Refrigerated / Heating and Heat only Circulating Baths


Refrigerated / Heating and Heat only Circulating Baths are equipped with a drain valve and port located either beneath the front access panel or on the right hand side of the unit.

To drain fluid from the bath, attach a short length of suitable 0.45 inch / 11.5 mm ID tubing to the drain port and secure it using a hose clamp with a minimum ID of 0.7 inch / 18 mm. Open the drain valve using a flat blade screwdriver. When closing the valve, do not over tighten.



WARNING: Be sure to close the drain valve before refilling the bath reservoir. Do not over tighten.

Drain Valve



Drain Port

Checking the Over-Temperature / Low Liquid Level Safety Systems

Your Circulator incorporates over-temperature and low liquid level protection according to IEC 61010-2-010. For optimum safety, these systems should be checked at least every six months for proper operation. These checks must be performed with the unit running.

Over-Temperature Protection


1. Press  until the Safety sub-menu is displayed.



Safety Set Temperature

2. Using a No.1 Philips head screwdriver, rotate the Safety Set Thermostat on the rear of the Temperature Controller until the unit shuts down. The Safety Set temperature at this point should be the same as the actual bath temperature.
3. Return the Safety Set temperature to the desired over-temperature value.
4. Press the  key to return to the Main Operational Display.

Safety Set
Thermostat



Low Liquid-Level Protection

1. Set the temperature set point to ambient and allow the Circulator to stabilize at that temperature.
2. Increase the temperature set point to about 5°C above ambient and slowly drain fluid from the bath.
3. Continue draining fluid until the unit shuts down. The fluid level at this point should be approximately 3.75 in. / 9.5 cm below the underside of the Circulator's top deck.
4. Replace the bath fluid and return to normal operation.

Cleaning Your Circulator



WARNING: It is the user's responsibility to properly decontaminate the unit in the event hazardous materials are spilled on exterior or interior surfaces. Consult the manufacturer if there is any doubt regarding the compatibility of decontamination or cleaning agents.

Temperature Controller

Turn the Temperature Controller OFF by pressing  and unplug power cord from the electrical outlet.

Wipe the housing with a clean cloth dampened with a mild detergent and water or mild all-purpose cleaner.



CAUTION: Do not spray cleaning liquids directly onto the Temperature Controller or allow them to enter the Controller's vents. Do not use abrasives as these could scratch the housing or the digital display.

Bath Reservoir

Bath Reservoir and Wetted Components — A vinegar solution or calcium/lime/rust remover can be used to clean mineral deposits from the Temperature Controller's wetted parts. The cleaner should be added to the bath reservoir at the prescribed dosage and circulated at 60°C / 140°F until the scale is removed.



CAUTION: Do not use steel wool to clean your Circulator's bath reservoir.

External Surfaces — Only mild detergents and water or an approved cleaner should be used on the top deck and other external surfaces of your Circulator. Do not allow cleaning liquids or sprays to enter the vents on the rear of the Temperature Controller.

Pump Impeller

In the unlikely event that debris becomes lodged in the pump impeller, a soft brush can be used to remove any lodged particles. If necessary, soak in a solution of distilled water and vinegar to soften before brushing.




CAUTION: Do not use hard utensils or abrasive pads to remove trapped debris.

Condenser, Air Vents, and Reusable Filter

(Refrigerating / Heating Circulators only)

To keep the refrigeration system operating at optimum cooling capacity, the condenser, removable air filter, and all air vents (front, side, back) should be kept free of dust and dirt. Be sure to check them on a regular basis and clean as required.

The reusable filter is easily accessed from the front of the unit by simply removing the access panel. Use a mild detergent and water solution to wash off any accumulated dust and dirt. Rinse and dry thoroughly before reinstalling.



Temperature Controller Removal and Re-Installation

Removal

The Temperature Controller on your Circulating Bath is designed to be easily removed from the top deck without the use of special tools. It is removed as follows:


1. Place the tip of a small flat blade screwdriver under the retaining ring locking tab and pry up gently.



2. Rotate the Temperature Controller clockwise until it stops (about 0.75 inch / 1.9 cm).




3. Lift the Controller straight up and out of the opening of the Circulator's top deck.



Re-Installation

The top deck of your Circulator incorporates four pins to facilitate positioning of the Temperature Controller when it is being reinstalled. These pins correspond to keyhole slots on the interior of the Circulator's retaining ring.

1. With the retaining ring locking tab oriented above one of the indents on the top deck, slowly lower the Temperature Controller into the top deck opening until it is resting on top of the positioning pins.



2. Gently rotate the Temperature Controller until it drops down on the positioning pins.



3. Rotate the Temperature Controller in the counter-clockwise until the Locking Tab engages the indent on the top deck.



Viewing Component Operating Time

Your Circulator tracks pump and (on Refrigerating/Heating Circulators) compressor operation, This run time is initially tracked in hours and minutes and then in days. You can view this information by pressing the  key until the pODO or cODO screen appears:



Pump Operating Time – (shown here in days)



Compressor Operating Time – (shown here in hours and minutes)

To Exit: Press , , , or allow the LCD to return to the main operational display (approximately 10 seconds).

Troubleshooting Chart

| Problem | Possible Causes | Corrective Action |
|---|--|--|
| Unit does not run (Digital Display is blank) | No power to unit | Check that the electrical cord is secure and connected to an operating electrical outlet. |
| Unit does not run (STANDBY appears on Digital Display) | Unit in Standby mode | Press Power Key on front panel. |
| No fluid circulation | Insufficient fluid in reservoir Pump impeller jammed | Add fluid to reservoir. Inspect pump and remove debris as required. |
| Insufficient circulation | Fluid viscosity too high External tubing diameter too small Low line voltage | Replace with lower viscosity bath fluid. Replace with larger diameter tubing. Check and correct as required. |
| Unit does not heat | Insufficient fluid in reservoir Temperature set point too low Safety Set Temperature too low | Add fluid to reservoir. Increase temperature set point. Increase Safety Set temperature. |
| Insufficient heating | Insufficient circulation Low line voltage Ambient temperature too cool Excessive heat loss | See Insufficient circulation, above. Check and correct as required. Increase ambient temperature or relocate unit. Check for heat loss from external tanks and hoses; Check for vapor/heat loss from internal reservoir. |
| Temperature unstable | Insufficient circulation Debris or mineral build-up on pump, heater, or temperature sensor. | Check pump flow and operation. Clean as required. |
| Unit does not cool | Dust build up on air filter or condenser Blocked air ventilation screens Temperature set point is too high Excessive heat load Ambient air temperature too high (>35°C / 95°F) Low or high line voltage | Clean air filter and/or condenser as required. Remove blockages as required. Decrease temperature set point. Check that heat load does not exceed capacity of bath; correct as required. Decrease ambient air temperature. Check and correct as required. |

| Problem | Possible Causes | Corrective Action |
|--|---|--|
| Insufficient cooling | Dust build up on air filter or condenser Blocked air ventilation screens Temperature set point is too high Excessive heat load Ambient air temperature too high (>35°C / 95°F) Low or high line voltage | Clean air filter and/or condenser as required. Remove blockages as required. Decrease temperature set point. Check that heat load does not exceed capacity of bath; correct as required. Decrease ambient air temperature. Check and correct as required. |
| Unable to achieve low end extreme temperatures | Pump speed too high Incorrect bath fluid Insufficient insulation on external fluid lines Ambient air temperature too high (>35°C / 95°F) Low or high line voltage Dust build up on air filter or condenser Blocked air ventilation screens Excessive heat load | Reduce pump speed. Check that the fluid being circulated is capable of reaching the required temperature. Check external fluid lines for proper insulation. Decrease ambient air temperature as required. Check and correct as required. Clean air filter or condenser as required. Remove blockages as required. Check that heat load does not exceed capacity of bath; correct as required. |

Technical Information

Performance Specifications

| | | |
|------------------------------|----------------------------------|--------------------|
| Operating Temperature Range: | Model dependent; see table below | |
| Temperature Stability: | ±0.01C (±0.02°F) | |
| Pump Type: | Variable speed pressure/suction | |
| | <u>60Hz models</u> | <u>50Hz models</u> |
| Maximum Pressure: | 4.3 psi (0.30 bar) | 3.6 psi (0.25 bar) |
| Maximum Pressure Flow Rate: | 5.3 gpm (20.1 lpm) | 4.4 gpm (16.7 lpm) |
| Maximum Suction Flow Rate: | 3.9 gpm (14.7 lpm) | 3.2 gpm (12.2 lpm) |
| Heater Wattage: | 1100 watts | 2200 watts |

| Model Type | Reservoir Capacity | Temperature Range | Electrical Requirements | |
|---------------------------------------|--------------------|--|---------------------------|---------------------------|
| | | | 60Hz Units | 50Hz Units |
| AD07R-20 Refrigerating / Heating Bath | 7 liters | -20° to 200°C -4° to 392°F | 120V, 60Hz, 12A | 240V, 50Hz, 12A |
| AD07R-40 Refrigerating / Heating Bath | 7 liters | -40° to 200°C -40° to 392°F | 120V, 60Hz, 12A | 240V, 50Hz, 12A |
| AD7LR-20 Refrigerating / Heating Bath | 7 liters | -20° to 200°C -4° to 392°F | 120V, 60Hz, 12A | 240V, 50Hz, 12A |
| AD15R-30 Refrigerating / Heating Bath | 15 liters | -30° to 200°C -22° to 392°F | 120V, 60Hz, 13A | 240V, 50Hz, 13A |
| AD15R-40 Refrigerating / Heating Bath | 15 liters | -40° to 200°C -40° to 392°F | 120V, 60Hz, 13A | 240V, 50Hz, 13A |
| AD20R-30 Refrigerating / Heating Bath | 20 liters | -30° to 200°C -22° to 392°F | 120V, 60Hz, 13A | 240V, 50Hz, 13A |
| AD28R-30 Refrigerating / Heating Bath | 28 liters | -30° to 200°C -22° to 392°F | 120V, 60Hz, 13A | 240V, 50Hz, 13A |
| AD45R-20 Refrigerating / Heating Bath | 45 liters | -25° to 135°C -13° to 275°F | 208-240V, 50/60Hz, 12A | 208-240V, 50/60Hz, 12A |
| AD07H200 Heating Only Bath | 7 liters | Ambient +10° to 200°C Ambient +20° to 392°F | 120V, 60Hz, 10A | 240V, 50Hz, 10A |
| AD15H200 Heating Only Bath | 15 liters | Ambient +10° to 200°C Ambient +20° to 392°F | 120V, 60Hz, 10A | 240V, 50Hz, 10A |
| AD20H200 Heating Only Bath | 20 liters | Ambient +10° to 200°C Ambient +20° to 392°F | 120V, 60Hz, 10A | 240V, 50Hz, 10A |
| AD28H200 Heating Only Bath | 28 liters | Ambient +10° to 200°C Ambient +20° to 392°F | 120V, 60Hz, 10A | 240V, 50Hz, 10A |

| | |
|--------------------------|--------------------------|
| Environmental Conditions | Indoor use only |
| Maximum Altitude: | 2000 meter |
| Operating Ambient: | 5° to 35°C (41° to 95°F) |
| Relative Humidity: | 80%, non-condensing |
| Installation Category: | II |
| Pollution Degree: | 2 |
| Ingress Protection: | IP 31 |
| Climate Class: | SN |
| Software Class: | B |
| Output Waveform: | Sinusoidal |

Specifications subject to change without notice.

Reservoir Fluids

Depending on your needs, a variety of fluids can be used with your Circulator. No matter what bath fluid is selected, it must be chemically compatible with the reservoir and the materials in your Circulator. It must also be suitable for the desired temperature range.

| | | | |
|--|---|--|---|
|  | <p>WARNING: When using Class III flammable fluids per DIN 12876-1, the user must attach the following warning labels to the front of the unit so that they are well visible:</p> | | |
| | <p>Warning Label W09 Colors: Yellow/black</p> |  | <p>Danger Area. Attention! Observe instructions (operating manual, safety data sheet)</p> |
| <p>Mandatory Label M018 Colors: Blue/white</p> <p>or</p> <p>Semi S1-0701 Table A1-2 #9 Colors: Blue/white</p> |   | <p>Carefully read the user information prior to beginning operation. Scope: EU</p> <p>Carefully read the user information prior to beginning operation. Scope: NAFTA</p> | |

| | |
|---|---|
|  | <p>WARNING: Always use fluids that satisfy safety, health, and equipment compatibility requirements. Be aware of the chemical hazards that may be associated with the bath fluid used. Observe all safety warnings for the fluids used as well as those contained in the material safety data sheet.</p> |
|---|---|

For optimum temperature stability, the fluid's viscosity should be 50 centistokes (cSt) or less at its lowest operating temperature. This permits good fluid circulation and minimizes heating from the pump.

For temperatures from 10°C to 90°C, distilled water is recommended. For temperatures below 10°C, a mixture of laboratory grade ethylene glycol and water should be used. Do not use deionized water.

| | |
|---|--|
|  | <p>WARNING: DO NOT USE THE FOLLOWING LIQUIDS</p> <ul style="list-style-type: none"> • Automotive antifreeze with additives** • Hard tap water** • Deionized water with a specific resistance > 1 meg ohm • Concentrations of acids or bases • Solutions with halides: chlorides, fluorides, bromides, iodides or sulfur • Bleach (Sodium Hypochlorite) • Solutions with chromates or chromium salts • Glycerine • Syltherm fluids <p>** At temperatures above 40°C, additives or mineral deposits can adhere to the heater. If deposits are allowed to build up, the heater may overheat and fail. Higher temperatures and higher concentrations of additives will hasten deposit build up.</p> |
|---|--|

Application Notes

At a fluid's low temperature extreme:

- The presence of ice or slush adversely affects temperature stability.
- A viscosity above 10 centistokes adversely affects temperature uniformity.
- A high fluid viscosity and high pump speed adds heat to the fluid being pumped.

At a fluid's temperature above ambient without refrigeration:

- If your set point temperature is less than 15°C above the ambient temperature, the viscosity of the fluid should be 10 centistokes or less to minimize friction heating of the fluid.
- Heat loss should be encouraged by uncovering the fluid and lowering the pump speed.

At fluid's high temperature extreme:


- Heat loss from vapor adversely affects temperature stability.
- To prevent the accumulation of vapors inside the room, the reservoir may need to be placed in a fume hood.
- Use a cover and/or floating hollow balls to help prevent heat and vapor loss.
- Replenish fluid lost from vapor frequently.

Tubing and Fitting Temperature Ranges

| Material | Temperature Range |
|------------------------------|-------------------|
| Buna N tubing | -40° to 120°C |
| Viton® tubing | -32° to 200°C |
| Braided Teflon® lined tubing | -50° to 225°C |
| Stainless steel fittings | -45° to 225°C |
| Nylon fittings | -40° to 90°C |
| Brass fittings | -40° to 80°C |

RS232/RS485 Configuration

Selecting RS232 or RS485 Communication — Your Circulator may be configured for either RS232 or addressable RS485 serial communication. The selection is made on the COM sub-menu (see *Normal Operation, Selecting the Remote Communication and Control Protocol*).



To Access: Press the  key until COM and the currently selected protocol appears on the display.

To Change: To configure the Circulator for RS232 communication, touch the down arrow; to configure the Circulator for RS485 communication, touch a point between the down arrow and the word SET.

To Accept: Press . You will be prompted to select the communications baud rate.

Selecting the Communications Baud Rate — This sub-menu allows you to select the speed at which your Circulator will transmit data. The baud rate setting may be 1200, 2400, 4800, 9600, 19200, 38400, 57600, or 115200; the factory default is 9600.



NOTE: When using RS232/RS485 communication, the Circulator's baud rate must match that of the connected device.



Touch up arrow for fastest baud rate
Touch down arrow for slowest baud rate

To Change: To select the 1200 baud rate setting, touch the bottom of the scroll bar; to select the 115200 setting, touch the top of the scroll bar. Rates between these two extremes are selected by touching the corresponding area of the scroll bar.

To Accept: Press . If you have selected RS232 communication, the Timer sub-menu will appear; if you have selected RS485 communication, you will be prompted to enter the desired RS485 address.

Entering the RS485 Address — When RS485 has been selected as the serial communications protocol, you will be prompted to enter the RS485 address after the baud rate has been selected. An address from 001 to 255 may be entered.



Touch and hold up/down arrows
to make large changes
Touch and release to make minor
changes

To Change: Touch the touch scroll bar or touch the arrow symbols until the desired RS485 address is displayed.

To Accept: Press , , , or allow the LCD to return to the main operational display (approximately 10 seconds).

RS232/RS485 Communications



CAUTION: Always turn electrical power to the Circulator OFF before making a connection to the serial (DB9) port.

Serial Connector — A DB9 connector is provided on the back panel of the Controller for RS232/RS485 data communication.

| Pin | RS232 | RS485 |
|-----|-------|-------|
| 1 | | |
| 2 | TX | A |
| 3 | RX | |
| 4 | DTR | |
| 5 | GND | GND |
| 6 | DSR | |
| 7 | RTS | |
| 8 | CTS | B |
| 9 | | |

Serial Communications Protocol — The Controller uses the following serial communications settings:

Data bits — 8

Parity — None

Stop bits — 1

Flow control — None

Baud rate — Selectable (Controller/PC baud rates must match). 57600 is recommended.

Communications Commands — RS232 commands must be entered using the command only; RS485 commands must be entered using @ and the RS485 address (E.g., @001) followed by the actual command. All commands must be entered in the exact format shown. Do not send a [LF] (line feed) after the [CR] (carriage return). Be sure to follow character case exactly.

A response followed by an exclamation point (!) indicates that a command was executed correctly. A question mark (?) indicates that the Controller could not execute the command (either because it was in an improper format or the values were outside the allowable range). A response must be received from the Controller before another command can be sent. All responses are terminated with a single [CR].



IMPORTANT INFORMATION REGARDING SERIAL COMMANDS ENTERED OVER ETHERNET VIA USER DATAGRAM PROTOCOL (UDP): Only Read commands will be allowed when entering serial commands over Ethernet unless (1) Ethernet is configured as Unlocked and (2) the proper password is included as part of the serial command. See *Configuring the Ethernet Connection*.

Set commands entered over Ethernet must include Pxxx immediately following the serial command. For example, to change set point, enter the following command: **SSiii.iPxxx [CR]** where:

SS = Set set point

iii.ii = set point value

P = password identifier

xxx = password

[CR] = carriage return

| Command | Format | Values | Return Message |
|---------------------------------|--------------|---|----------------|
| Set Command Echo | SEi[CR] | Echo: i = 1 No Echo: i = 0 | ![CR] |
| Set Set Point | SSiii.ii[CR] | i = any integer from 0-9 | ![CR] |
| Set On Off | SOi[CR] | On: i = 1 Off: i = 0 | ![CR] |
| Set High Alarm | SHiii[CR] | i = any integer from 0-9 | ![CR] |
| Set Low Alarm | SLiii[CR] | i = any integer from 0-9 | ![CR] |
| Set Pump Speed | SMi[CR] | i = any integer from 5-100 in increments of 5 | ![CR] |
| Set Restart Power Status | SWi[CR] | Restart: i = 1 Standby: i = 0 | ![CR] |
| Set Internal / External Control | SJi[CR] | External: i = 1 Internal: i = 0 | ![CR] |
| Read Set Point Temperature | RS[CR] | | iii.ii[CR] |
| Read Units of Operation | RU[CR] | | C[CR] or F[CR] |
| Read Internal Temperature | RT[CR] | | iii.ii[CR] |
| Read External Temperature | RR[CR] | | iii.ii[CR] |
| Read Operating Status | RO[CR] | Running: i = 1 Standby: i = 0 | i[CR] |
| Read High Alarm Setting | RH[CR] | | iii[CR] |
| Read Low Alarm Setting | RL[CR] | | iii[CR] |
| Read Pump Speed | RM[CR] | i = any integer from 5-100 | i[CR] |
| Read Power Status | RW[CR] | | i[CR] |
| Read Alarm Status | RF[CR] | No Faults: i = 0 Fault: i = 1 | i[CR] |
| Read the Auto-Cool Set Point | RA[CR] | | ii[CR] |
| Read the Firmware Version | RB[CR] | | viii[CR] |

USB Data Logging

Data is stored in a CSV file that can be read in spreadsheet programs such as Microsoft Excel[®]. The data is output in the following format:

| MIN.SEC | [V] | UNIT | SET | P1 | P2 |
|---------|---------|------|-----|-------|-------|
| 4 | [v0013] | C | 10 | 15.19 | 14.89 |
| 4 | [v0013] | C | 10 | 15.18 | 14.88 |
| 4.22 | [v0013] | C | 10 | 15.03 | 14.71 |
| 4.23 | [v0013] | C | 10 | 15.02 | 14.7 |
| 4.24 | [v0013] | C | 10 | 15 | 14.69 |
| 4.25 | [v0013] | C | 10 | 15 | 14.69 |
| 4.26 | [v0013] | C | 10 | 14.99 | 14.68 |
| 4.27 | [v0013] | C | 10 | 14.98 | 14.67 |

MIN.SEC = Minutes and seconds.

V = Firmware version.

UNIT = Selected temperature unit.

SET = Temperature set point

P1 = Fluid temperature as measured by the Internal probe

P2 = Fluid temperature as measured by the External probe (if connected)



NOTE: The time stamp on the first line of data represents the time when data began being logged. The time stamp will indicate 0.01 if data logging was enabled after a flash drive or computer was connected to the USB port. It will show a later time if data logging was enabled before the storage device was connected. In the example shown above, the flash drive was connected to the USB port 4 minutes after data logging was enabled.

USB B Setup, Monitoring, and Control

The Advanced Digital Temperature Controller can be monitored and controlled using a personal computer connected to its USB B port.

Initial Setup

1. Plug the USB A connector of a USB A to USB B cable into an open USB port on the personal computer and plug the USB-B connector into the corresponding port on the Temperature Controller.
2. The computer operating system will automatically detect a new device and ask to install the drivers for the device. The drivers are located on the Resource Disk that came with your Temperature Controller.
3. Place the Resource Disk in the computers CD drive and install the drivers.
4. Determine the identity of the communications port your computer assigned to the Temperature Controller.

Monitoring and Control

When all connections have been made and drivers installed, you can begin monitoring and control the Temperature Controller using a terminal program and the active serial communications commands.

Configuring the Ethernet Connection



NOTE: The Advanced Digital Controller's supports Ethernet communication using User Datagram Protocol (UDP).

The Advanced Digital Controller may be connected directly to a laptop or desktop computer via its Ethernet connection or indirectly via your facility's wired or wireless network. Either type of connection enables you to control and/or monitor the operation of your circulating bath using User Datagram Protocol (UDP). The Controller's active serial communications commands are used to retrieve and/or change operational information.

Direct Computer to Controller Configuration



NOTE: If your computer is already connected to a wired network, you will not be able to directly connect it to your circulating bath's Temperature Controller.

1. Open the utility that displays your computer's available network connections.
2. Select an available wired connection.
3. Enter the following properties for the selected connection:
 - A. Internet protocol: TCP/IP
 - B. IP address: any three sets of identical numbers (e.g., 111.111.111) followed by a single digit (e.g., 111.111.111.5).
 - C. Subnet mask: 255.255.255.0
4. Enter the IP address into your Controller as outlined in *Fixed IP Addressing*, below. Please note that the first three series of numbers in IP address on the Controller must match those entered in Step 3B, above. The last number does not have to match.

Wired or Wireless Network Configuration

The Advanced Digital Controller supports both dynamic (DHCP) and fixed IP configuration. With either type of IP configuration, it may be necessary to work with your IT department to ensure that both the network and your Controller are configured properly.



Selects Ethernet

To Access: Press the  key until COM and the currently selected protocol appears on the display.

To Change: Touch the up arrow to select Ethernet.

To Accept: Press . A screen allowing you to select whether only Read commands will be allowed or if both Read and Set commands will be allowed will appear.



SET LOCKED – Read Commands Only Permitted



UNLOCK – Set and Read Commands Permitted

| | |
|---|---|
|  | <p>IMPORTANT: When Unlock is selected, a three digit numeric field will also appear. This field serves as a password that must be entered along with Set commands in order for the command to be executed.</p> <p>For example, the serial command required to change set point via Ethernet is: SSiii.iiPxxx [CR] where:</p> <ul style="list-style-type: none"> SS = Set set point iii.ii = set point value P = password identifier xxx = password [CR] = carriage return |
|---|---|

To Change: Use the touch scroll bar to change the displayed selection and, if Unlock is selected, also enter the password that will be used when executing Set commands over Ethernet. You may enter any numerical value from 100 to 255 as your password.

To Accept: Press . A screen displaying Fixed or dynamic (DHCP) Ethernet will appear.



To Change: Touch the down arrow to select Fixed IP; touch the up arrow to select DHCP IP.

To Accept: Press . The IP address for the selection will appear. If Fixed is selected, you will be allowed to change the IP address (see *Fixed IP Addressing*). If DHCP is selected, you will be allowed to view, but not change, the IP address (see *Dynamic IP Addressing*).

Fixed IP Addressing — If your Controller will be connected directly to a computer or your network uses fixed IP addressing, you must manually enter the IP address that will be used into your Controller. This is done as follows:

1. Select Fixed IP, select SET LOCKED or UNLOCKED, and then press . The IP address screen will be displayed.



To Change: Press **SET**. The last digit and decimal point in the first field of the IP address will begin flashing. Use the scroll bar to enter the correct IP address for that field and press **SET** to advance to the next address field. Repeat until the desired values for all IP address fields have been entered.

To Accept: Press **Menu**. The network mask screen will appear.

2. Enter the network mask address.



To Change: Press **SET**. The last digit and decimal point in the first field of the network mask will begin flashing. Use the scroll bar to enter the correct network mask for that field and press **SET** to advance to the next field. Repeat until the desired values for all network mask fields have been entered.

To Accept: Press **Menu**. The router address screen will appear.

3. Enter the router address.



To Change: Press **SET**. The last digit and decimal point in the first field of the router address will begin flashing. Use the scroll bar to enter the correct router address for that field and press **SET** to advance to the next address field. Repeat until the desired values for all router address fields have been entered.

To Accept: Press **☰**.

Dynamic IP Addressing — If your network uses DHCP (Dynamic Host Configuration Protocol), the IP address, network mask address, and router address will be assigned automatically. Addressing is assigned as follows:

1. Configure Ethernet communication as outlined in Wired or Wireless previously.
2. Turn power to the Controller OFF using the power switch / circuit breaker on the rear of the unit.
3. Turn power to the Controller back ON.




NOTE: This power down sequence is only required the first time the Controller is configured for DHCP.

To view the assigned addresses, select ETH DHCP, select SET LOCKED or UNLOCKED, and then press **☰**. The IP address screen will appear. Continue pressing **☰** to scroll to the network mask and router address screens.



Press the SET button to advance to subsequent address screens.

Equipment Disposal (WEEE Directive)



This equipment is marked with the crossed out wheeled bin symbol to indicate that this equipment must not be disposed of with unsorted waste.

Instead, it is your responsibility to correctly dispose of your equipment at lifecycle-end by handing it over to an authorized facility for separate collection and recycling. It's also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect from health hazards the persons involved in the disposal and recycling of the equipment

For more information about where you can drop off your waste equipment, please contact your local dealer from whom you originally purchased this equipment.

By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

Service and Technical Support

Web Resources

Visit the VWR's website at www.vwr.com for:

- Complete technical service contact information
- Access to VWR's Online Catalogue, and information about accessories and related products
- Additional product information and special offers

Contact us For information or technical assistance contact your local VWR representative or visit. www.vwr.com.

Warranty

VWR International warrants that this product will be free from defects in material and workmanship for a period of two (2) years from date of purchase. If a defect is present, VWR will, at its option, repair, replace, or refund the purchase price of this product at no charge to you, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear.

For your protection, items being returned must be insured against possible damage or loss. This warranty shall be limited to the replacement of defective products. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

Local VWR Offices in Europe and Asia Pacific

Austria

VWR International GmbH
Graumanngasse 7
1150 Wien
Tel.: 01 97 002 0
Fax: 01 97 002 600
E-mail: info@at.vwr.com

Belgium

VWR International bvba
Researchpark Haasrode 2020
Geldenaaksebaan 464
3001 Leuven
Tel.: 016 385 011
Fax: 016 385 385
E-mail: customerservice@be.vwr.com

China

VWR International China Co., Ltd
Suite 1802 - 1803,
Xing Ye Bank Mansion, No 168,
168 Jiangning Road
Shanghai 200041, China
Tel.: +86- 21 521 388 22
Fax: +86- 21 521 33 933
E-mail: sales_china@vwr.com

Denmark

VWR - Bie & Berntsen
Transformervej 8
2730 Herlev
Tel.: 43 86 87 88
Fax: 43 86 87 90
E-mail: info@dk.vwr.com

Finland

VWR International Oy
Valimotie 9
00380 Helsinki
Tel.: +358 9 80 45 51
Fax: +358 9 80 45 52 00
E-mail: info@fi.vwr.com

France

VWR International S.A.S.
Le Périgares – Bâtiment B
201, rue Carnot
94126 Fontenay-sous-Bois cedex
Tel.: 0 825 02 30 30 (0,15 EUR
TTC/min)
Fax: 0 825 02 30 35 (0,15 EUR
TTC/min)
E-mail: info@fr.vwr.com

Germany

VWR International GmbH
Hilpertstrasse 20a
D - 64295 Darmstadt
Tel.: 0180 570 20 00*
Fax: 0180 570 22 22*
E-mail: info@de.vwr.com
*0,14 €/Min. aus d. dt. Festnetz,
Mobilfunk max. 0,42 €/Min.

Hungary

VWR Spektrum-3D Kft.
Simon László u. 4.
4034 Debrecen
Tel.: (52) 521-131
Fax: (52) 470-069
E-mail: info@spektrum-3d.hu

India

VWR Lab Products Pte Ltd
2nd Floor, Front Wing, 135/12,
Brigade Towers
Brigade Road
Bangaluru 560025 India
Tel: +91-2522-277876 (Mumbai)
Tel: +91-80-41117124 (Bangalore)
Fax +91-80-41117120
E-mail: vwr_india@vwr.com

Ireland / Northern Ireland

VWR International Ltd / VWR
International (Northern Ireland) Ltd
Orion Business Campus
Northwest Business Park
Ballycoolin
Dublin 15
Tel.: 01 88 22 222
Fax: 01 88 22 333
E-mail sales@ie.vwr.com

Italy

VWR International s.r.l.
Via Stephenson 94
20157 Milano (MI)
Tel.: 02 332 03 11
Fax: 800 152 999
E-mail: info@it.vwr.com

The Netherlands

VWR International B.V.
Postbus 8198
1005 AD Amsterdam
Tel.: 020 4808 400
Fax: 020 4808 480
E-mail: info@nl.vwr.com

Norway

VWR International AS
Haavard Martinsens vei 30
0978 Oslo
Tel.: 02290
Fax: 815 00 940
E-mail: info@no.vwr.com

Poland

Labart Sp. z o.o.
A VWR International Company
Limbowa 5
80-175 Gdansk
Tel.: 058 32 38 210
Fax. 058 32 38 205
E-mail: labart@labart.pl

Portugal

VWR International - Material de
Laboratório, Lda
Edifício Neopark
Av. Tomás Ribeiro, 43- 3 D
2790-221 Carnaxide
Tel.: 21 3600 770
Fax: 21 3600 798/9
E-mail: info@pt.vwr.com

Singapore

VWR Singapore Pte Ltd
18 Gul Drive
Singapore 629468
Tel: +65 6505 0760
Fax: +65 6264 3780
E-mail: sales@sg.vwr.com

Spain

VWR International Eurolab S.L.
C/ Tecnología 5-17
A-7 Llinars Park
08450 - Llinars del Vallès
Barcelona
Tel.: 902 222 897
Fax: 902 430 657
E-mail: info@es.vwr.com

Sweden

VWR International AB
Fagerstagatan 18a
163 94 Stockholm
Tel.: 08 621 34 00
Fax: 08 621 34 66
E-mail: info@se.vwr.com

Switzerland

VWR International AG
Lerzenstrasse 16/18
8953 Dietikon
Tel.: 044 745 13 13
Fax: 044 745 13 10
E-mail: info@ch.vwr.com

UK

VWR International Ltd
Customer Service Centre
Hunter Boulevard
Magna Park
Lutterworth
Leicestershire
LE17 4XN
Tel.: 0800 22 33 44
Fax: 01455 55 85 86
E-mail: uksales@uk.vwr.com