

## Basket Hangers

### Basket Hangers

The basket hangers are made of stainless steel. The assembly is notched for holding USP disintegration baskets. These basket hangers are for use with the Standard Disintegration Apparatus (single-basket, three-basket, six-basket, or 12-basket systems).



Basket Hanger

### Basket Hangers

Item Description	Part Number
Basket Hanger, 1-unit	<b>37-1000</b>
Basket Hanger, 3-unit	<b>37-1005</b>
Basket Hanger, 6-unit	<b>37-1015</b>



Basket Suspension Rod

## Basket Suspension Rod

The basket suspension rod is made of type 316 stainless steel. It has an adjustable clip for easy setup.

### Basket Suspension Rod

Item Description	Part Number
Suspension Rod for Standard Disintegration Apparatus	<b>37-1040</b>

## Disintegration Accessories



### VK Verified

- Individual certificate
- Critical dimensions measured
- Results are traceable

Varian offers several individually verified disintegration accessories shown below. For a nominal charge, each of the accessories listed is measured to ensure conformance. A certificate is included to provide the actual measured values, the type of measuring device used, and its traceability.

Should you require verification of another component—give us a call!

### Basket Rack Assembly

Item Description	Part Number
Basket Rack Assembly for VK 100, 6-tube, 10 mesh, 1905 $\mu\text{m}$ stainless steel screen, fits VK 100, includes individual certificate	<b>37-3001V</b>
Basket, 6-tube, 10 mesh, 1905 $\mu\text{m}$ stainless steel screen for standard disintegration unit, includes individual certificate	<b>37-2001V</b>



Standard 6-tube Basket  
for Automated Unit

## Basket Rack Assemblies

### Standard 6-tube

The standard 6-tube basket rack assembly is designed to meet or exceed current USP specifications. This assembly has 10 mesh type 316 stainless steel screens and 80 mm length tubes. The tubes have a 25 mm outside diameter.

### Standard 6-tube with 10 mesh, 1905 $\mu\text{m}$ Stainless Steel Screen

Item Description	Part Number
For use with Standard Disintegration Apparatus	<b>37-2001</b>
For use with Automated Disintegration Tester, VK 100	<b>37-3001</b>



Standard 6-tube Basket  
for Automated Unit

### Special 6-tube Assembly

The special 6-tube assembly uses 20 mesh screens made of type 316 stainless steel. The tubes are 80 mm in length and have a 25 mm outside diameter.

### Special 6-tube with 20 mesh, 864 $\mu\text{m}$ Stainless Steel Screen

Item Description	Part Number
For use with Standard Disintegration Apparatus	<b>37-2010</b>
For use with Automated Disintegration Apparatus, VK 100	<b>37-3010</b>



Nylon 6-tube Basket  
for Automated Unit

### Nylon 6-tube Assembly

The Nylon 6-tube basket assembly is designed to meet current USP specifications. All components are Nylon except the 10 mesh type 316 stainless steel screen. The tubes are 80 mm in length and have a 25 mm outside diameter.

### Nylon 6-tube Basket, with 10 mesh, 1905 $\mu\text{m}$ Stainless Steel Screen

Item Description	Part Number
For use with Standard Disintegration Apparatus	<b>37-2020</b>
For use with Automated Disintegration Apparatus, VK 100	<b>37-3020</b>

**3-tube Basket for Bolus and Large Tablets**

The 3-tube basket for bolus and large tablets has 10 mesh type 316 stainless steel screens. The tubes are 80 mm in length and have a 38 mm outside diameter.

**3-tube Basket with 10 mesh, 1905 µm Stainless Steel Screen**

Item Description	Part Number
For use with Standard Disintegration Apparatus	<b>37-2030</b>
For use with Automated Disintegration Apparatus, VK 100	<b>37-3030</b>



3-tube Basket  
for Automated Unit

**Single-tube Basket for Bolus and Large Tablets**

The single-tube basket for bolus and large tablets has a 10 mesh type 316 stainless steel screen. The tube is 92 mm in length and has a 44 mm outside diameter.

**Single-tube Basket, with 10 mesh, 1905 µm Stainless Steel Screen**

Item Description	Part Number
For use with Standard Disintegration Apparatus	<b>37-2040</b>
For use with Automated Disintegration Apparatus, VK 100	<b>37-3040</b>



Single-tube Basket for  
Standard Unit

**6-tube Basket Assembly with 40 mesh screens and 40 mesh cover**

This basket assembly is designed to meet current USP specifications. It includes a basket and special 6-hole stainless steel cover with 40 mesh stainless steel screen. The tubes are 80 mm in length and have a 25 mm outside diameter.

**6-tube Basket, with 40 mesh Stainless Steel Screen**

Item Description	Part Number
For use with Standard Disintegration Apparatus	<b>37-2050</b>
For use with Automated Disintegration Apparatus, VK 100	<b>37-3050</b>




6-tube Basket  
for Standard Unit

## Basket Rack Assembly Accessories

### Covers for Basket Assemblies

Covers are designed to meet current USP specifications. One stainless steel cover is required per basket assembly. These covers fit both the standard and automated (VK 100) Apparatus basket assemblies.



### Stainless Steel Cover with Stainless Steel Screen

Item Description	Part Number
UPS 10 mesh, 1905 µm for 6-tube basket (screen type) (A)	<b>37-2100</b>
10 mesh, 1905 µm for 3-tube basket (3-hole type) (B)	<b>37-2105</b>
40 mesh for 6-tube basket (six-hole type) (C)	<b>37-2110</b>

### Replacement Screens for Basket Assembly

Replacement screens fit both standard and automated (VK 100) apparatus basket assemblies. Screens are sold in sets of six.

### Stainless Steel Screen

Item Description	Part Number
USP 10 mesh, 1905 µm for 6-tube basket	<b>37-2200</b>
20 mesh, 864 µm for 6-tube basket	<b>37-2205</b>
40 mesh for 6-tube basket	<b>37-2210</b>

### Replacement Glass Tubes

Replacement tubes fit both standard and automated (VK 100) apparatus basket assemblies.

### Replacement Glass Tubes

Item Description	Part Number
USP-type Tube for 6-tube basket rack 80 mm length, 25 mm outside diameter, set of 6	<b>37-5000</b>
Tube for 3-tube basket rack 80 mm length, 38 mm outside diameter, set of 3	<b>37-5005</b>
Tube for single-tube basket rack 92 mm length, 44 mm outside diameter	<b>37-5010</b>

### Depth Gauge

This depth gauge is designed to verify that the disintegration baskets do not dip lower than within 25 mm of the bottom of the disintegration water bath.

### Depth Gauge

Item Description	Part Number
Depth gauge, 25 mm	<b>37-9000</b>

### Plastic Disks

Plastic disks fit both standard and automated (VK 100) apparatus basket assemblies.



### Plastic Disks

Item Description	Part Number
USP Fluted Disks for 6-tube assembly, set of 6	<b>37-4000</b>
Disks for 3-tube assembly, set of	<b>37-4010</b>
Disks for single-tube assembly	<b>37-4020</b>

### Plunger Assemblies

Plunger assemblies fit both standard and automated (VK 100) apparatus basket assemblies.



### Plunger Assemblies

Item Description	Part Number
Type for 6-tube assembly	<b>37-4100</b>
Special type for 3-tube assembly	<b>37-4110</b>
Canadian type for 6-tube assembly	<b>37-4120</b>